

Filter & Blotting Papers



Ahlstrom - Munksjö
Filter Paper

Qualitative, Standard Grades (for gravity filtration)

Slow flow (4 mL/min), 1.5 µm retention, Grade 610

(Whatman 5; S&S 602)

9 cm, Pk/100. 6100-0900	\$9.85
12.5 cm, Pk/100. 6100-1250	\$12.40
15 cm, Pk/100. 6100-1500	\$14.95
20.5 cm, Pk/100. 6100-2050	\$23.60
24 cm, Pk/100. 6100-2400	\$33.55

Medium flow (30 mL/min), 2 µm retention, Grade 642

(Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)

9 cm, Pk/100. 6420-0900	\$7.65
15 cm, Pk/100. 6420-1500	\$13.90

Medium flow (55 mL/min), 2.5 µm retention, Grade 601

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

4.25 cm, Pk/100. 6010-0425	\$5.90
9 cm, Pk/100. 6010-0900	\$7.90
12.5 cm, Pk/100. 6010-1250	\$11.45
15 cm, Pk/100. 6010-1500	\$13.95
18.5 cm, Pk/100. 6010-1850	\$20.20
20.5 cm, Pk/100. 6010-2050	\$21.80
33 cm, Pk/100. 6010-3300	\$58.85

Medium flow (60 mL/min), 6 µm retention, Grade 613

(Whatman Student; S&S 595, 0903)

9 cm, Pk/100. 6130-0900	\$6.75
15 cm, Pk/100. 6130-1500	\$11.80
20.5 cm, Pk/100. 6130-2050	\$17.75
33 cm, Pk/100. 6130-3300	\$39.40

Medium flow (200 mL/min), 10 µm retention, Grade 631

(Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)

20.5 cm, Pk/100. 6310-2050	\$22.00
33 cm, Pk/100. 6310-3300	\$59.45

Fast flow (235 mL/min), 25 µm retention, creped, Grade 615

9 cm, Pk/100. 6150-0900	\$6.95
12.5 cm, Pk/100. 6150-1250	\$10.15
15 cm, Pk/100. 6150-1500	\$12.20
20.5 cm, Pk/100. 6150-2050	\$18.30
33 cm, Pk/100. 6150-3300	\$40.60

Very Fast flow (360 mL/min), 35 µm retention, creped, Grade 617

9 cm, Pk/50. 6170-0900	\$6.30
12.5 cm, Pk/50. 6170-1250	\$8.50
15 cm, Pk/50. 6170-1500	\$9.90
20.5 cm, Pk/50. 6170-2050	\$15.10
33 cm, Pk/50. 6170-3300	\$32.35

Qualitative, Pleated (Fluted)

Slow flow (30 mL/min), 2 µm retention, Grade 542

(Whatman 2V)

12.5 cm, Pk/100. 5420-1250	\$68.05
18.5 cm, Pk/100. 5420-1850	\$77.05

Qualitative, Wet-Strengthened Grades (for vacuum filtration)

Medium flow (65 mL/min), 3 µm retention, Grade 909

(Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)

32 cm, Pk/100. 9090-3200	\$64.10
---------------------------------------	---------

Fast flow (245 mL/min), 27 µm retention, creped, Grade 962

(Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)

11 cm, Pk/100. 9620-1100	\$9.55
---------------------------------------	--------

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

Slow flow (4 mL/min), 1.5 µm retention, Grade 94

(Whatman 42; S&S 589/3)

9 cm, Pk/100. 0940-0900	\$27.50
15 cm, Pk/100. 0940-1500	\$47.40

Medium flow (30 mL/min), 2 µm retention, Grade 74

(Whatman 40; S&S 589/2, 589/4)

9 cm, Pk/100. 0740-0900	\$27.50
15 cm, Pk/100. 0740-1500	\$47.40

Fast flow (220 mL/min), 15 µm retention, Grade 55

(Whatman 54, 541; S&S 1505)

9 cm, Pk/100. 0550-0900	\$31.10
15 cm, Pk/100. 0550-1500	\$49.90

Blotting/Chromatography/Electrophoresis Papers

0.19 mm thick, 125 mm/30 min absorption, Grade 601

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

46 x 57 cm, Pk/100. 6018-4657	\$114.00
--	----------

0.34 mm thick, 140 mm/30 min absorption, Grade 238

(Whatman 3MM, S&S 598, 2316)

7 x 10 cm, Pk/100. 2388-0710	\$26.90
20 x 20 cm, Pk/100. 2388-2020	\$59.35
46 x 57 cm, Pk/100. 2388-4657	\$218.00

0.95 mm thick, 180 mm/30 min absorption, Grade 243

(Whatman 17CHR; S&S 2668)

11 x 14 cm, Pk/100. 2438-1114	\$25.65
20 x 20 cm, Pk/100. 2438-2020	\$47.55

2.48 mm thick, 240 mm/30 min absorption, Grade 320

(S&S 2727)

11 x 14 cm, Pk/50. 3208-1114	\$40.55
20 x 20 cm, Pk/50. 3208-2020	\$69.00



Smooth tips &
polypropylene grips

Forceps, Filter, 11.5 cm. 300-025	\$4.40
--	--------



Vinyl-coated lead
wire is 300 mm long.

Flask Weight, lead wire, 114 g. 183240000	\$33.50
--	---------