

Stereo Microscopes



420T-430PLL-10

National

420-430PLL-10 Stereo, Binocular, 10-40X Zoom..... \$568.00 ↓

- Widefield 10X/20 eyepieces; accept 23 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 10X-40X total magnification; working distance of 79 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual illuminators offer light from top only (LED with rheostat control and directional beam), bottom only (LED), or both top and bottom
- Post stand with locking support collar allows viewing specimens up to 93 mm tall
- Frosted glass and reversible black/white plastic stage plates included
- 24.2 cm(D) x 17.1 cm(W) x 43.3 cm(H)
- Lifetime mechanical and optical/3-year electrical warranty

420T-430PLL-10 Stereo, Trinocular, 10-40X Zoom..... \$657.00 ↓

- Trinocular version of 420-430PLL-10
- 2X SLR adapter & 0.4X C-mount adapter included

425-445PLL-10



National

425-445PLL-10 Stereo, Binocular, 7.5-50X Zoom..... \$912.00 ↓

- Widefield 10X/23 eyepieces; accept 24.95 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 7.5-50X total magnification; working distance of 113 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual 3W LED illuminators with rheostat control offer cool light from top, bottom, or both top and bottom
- Post stand has a large ergonomic base and locking support collar
- Frosted glass and reversible black/white plastic stage plates included
- Lifetime mechanical and optical/3-year electrical warranty

425T-445PLL-10 Stereo, Trinocular, 7.5-50X Zoom..... \$972.00 ↓

- Trinocular version of 425-445PLL-10




Meiji Access www.meijitechno.com
and call us to place an order


Stands



SMS16B



SMS20



SMS25

Please call regarding availability of items not listed

Stereo Boom Stands, Diagnostic Instruments

All boom stands require a focus mount

Standard, nylon bearings. **SMS16B** \$599.00

Ball Bearing

Medium-Duty. **SMS6B** \$759.00

Heavy-Duty. **SMS20-*** See note below

Articulating Arm. **SMS25** \$1,409.00

**Access www.spotimaging.com and call us with microscope model and boom stand version to place an order*



Meiji PBH Stand
with Adapter Sleeve

Stereo Dual Illumination Stands, 6V/10W, Meiji

Pole-Type Base, 84 mm focus mount. **PBH**..... \$563.00

Some stereo brands require an adapter sleeve

Stereo Stand Sleeve for 76 mm dia. bodies \$30.00

Stereo Stage Plates, 94.5 mm

Clear Glass. **MA569** \$45.00

Filter, Blue Frosted, 40 mm. **MA563/05** \$50.00

Black & White Plastic. **MA568** \$45.00