

Funnels (cont)

Standard Stem, Premium, polypropylene, Bel-Art
 Exclusive filter-step internal design provides speedier filtration; external ribs prevent air-lock.

55 mm Top ID, takes 9 cm filter, Pk/12. 146760055	\$47.75
65 mm Top ID, takes 11 cm filter, Pk/12. 146760065	\$49.60
75 mm Top ID, takes 12.5 cm filter, Pk/6. 146760075	\$25.65
100 mm Top ID, takes 18.5 cm filter, Pk/4. 146760100	\$32.60

Standard Stem, Economy, polypropylene

50 mm Top ID, each. 310FUNPP50	\$0.80
62 mm Top ID, each. 310FUNPP62	\$0.95
75 mm Top ID, each. 310FUNPP75	\$1.15
100 mm Top ID, each. 310FUNPP100	\$1.65

These 4 x 4" nonsterile gauze sponges are 12-ply.



Gauze Pads, 10 pk/200, Cs/2000. 44122000 \$71.00

Sani-Cloth Plus



Germicidal Wipes, Tub/160. 440Q89072 \$9.50

Gloves & Hand Protection


Cotton Gloves



Heat-Resistant Gloves



Nitrile Gloves



Gloves

Cotton

Protect optics and photographs from fingerprints.
 Can also be used as glove liners.

Small, Pk/12. GLV-190.10	\$4.95
Large, Pk/12. GLV-190.20	\$5.25

Heat-Resistant

Terry cloth gloves protect hands from high temperatures.

33 cm (13") Long, Pair. 132010000	\$34.20
47 cm (18.5") Long, Pair. 132010001	\$46.55

Latex, powder-free, 5 mil, fully textured	<u>1 bx</u>	<u>cs/10 bx</u>
X-Small, Bx/100. 73101	\$9.70	87.00
Small, Bx/100. 73103	\$9.70	87.00
Medium, Bx/100. 73105	\$9.70	87.00
Large, Bx/100. 73107	\$9.70	87.00
X-Large, Bx/100. 73109	—	87.00
XX-Large, Bx/100. 73111	—	87.00

Nitrile, powder-free, 5.9 mil, textured fingers, Halyard

Nitrile gloves are hypoallergenic and provide greater puncture and chemical resistance than latex or vinyl gloves.




	<u>1 bx</u>	<u>cs/10 bx</u>
X-Small, Bx/100. 73121	\$13.55	120.00
Small, Bx/100. 73123	\$13.55	120.00
Medium, Bx/100. 73125	\$13.55	120.00
Large, Bx/100. 73127	\$13.55	120.00
X-Large, Bx/90. 73129	\$13.55	120.00
XX-Large, Bx/90. 73131	—	120.00



Safely grip hot or cold labware with this silicone rubber hand protector partially lined with mini suction cups. Withstands temperatures from 260°C (500°F) to -57°C (-70°F).

Hand Protector, red, Hot Hand®. 380000000 \$26.50

Homogenizers & Tissue Grinder



Rapidly homogenize, emulsify, suspend and disrupt 0.1 mL to 250 mL samples. Includes a 5 x 50 mm flat-bottom probe & 7 x 50 mm saw tooth probe for use with micro tubes.

Homogenizer, Handheld. D1000..... \$978.00



Homogenize a variety of sample types, from soft tissue to bone, in 35 seconds or less.
 Benchmark BeadBlaster™ 24

Homogenizer, Benchtop, Microtube. D2400..... \$7,049.00

- Holds (24) 2.0 mL tubes, or (6) 5 mL tubes w/optional tube holder
- Speed adjusts from 4.00 to 7.00 m/s in 0.05 m/s increments
- Stores 50 programs in memory, each with a maximum of 10 cycles
- 2-year warranty

Homogenizer Prefilled Tubes

2.0 mL screw-capped tubes are prefilled with your choice of the following beads:

Glass, acid washed	
0.1 mm, Pk/50. D1031-01	\$105.00
0.5 mm, Pk/50. D1031-05	\$105.00
Zirconium Beads, acid washed/heat treated/nuclease-free	
1.5 mm, Pk/50. D1032-15	\$147.00
3.0 mm, Pk/50. D1032-30	\$147.00
Stainless Steel Beads, acid washed	
2.8 mm, Pk/50. D1033-28	\$147.00

Continued

Continued