

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



Wax Pencil

Markers, USA

Wax Pencil (China Marker), Bx/12. **02089** \$11.50



Sharpie Retractable, permanent

Fine Tip	<u>ea</u>	<u>bx/12</u>
Black. 36701	---	\$29.65
Blue. 36703	\$2.85	29.65
Red. 36702	\$2.85	29.65
Ultra-Fine Tip	<u>ea</u>	<u>bx/12</u>
Black. 1735790	\$2.85	29.65
Blue. 1735792	---	29.65

Mortar and Pestle Sets

Glass, **USA**

59 mL (2 oz) Capacity. 826-2	\$16.05
118 mL (4 oz) Capacity. 826-3	\$16.90
237 mL (8 oz) Capacity. 826-4	\$20.00



Porcelain, outside wall and rim glazed

30 mL Capacity, 63 mm top OD. LPC-154	\$4.70
80 mL Capacity, 80 mm top OD. LPC-154-1	\$6.30
100 mL Capacity, 105 mm top OD. LPC-155	\$7.75
320 mL Capacity, 135 mm top OD. LPC-156	\$10.55
500 mL Capacity, 170 mm top OD. LPC-157	\$15.30

Petri Dishes

PS (polystyrene)

90 x 15 mm, economy Pk/20, nonsterile. 310150010	\$5.00
35 x 10 mm, vented lid, stacking ring, Simport Cs/500, sterile sleeve/20. D210-15	\$136.00
60 x 15 mm, stackable, Kord-Valmark Cs/500, sterile sleeve/20. 2901	\$96.00
100 x 15 mm, Kord-Valmark Slippable (rimless lid for automation)	
Standard, Cs/500, sterile sleeve/25. 2917	\$105.00
ISO Mark, Cs/500, sterile sleeve/25. 2910P	\$105.00
★ Stackable	
Standard, Cs/500, sterile sleeve/25. 2900	\$109.00
100 x 25 mm, stackable, Kord-Valmark Cs/340, sterile sleeve/17. 2906	\$109.00
Glass (borosilicate), reinforced rims, economy	
90 x 15 mm, each. 31015003	\$2.50 ↓
100 x 15 mm, Pk/10. 31015004	\$29.20



Measure pH with 1 pH unit accuracy from 1-14. Package contains 100 paper strips, 2 x 1/4", in a plastic vial with color chart.

pH Paper, pH 1-14, Pk/100. **600PHP100** \$3.75

Pipetting



Amber Latex



Natural Rubber

Natural Rubber w/Flange

Dropper Bulbs

Amber Latex


1 mL Capacity, 43 mm long, Pk/12. 20621A	\$6.50
1 mL Capacity, 43 mm long, Pk/72. HS20622A	\$33.95
2 mL Capacity, 52 mm long, Pk/12. 20621B	\$6.50
2 mL Capacity, 52 mm long, Pk/72. HS20622B	\$33.95

Natural Rubber

1 mL Capacity, 30 mm long, black, Pk/12. 200-1B	\$3.55
2 mL Capacity, 32 mm long, red, Pk/12. 200-2B	\$5.95


Natural Rubber w/Flange, blue

1 mL Capacity, 29 mm long, Pk/24. 200-1C	\$12.15
---	---------



glass w/rubber bulb; 75 mm (3")

Medicine (Eye) Dropper, Pk/12. **97-4301** \$3.60



Pipet Fillers

3-Valve, rubber, black. 207B	\$14.70
Wheel-Type, Pipette Pump™	
Fits Pipets Up To 2 mL, blue. 378970000	\$25.65
Pipette Pump Chuck (replacement). 937890070	\$8.90
Fits Pipets Up To 10 mL, green. 378980000	\$26.75
Pipette Pump Chuck (replacement). 937890080	\$9.75
Fits Pipets Up To 25 mL, red. 378990000	\$36.85
Pipette Pump Chuck (replacement). 937890090	\$9.75
Wheel-Type, Pipette Pump™ II, quick-release	<u>ea</u> <u>pk/12</u>
Fits Pipets Up To 2 mL, blue. 379111002	\$32.80 353.00
Fits Pipets Up To 10 mL, green. 379111010	\$32.80 353.00
Fits Pipets Up To 25 mL, red. 379111025	\$43.55 469.00

Continued



Package includes: Diamond Serologic portable pipette controller; two rechargeable batteries; battery charger; three PTFE membrane filters – (1) 0.20 µm & (2) 0.45 µm (one installed); desktop stand and wall mounting bracket.

Pipet Fillers (cont)

Motorized, rechargeable
Diamond Serologic (1-year warranty). **3380**..... \$262.00



Dimensions: 29 x 8.8 x 40 cm

Pipet Rack, magnetized, blue, Heathrow. **HS20612M** \$75.15

Pipets, Mohr (Measuring), Class A, individually certified


5 mL, (±0.020 mL), Pk/5. **600PM7058H**..... \$66.35
10 mL, (±0.030 mL), Pk/5. **600PM7058I** \$66.35
25 mL, (±0.050 mL), Pk/5. **600PM7058J** \$98.80



Dispense fluids accurately and safely with these color-coded and plugged pipets. Negative graduations provide extra volume (not available on 25 mL). RNase-/DNase-/Pyrogen-free, and noncytotoxic. Calibrated “To Deliver” (TD).

Pipets, Serological, polystyrene, UniPlast™

Nonsterile
5 mL, blue band, pk/25, Cs/375. **1755**..... \$81.15
10 mL, orange band, pk/25, Cs/250. **1775** \$55.60
25 mL, red band, pk/10, Cs/100. **1784** \$51.70
Sterile Pk/1 (individually paper/plastic wrapped)
1 mL, yellow band, Cs/500. **1700** \$81.45
2 mL, green band, Cs/500. **1720** \$87.30
5 mL, blue band, reverse grad., Cs/250. **1740** \$64.45
10 mL, orange band, reverse grad., Cs/200. **1760** \$56.20
25 mL, red band, reverse grad., Cs/100. **1780** \$60.75
50 mL, black band, reverse grad., Cs/100. **1820** \$152.00
Sterile Bg/25
1 mL, yellow band, Cs/1000. **1710** \$149.00
10 mL, orange band, Cs/250. **1770** \$66.15
25 mL, red band, Cs/100. **1782** \$53.80



Pipets, Transfer, LDPE (low-density polyethylene)

Graduated, Nonsterile, 15 cm long
5 mL capacity, 1 mL graduated, 25 drops/mL
#1 Bx/500. **PTP-02** \$17.50
Cs/10 Bx. **PTP-02** \$168.00
7 mL capacity, 3 mL graduated, 14 drops/mL
Bx/500. **PTP-01** \$16.45
Cs/10 Bx. **PTP-01** \$153.00
Graduated, Sterile Pk/1 (individually wrapped), 15 cm long
5 mL capacity, 1 mL graduated, Simport
Bx/400. **P200-521S** \$42.75
Cs/10 Bx. **P200-521S** \$392.00
Graduated, Sterile Pk/20, 15 cm long
5 mL capacity, 1 mL graduated, Bx/500. **HS206371B** \$37.50
Long Stem, Sterile Pk/1 (individually wrapped)
22.5 cm (9") long, 6 mL capacity, Cs/400. **139038** \$56.20 ↓
Long Stem, Nonsterile
30 cm (12") long, 23 mL capacity, Bx/100. **139050** \$20.50

Made of borosilicate glass in compliance with ASTM specification E-969, Class A requirements. Individually serialized and certified. Calibrated “To Deliver” (TD).



Pipets, Volumetric, Class A, individually certified

2 mL (±0.006 mL), orange band, Pk/5. **600PT7100B** \$67.40
5 mL (±0.010 mL), white band, Pk/5. **600PT7100C** \$70.80
10 mL (±0.020 mL), red band, Pk/5. **600PT7100D** \$80.90
20 mL (±0.030 mL), yellow band, Pk/5. **600PT7100E** \$114.00



These mini, air displacement pipettors provide an economical and accurate method for pipetting fixed volumes. Pipettors feature a two-step plunger and are fully autoclavable.

Pipettors, Single Channel

Fixed Volume, pocket-size, Diamond Jr.
5 µL, red. **3322-5** \$18.20
10 µL, pink. **3322-10** \$18.20
20 µL, green. **3322-20** \$18.20
50 µL, purple. **3322-50** \$18.20
100 µL, orange. **3322-100** \$18.20



Adjustable Volume, Diamond Advance™ (3-year warranty)
0.5-10 µL, 0.02 µL increments, red. **3341-10** \$219.00
5-50 µL, 0.1 µL increments, yellow. **3341-50** \$219.00
10-100 µL, 0.1 µL increments, blue. **3341-100** \$219.00
100-1000 µL, 1 µL increments, orange. **3341-1000** \$219.00

Continued