


Rulers, Meter Sticks & Tape Measure

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. 4461	\$0.95
White. 4460	\$0.55
Stainless Steel, 12 mm wide, USA. 605-065	\$5.60

Meter Sticks, lacquered wood

One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

Plain Ends. 605-050	\$4.90
Metal Ends. 605-051	\$6.15

Tape Measure, 152.4 cm/60", reinforced fiberglass. **11036**

One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.



Sharps Containers, red polypropylene

1 Quart, translucent locking lid. 440065	\$3.80
Dimensions: 4.25" x 4.25" x 6.25"(H)	
5.4 L, horizontal entry lid, Prevent™, each. 4402269	\$8.15
Dimensions: 12" x 4.75" x 11"(H)	
12 Gallon, hinged vertical entry lid, each. 440321391	\$28.40
Dimensions: 18.25" x 12.75" x 18.75"(H)	

Screw-Cap, 90 mL



Specimen Containers, polypropylene (PP)

Screw-Cap

90 mL, sterile, Cs/400. 840175	\$110.00
480 mL (16 oz), 83 mm (3 1/4") opening, Cs/24. 14510	\$71.95
960 mL (32 oz), 111 mm (4 3/8") opening, Cs/12. 14511	\$59.05

Snap-Cap, white HDPE bottom/white LDPE lid

480 mL (16 oz), Cs/100. 271016W	\$61.35
9.9 cm opening x 11.1 cm high	
1000 mL (34 oz), Cs/100. 271032W	\$90.20
11.7 cm opening x 13.7 cm high.	



Spot Plates, glazed white porcelain, 21 mm diameter wells

6-Well, each. LPC-161	\$2.35
12-Well, 115 x 95 x 7 mm, each. 3104187701	\$3.85

Stirring Bars



FDA approved Teflon® (PTFE) coating is a minimum of 1.6 mm (0.062") thick, except for micro bars which have a 1.32 mm (0.050") thick cover.

Stirring Bars, Teflon®, Bel-Art

Centered Sphere, Saturn Spinbar®

40 mm (1 1/2") Long. 371160002	\$25.70
50 mm (2") Long. 37116003	\$38.40

Micro, Spinbar®

White	
5 x 2 mm. 371190005	\$5.10
7 x 2 mm. 371190007	\$5.10
10 x 3 mm. 371190010	\$5.65
12.7 x 3 mm. 371190127	\$5.65
15 x 1.5 mm. 371190015	\$5.00

Red

7 x 2 mm. 371210013	\$5.65
8 x 1.5 mm. 371210001	\$5.65

Octagonal, molded-on pivot ring, Spinbar®

White	
12.7 x 8 mm (1/2 x 5/16"). 371101256	\$4.80
15.9 x 8 mm (5/8 x 5/16"). 371100058	\$4.80
25.4 x 9.5 mm (1 x 3/8"). 371100138	\$6.35
38.1 x 9.5 mm (1 1/2 x 3/8"). 371101128	\$7.50
50.8 x 9.5 mm (2 x 3/8"). 371100238	\$10.00
63.5 x 9.5 mm (2 1/2 x 3/8"). 371102128	\$11.55
76.2 x 12.7 mm (3 x 1/2"). 371100003	\$20.75

Red, 12.7 x 3.2 mm (1/2 x 1/8"), no pivot ring. **371090034**

Stirring Bar Assortment, octagonal, Pk/6. **371600000**

Includes two of each size: 25.4 x 8 mm (1 x 5/16"), 38.1 x 8 mm (1 1/2 x 5/16"), and 50.8 x 8 mm (2 x 5/16").

Polygonal, molded-on pivot ring, Spinbar®

12 x 4.5 mm (7/16 x 3/16"). **371220012**

Triangular, vanes w/cross-mounted magnets, Spinvane®

Fits 1 mL Vials. 371340000	\$9.95
Fits 3/5 mL Vials. 371350000	\$9.95






Stirring Bar Retainer, Spinbar® Restrainer. 377861000..... \$14.10


Stirring Bar Retrievers, magnetic, Spinbar®

30.5 cm (12") Long. **377720000** \$15.55
 45.7 cm (18") Long. **377730000** \$18.00



Prevent stirring bars from going down the sink. Strainer's four strong magnets capture stirring bars, but fluids drain easily. Blue.

Stirring Bar Sink Strainer, magnetic, Spinbar®. 377870000 \$30.30



Stirring Rods, glass, annealed ends

150 x 5 mm (6 x 3/16"), Pk/12. **600GSR006** \$4.40
 Rubber Policeman, fits 5 mm rod, Pk/12. **201-BB** \$6.00
 254 x 6 mm (10 x 1/4"), Pk/12. **600GSR010** \$6.25
 305 x 10 mm (12 x 3/8"), Pk/12. **600GSR012** \$23.40

Stoppers



Stoppers, rubber, black

Stopper Sizing Guide

Size #	Top (mm)	Bottom (mm)	Fits ID Openings (mm)
1	19	14	15-17
2	20	16	16-18.5
3	24	18	18-21
4	26	20	20-23
5	27	23	23-25
6	32	26	26-27
6½	34	27	27-31.5
7	37	30	30-34
8	41	33	33-37
9	45	37	37-41
10	50	42	42-46
12	64	54	54-59

All stoppers are 25 mm long. Stopper size #1 has a 4 mm hole; sizes #2 and larger have 5 mm holes.

Solid		<u>pk/5</u>	<u>lb</u>
Size 2, ~39 pcs/lb. 50020	\$3.15		13.35
Size 3, ~29 pcs/lb. 50030	\$4.25		13.35
Size 4, ~23 pcs/lb. 50040	\$5.35		13.35
Size 5, ~20 pcs/lb. 50050	—		13.35
Size 6, ~15 pcs/lb. 50060	—		13.35
Size 6.5, ~14 pcs/lb. 500650	—		13.35
Size 7, ~11 pcs/lb. 50070	—		13.35
Size 8, ~9 pcs/lb. 50080	—		13.35
1-Hole		<u>pk/5</u>	<u>lb</u>
Size 1, ~51 pcs/lb. 50011	\$2.40		13.35
Size 2, ~43 pcs/lb. 50021	\$2.85		13.35
★ Size 3, ~30 pcs/lb. 50031	\$4.10		13.35
Size 4, ~24 pcs/lb. 50041	\$5.10		13.35
Size 5, ~21 pcs/lb. 50051	—		13.35
Size 6, ~18 pcs/lb. 50061	—		13.35
Size 9, ~8 pcs/lb. 50091	—		13.35
Size 12, ~3 pcs/lb. 50121	—		13.35
2-Hole			<u>lb</u>
Size 6, ~17 pcs/lb. 50062			\$13.35
Size 7, ~12 pcs/lb. 50072			\$13.35
Size 8, ~9 pcs/lb. 50082			\$13.35
Size 9, ~7 pcs/lb. 50092			\$13.35
Size 10, ~5 pcs/lb. 50102			\$13.35
Assortment, 1 lb. 50012			\$17.35

Solid, 1-hole, and 2-hole stoppers in sizes 1-6.

Continued