




Stirring Bar Retainer, Spinbar® Restrainer. 377861000..... \$14.10


Stirring Bar Retrievers, magnetic, Spinbar®

30.5 cm (12") Long. **377720000** \$15.55
 45.7 cm (18") Long. **377730000** \$18.00



Prevent stirring bars from going down the sink. Strainer's four strong magnets capture stirring bars, but fluids drain easily. Blue.

Stirring Bar Sink Strainer, magnetic, Spinbar®. 377870000 \$30.30



Stirring Rods, glass, annealed ends

150 x 5 mm (6 x 3/16"), Pk/12. **600GSR006** \$4.40
 Rubber Policeman, fits 5 mm rod, Pk/12. **201-BB** \$6.00
 254 x 6 mm (10 x 1/4"), Pk/12. **600GSR010** \$6.25
 305 x 10 mm (12 x 3/8"), Pk/12. **600GSR012** \$23.40

Stoppers



Stoppers, rubber, black

Stopper Sizing Guide

Size #	Top (mm)	Bottom (mm)	Fits ID Openings (mm)
1	19	14	15-17
2	20	16	16-18.5
3	24	18	18-21
4	26	20	20-23
5	27	23	23-25
6	32	26	26-27
6½	34	27	27-31.5
7	37	30	30-34
8	41	33	33-37
9	45	37	37-41
10	50	42	42-46
12	64	54	54-59

All stoppers are 25 mm long. Stopper size #1 has a 4 mm hole; sizes #2 and larger have 5 mm holes.

Solid	pk/5	lb
Size 2, ~39 pcs/lb. 50020	\$3.15	13.35
Size 3, ~29 pcs/lb. 50030	\$4.25	13.35
Size 4, ~23 pcs/lb. 50040	\$5.35	13.35
Size 5, ~20 pcs/lb. 50050	—	13.35
Size 6, ~15 pcs/lb. 50060	—	13.35
Size 6.5, ~14 pcs/lb. 500650	—	13.35
Size 7, ~11 pcs/lb. 50070	—	13.35
Size 8, ~9 pcs/lb. 50080	—	13.35
1-Hole	pk/5	lb
Size 1, ~51 pcs/lb. 50011	\$2.40	13.35
Size 2, ~43 pcs/lb. 50021	\$2.85	13.35
★ Size 3, ~30 pcs/lb. 50031	\$4.10	13.35
Size 4, ~24 pcs/lb. 50041	\$5.10	13.35
Size 5, ~21 pcs/lb. 50051	—	13.35
Size 6, ~18 pcs/lb. 50061	—	13.35
Size 9, ~8 pcs/lb. 50091	—	13.35
Size 12, ~3 pcs/lb. 50121	—	13.35
2-Hole	lb	
Size 6, ~17 pcs/lb. 50062	\$13.35	
Size 7, ~12 pcs/lb. 50072	\$13.35	
Size 8, ~9 pcs/lb. 50082	\$13.35	
Size 9, ~7 pcs/lb. 50092	\$13.35	
Size 10, ~5 pcs/lb. 50102	\$13.35	
Assortment, 1 lb. 50012	\$17.35	

Solid, 1-hole, and 2-hole stoppers in sizes 1-6.

Continued