



Stirring Bar Retainer, Spinbar® Restrainer. 377861000..... \$14.10


Stirring Bar Retrievers, magnetic, Spinbar®

30.5 cm (12") Long. **377720000** \$15.55
 45.7 cm (18") Long. **377730000** \$18.00



Prevent stirring bars from going down the sink. Strainer's four strong magnets capture stirring bars, but fluids drain easily. Blue.

Stirring Bar Sink Strainer, magnetic, Spinbar®. 377870000 \$30.30



Stirring Rods, glass, annealed ends

150 x 5 mm (6 x 3/16"), Pk/12. **600GSR006** \$4.40
 Rubber Policeman, fits 5 mm rod, Pk/12. **201-BB** \$6.00
 254 x 6 mm (10 x 1/4"), Pk/12. **600GSR010** \$6.25
 305 x 10 mm (12 x 3/8"), Pk/12. **600GSR012** \$23.40

Stoppers



Stoppers, rubber, black

Stopper Sizing Guide			
Size #	Top (mm)	Bottom (mm)	Fits ID Openings (mm)
1	19	14	15-17
2	20	16	16-18.5
3	24	18	18-21
4	26	20	20-23
5	27	23	23-25
6	32	26	26-27
6½	34	27	27-31.5
7	37	30	30-34
8	41	33	33-37
9	45	37	37-41
10	50	42	42-46
12	64	54	54-59

All stoppers are 25 mm long. Stopper size #1 has a 4 mm hole; sizes #2 and larger have 5 mm holes.

Solid		pk/5	lb
Size 2, ~39 pcs/lb. 50020	\$3.15		13.35
Size 3, ~29 pcs/lb. 50030	\$4.25		13.35
Size 4, ~23 pcs/lb. 50040	\$5.35		13.35
Size 5, ~20 pcs/lb. 50050	—		13.35
Size 6, ~15 pcs/lb. 50060	—		13.35
Size 6.5, ~14 pcs/lb. 500650	—		13.35
Size 7, ~11 pcs/lb. 50070	—		13.35
Size 8, ~9 pcs/lb. 50080	—		13.35
1-Hole		pk/5	lb
Size 1, ~51 pcs/lb. 50011	\$2.40		13.35
Size 2, ~43 pcs/lb. 50021	\$2.85		13.35
★ Size 3, ~30 pcs/lb. 50031	\$4.10		13.35
Size 4, ~24 pcs/lb. 50041	\$5.10		13.35
Size 5, ~21 pcs/lb. 50051	—		13.35
Size 6, ~18 pcs/lb. 50061	—		13.35
Size 9, ~8 pcs/lb. 50091	—		13.35
Size 12, ~3 pcs/lb. 50121	—		13.35
2-Hole			lb
Size 6, ~17 pcs/lb. 50062			\$13.35
Size 7, ~12 pcs/lb. 50072			\$13.35
Size 8, ~9 pcs/lb. 50082			\$13.35
Size 9, ~7 pcs/lb. 50092			\$13.35
Size 10, ~5 pcs/lb. 50102			\$13.35
Assortment, 1 lb. 50012			\$17.35

Solid, 1-hole, and 2-hole stoppers in sizes 1-6.

Continued



Accurately determine the correct size of rubber stopper, size 00-12, or cork, size 0-26, for any vessel using this durable plastic tool. Reverse side has a circular scale for sizing 4-52 mm OD tubing.

Stopper/Tubing Sizer. 133920000 \$13.80

Non-fenestrated towels are sterile, latex-free and poly-lined.



Surgical Drape, 18 x 26", Bx/50. 44025-515..... \$25.50

Tape



Tape

Autoclave, steam


$\frac{3}{4}$ " x 60 Yd Roll. **44073ST048** \$8.50
 1" x 60 Yd Roll. **44073ST036** \$10.35

Kraft paper with rubber resin adhesive.



Label, Write-On™, Bel-Art
 1" Core, 15 yd (540") roll

$\frac{1}{2}$ " Wide, each \$5.65
 Green. **134623005** Yellow. **134622005**
 $\frac{3}{4}$ " Wide, Pk/4 \$24.35
 Green. **134820075** Yellow. **134810075**
 White. **134800075**
 3" Core, 40 yd roll
 $\frac{1}{2}$ " Wide, each \$10.30
 Green. **134633005** White. **134630005**
 $\frac{3}{4}$ " Wide, each \$13.70
 Green. **134633075** Yellow. **134632075**
 $\frac{3}{4}$ " Wide, White, Pk/4. **134830075**..... \$54.70
 $\frac{3}{4}$ " Wide, Mixed, Pk/6. **134630600** \$83.25
 blue, green, orange, red, white & yellow



Label Tape Dispenser, multi-roll. 134730000..... \$43.30



Versatile multi-format block holds (96) 0.2 mL tubes/strips, (39) 0.5 mL tubes, or (1) nonskirted 96-well PCR plate. In-situ adapter sold separately. Touchscreen control and height adjustable heated lid. 2-year warranty.



0.2 mL Tubes



96-Well Plate



0.5 mL Tubes



In-Situ Adapter

Thermal Cycler, multi-format block, Benchmark. **T5000-96** .. \$4,299.00
 In-Situ Adapter, holds 4 slides. **T5000-A-IS** \$422.00