



Accurately determine the correct size of rubber stopper, size 00-12, or cork, size 0-26, for any vessel using this durable plastic tool. Reverse side has a circular scale for sizing 4-52 mm OD tubing.

Stopper/Tubing Sizer. 133920000 \$13.80

Non-fenestrated towels are sterile, latex-free and poly-lined.



Surgical Drape, 18 x 26", Bx/50. 44025-515..... \$25.50

Tape



Tape

Autoclave, steam


¾" x 60 Yd Roll. **44073ST048** \$8.50
 1" x 60 Yd Roll. **44073ST036** \$10.35

Kraft paper with rubber resin adhesive.



Label, Write-On™, Bel-Art
 1" Core, 15 yd (540") roll

½" Wide, each \$5.65
 Green. **134623005** Yellow. **134622005**
 ¾" Wide, Pk/4 \$24.35
 Green. **134820075** Yellow. **134810075**
 White. **134800075**
 3" Core, 40 yd roll
 ½" Wide, each \$10.30
 Green. **134633005** White. **134630005**
 ¾" Wide, each \$13.70
 Green. **134633075** Yellow. **134632075**
 ¾" Wide, White, Pk/4. **134830075**..... \$54.70
 ¾" Wide, Mixed, Pk/6. **134630600** \$83.25
 blue, green, orange, red, white & yellow



Label Tape Dispenser, multi-roll. 134730000..... \$43.30



Versatile multi-format block holds (96) 0.2 mL tubes/strips, (39) 0.5 mL tubes, or (1) nonskirted 96-well PCR plate. In-situ adapter sold separately. Touchscreen control and height adjustable heated lid. 2-year warranty.



0.2 mL Tubes



96-Well Plate



0.5 mL Tubes



In-Situ Adapter

Thermal Cycler, multi-format block, Benchmark. T5000-96 .. \$4,299.00
 In-Situ Adapter, holds 4 slides. **T5000-A-IS** \$422.00



Fuses

¼ x 1¼"

Fast-Acting, AGC

0.1A/250V. 312-100	\$2.00
0.2A/250V. 312-200	\$2.00
1A/250V. 312-001	\$2.00
2A/250V. 312-002	\$2.00
3A/250V. 312-003	\$2.00
4A/250V. 312-004	\$2.00

Time-Delay (Slow-Blow), MDL

0.7A/250V. 313-700	\$2.00
1A/250V. 313-001	\$2.00
2.5A/250V	\$2.00

5 x 20 mm

Fast-Acting, GMA

0.5A/250V. 235-500	\$2.00
1A/250V. 235-001	\$2.00
3A/250V. 235-003	\$2.00
5A/125V. 235-005	\$2.00

Time-Delay (Slow-Blow), Glass tube

0.5A/250V. 239-500	\$2.00
1A/250V. 239-001	\$2.00
2A/250V. 239-002	\$2.00
3A/250V. 239-003	\$2.00
5A/250V. 239-005	\$2.00

Time-Delay (Slow-Blow), Ceramic

1A/250V. 215-001	\$2.00
-------------------------------	--------



Cone Point



Cup Point

Set Screws, steel

Cone Point, DIN 914

M4 x 0.7 x 5 mm length, 2 mm hex key size. 77702	\$0.50
M4 x 0.7 x 10 mm length, 2 mm hex key size. 77704	\$0.50
M4 x 0.7 x 12 mm length, 2.5 mm hex key size. 77706	\$0.50
M5 x 0.7 x 10 mm length, 2.5 mm hex key size. 77708	\$0.50

Cup Point, DIN 916

M3 x 0.5 x 3 mm length, 1.5 mm hex key size. 77722	\$0.50
M3 x 0.5 x 5 mm length, 1.5 mm hex key size. 77724	\$0.50

Tools




Cutter, diagonal, 5". PLR-766.00 \$10.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1½ x 3" opening. DEM-345.00 \$80.95

Gauges



Vernier Caliper measures in millimeters (0-150 mm with 0.01 mm precision) and inches (0-6"). Provided with a plastic case.

Vernier Caliper, digital, 3-button. 600VCD001 \$30.00
 Battery (replacement). **357** \$4.00



Compasses, replaceable tips, Germany

3" Long. DIV-320.00	\$11.95
6" Long. DIV-326.00	\$16.95
Tips (replacement), Pair. DIV-320.05	\$1.95



Reproducibility enhanced by rapid stop. Provided with padded, plastic case.

Micrometer, 0-25 mm in 0.01 mm div, Germany. 601A \$38.95



Measure inside and outside diameters of glass and plastic tubing, connecting pieces, rods, test tubes, etc. between 1.5 and 16 mm.

Tubing Sizer, metric diameter gauge, PP. 133900000 \$36.55

Continued