

Simport plastic culture tubes with dual-position caps.



**Tubes, Culture,**


- 10 x 75 mm (3 mL)  
Glass (Type 1, Class A borosilicate), rimless  
Cs/1000. **1503** ..... \$45.90
- 12 x 75 mm (5 mL)  
Polystyrene, Sterile, dual-position caps, Simport  
Cs/500, 20 sterile bags/25. **T415-2** ..... \$72.10  
Cs/1000, 8 sterile bags/125. **T415-3** ..... \$102.00
- Polystyrene, Nonsterile, Cs/1000. **110410** ..... \$36.40
- Polypropylene, Sterile, dual-position caps, Simport  
Cs/500, 20 sterile bags/25. **T415-2A** ..... \$72.10
- 13 x 100 mm (7 mL)  
Glass (Type 1, Class A borosilicate), rimless  
Cs/1000. **1510** ..... \$69.10
- 17 x 95 mm (14 mL)  
Polypropylene, Sterile, dual-position caps, Simport  
Cs/500, 20 sterile bags/25. **T416-2A** ..... \$95.95
- 16 x 125 mm (12 mL)  
Glass (Type 1, Class A borosilicate), rimless  
Cs/1000. **1515** ..... \$119.00
- 16 x 150 mm (15 mL)  
Glass (Type 1, Class A borosilicate), rimless  
Cs/1000. **1517** ..... \$113.00



RNase/DNase-free

**Tubes, Microcentrifuge, Economy, polypropylene**

- 0.5 mL, 0.1 mL graduations, natural, Pk/1000. **HS4422** ..... \$35.00
- 1.5 mL, 0.25 mL graduations, Heathrow  
Natural, max RCF is 16,400 x g  
Pk/500. **HS4323** ..... \$17.00  
Pk/1000. **HS4325** ..... \$33.00
- Black, max RCF is 18,000 x g, Pk/500. **HS4323K** ..... \$19.00
- 5 mL, natural, max RCF is 25,000 x g, Benchmark  
Nonsterile, 4 bags/50, Pk/200. **C1005T5** ..... \$53.00  
Sterile, 10 bags/20, Pk/200. **C1005T5ST** ..... \$55.75



Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable. Free of RNase/DNase.

**Tubes, Microcentrifuge, Premium, PP, G-Tube™, Bio Plas, USA**

- Amber color protects light sensitive specimens.
- 0.6 mL, Pk/500 ..... \$21.70  
Amber. **4048** Lavender. **4047** Natural. **4040**  
Mixed (blue, green, orange, pink, & yellow). **4049** ..... \$26.15
- 1.5 mL, 0.5 mL graduations, frosted area on side, Pk/500 ..... \$21.70  
Amber. **4038** Blue. **4033** Lavender. **4037** Natural. **4030**  
Mixed (blue, green, orange, pink, & yellow). **4039** ..... \$26.15
- 2.0 mL, 0.5 mL graduations, frosted area on side, Pk/500 ..... \$23.35  
Amber. **4058** Blue. **4053** Natural. **4050**  
Mixed (blue, green, orange, pink, & yellow). **4059** ..... \$27.80

**Tubes, Microcentrifuge, Siliconized, G-Tube™, Bio Plas, USA**

- 1.5 mL, natural color, Pk/500. **4165SL** ..... \$44.40  
*Siliconized products should NOT be autoclaved*



**Tubes, PCR, 0.2 mL thin-wall, polypropylene**

- Individual, flat caps, Simport  
Blue, Pk/1000. **T325-2VB** ..... \$45.75  
Natural, Pk/1000. **T325-2VN** ..... \$45.75
- Strips, Amplitube™, Low Profile, individually attached caps  
8-Tube, Pk/125 (1000 tubes). **T320-2LPN** ..... \$93.00  
Holds 0.1 mL. Tubes have flat caps & natural color.
- Strips, Bio Plas, *Caps are sold separately; available in natural only.*  
8-Tube, Pk/125 (1000 tubes) ..... \$72.10  
Blue. **5010-4** Natural. **5010-1**  
8-Tube Cap Strips, domed, Pk/300 (2400 caps). **5030-1**.. \$75.95  
12-Tube, Pk/100 (1200 tubes) ..... \$81.35  
Blue. **5020-4** Natural. **5020-1**  
12-Tube Cap Strips, domed, Pk/200 (2400 caps). **5040-1** \$75.95



**Vial Decapping Pliers, vinyl grip, USA**

- For 13 mm Crimp Seals. **D-13** ..... \$72.00
- For 20 mm Crimp Seals. **D-20** ..... \$72.00
- For 30 mm Crimp Seals. **D-30** ..... \$72.00