


Tube Storage



Baskets, test tube, aluminum, Heathrow

- Small, 15.2 x 15.2 x 15.2 cm(H). **HS20341A** \$27.00
 Medium, 25.6 x 15.2 x 15.2 cm(H). **HS20341B** \$30.00



Each cube holds (36) 1.5/2.0 mL tubes, or (96) 0.2 mL tubes, or (12) 0.2 mL strips, or (1) 96-well PCR plate and (10) 1.5/2.0 mL tubes. Maintains 0°C for 4 hours.

Tube Cooler, benchtop, 0°C, CoolCube™, Benchmark

- Microtube and PCR Plate. **R1000** \$112.00

SB2CC-64 Cryogenic
Tube Storage Box



Tube Storage Boxes, Cardboard, vented

- 2" Box Only, 5¼ x 5¼ x 2", each. **SB2CC** \$1.65
 2" box w/dividers, white

- ★ 100-Cell (10 x 10), Pk/12. **SB2CC-100** \$26.95
 81-Cell (9 x 9), Pk/12. **SB2CC-81** \$26.95
 64-Cell (8 x 8), Pk/12. **SB2CC-64** \$26.95

- 2" Mini box w/dividers, 3 x 3 x 2"
 25-Cell (5 x 5), Pk/12. **MB2C-25** \$29.55
 Holds 13 mm vials and 4.0-5.0 mL cryovials.

- 3" Box Only, 5¼ x 5¼ x 3", each. **SB3C** \$1.75
 3" box w/dividers

- 100-Cell (10 x 10), Pk/12. **SB3C-100** \$28.00
 81-Cell (9 x 9), Pk/12. **SB3C-81** \$28.00
 64-Cell (8 x 8), Pk/12. **SB3C-64** \$28.00

Dividers only, Fit 2 & 3" Boxes

- 100-Cell (10 x 10), Pk/12. **CD-100** \$11.55
 Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.
 81-Cell (9 x 9), Pk/12. **CD-81** \$11.55
 Holds 13 mm vials and 4.0-5.0 mL cryovials.
 64-Cell (8 x 8), Pk/12. **CD-64** \$11.55
 Holds 14 mm vials.
 49-Cell (7 x 7), Pk/12. **CD-49** \$11.55
 Holds 14 mm vials.
 25-Cell (5 x 5), 25 mm cell, Pk/12. **CD-25** \$11.55
 16-Cell (4 x 4), 31 mm cell, Pk/12. **CD-16** \$11.55

- 5" Box Only, 5¼ x 5¼ x 5", each. **SB4C** \$3.20
 5" Box w/dividers
 64-Cell (8 x 8), Pk/12. **SB4C-64** \$46.35
 36-Cell (6 x 6), 15 mL tubes, Pk/12. **SB4C-36** \$46.35
 16-Cell (4 x 4), 50 mL tubes, Pk/12. **SB4C-16** \$46.35
 Dividers only, Fit 5" Box
 64-Cell (8 x 8), Pk/12. **CD4-64** \$11.55
 36-Cell (6 x 6), Pk/12. **CD4-36** \$11.55
 16-Cell (4 x 4), Pk/12. **CD4-16** \$11.55



Holds 1.5-2.0 mL tubes
in a two-piece rack.


Tube Storage Boxes, Polystyrene Foam (styrofoam), USA

- 50-Well, Pk/10. **HS4283P** \$42.50
 100-Well, Pk/10. **HS4282** \$53.20



Tube Racks, Grip, polypropylene, No-Wire™, Bel-Art

- 15-16 mm Tubes, 60-place, each. **187490001** \$26.80
 18-20 mm Tubes, 40-place, each. **187490002** \$26.80



Tube Racks, Interlocking, polypropylene, Heathrow

- Cube, mixed colors, Pk/5. **HS29050A** \$42.35
 Holds (4) 50 mL tubes, (10) 15 mL tubes,
 (12) 12 x 75 or 12 x 100 tubes, or (16) 1.5-2.0 mL tubes.
 4-Way, mixed colors, Pk/5. **HS29022G** \$39.05
 Holds (4) 50 mL tubes, (12) 15 mL tubes,
 (32) 1.5 mL tubes, or (32) 0.5 mL tubes.
 4-Way Rotatable, Rota-Rack™, four connected modules
 Large Tubes, 22.4 x 11 x 11 cm, each. **HS29040A** \$41.20
 Each module holds (2) 50 mL tubes, (8) 15 mL tubes,
 (6) 20 mm tubes, or (10) 10-13 mm tubes. Mixed colors.

Continued

Continued