

**Tube Racks, Mega Rack**, polypropylene, Heathrow  
10-13 mm Tubes

- 216-Place, 33.7 x 26.7 x 9.1 cm(H)  
Blue, each. **HS24322A** ..... \$38.00
- White, each. **HS24322E** ..... \$38.00
- 432-Place, 67.4 x 26.7 x 9.1 cm(H)  
Blue, each. **HS120113** ..... \$69.00
- White, each. **HS120112** ..... \$69.00



80-Well



96-Well Reversible

**Tube Racks, Microtube**, polypropylene, Heathrow


- 80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm  
Holds 1.5-2.0 mL tubes in conical wells.
- Blue, Pk/5. **HS29025B** ..... \$25.15
- Pink, Pk/5. **HS29025D** ..... \$25.15
- Mixed, Pk/5. **HS29025G** ..... \$25.15
- 96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm  
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.
- Blue, Pk/5. **HS2345C** ..... \$48.90
- Green, Pk/5. **HS2345D** ..... \$48.90
- Pink, Pk/5. **HS2345E** ..... \$48.90

Available in blue,  
green, orange,  
purple & yellow.




**Tube Rack, PCR**, polypropylene, mixed colors, Heathrow

- Reversible, 96-Well, Pk/5. **HS2344A** ..... \$52.90
- Holds 0.2, 0.5 and 1.5 mL tubes.



**Tube Racks, Peg**, 14-17 mm Tubes, 50-place, polypropylene

- Blue, Pk/2. **HS24311B** ..... \$22.60
- White, Pk/2. **HS24311A** ..... \$22.60




Tilt your test tube rack 25°, 35°, or 45°. Create agar slant tubes, or adjust the angle of tube racks when pipetting to reduce hand strain.

- Tube Rack Slanter**, Bel-Art. **188900100** ..... \$21.50
- Platform size: 19.6 cm(W) x 12.7 cm(D)



Half-Size



Full-Size


**Tube Racks, Teardrop, OneRack™**, POM, Heathrow

- Half-Size  
16 mm Tubes, 36-place, 12.7 x 12.7 x 7 cm(H)  
Magenta, Pk/8. **HS27502D** ..... \$99.40
- Full-Size  
13 mm Tubes, 72-place, 20.2 x 10.4 x 5.9 cm(H)  
Blue, Pk/8. **HS27511B** ..... \$129.00
- Magenta, Pk/8. **HS27511D** ..... \$129.00
- 16 mm Tubes, 72-place, 25 x 12.7 x 7 cm(H)  
Blue, Pk/8. **HS27512B** ..... \$129.00
- Magenta, Pk/8. **HS27512D** ..... \$129.00



Polypropylene cube and lid has two snap-in, gridded inserts to hold either (16) 15 mL tubes or (9) 50 mL tubes.

- Tube Racks, Transformer Cube**, Pk/5. **HS23271A** ..... \$53.45



120-Place  
Wire Rack

**Tube Racks, Wire**, epoxy-coated, blue, Poxgrid®, Bel-Art

- 10-13 mm Tubes, 120-place (6 x 20), each. **187860750** ..... \$45.80
- Dimensions: 33.1 x 11 x 6 cm(H)
- 15-16 mm Tubes  
100-Place (5 x 20), each. **187860780** ..... \$46.55
- Dimensions: 38.6 x 10.8 x 6.4 cm(H)
- 180-Place (9 x 20), each. **187860770** ..... \$46.80
- Dimensions: 38.6 x 18.3 x 6.4 cm(H)

Continued