

Spade Tips



- Angled (Kuehne), 11 cm (4¼")
- Stainless Steel, economy. **33-1041** \$4.15
- Teflon® Coated, Bel-Art. **379420000** \$22.50
- Straight, stainless steel, 10.5 cm, economy. **57.926** \$8.45

Square Point Tips, straight



- Stainless Steel, 12 cm (4¾"), economy. **1244** \$2.25
- Teflon® Coated, 11 cm (4¼"), Bel-Art. **379430000** \$24.95

Locking



Forceps, stainless steel

- Self-Closing (cross-locking), squeeze handles to open tips
- Fine Point, smooth tips
- Economy, 13 cm (5⅞"). **300-066** \$4.95
- Peer-Vigor, 11 cm (4¼"), Swiss. **57.324** \$42.40
- Medium Point, serrated tips, 17 cm (6⅝"). **TWZ-948.00** \$3.95



Slide-Locking, grid-serrated tips


- Fine Point
- Economy, 14 cm (5½"). **57.0357A** \$9.00
- Dumont, 16 cm (6¼"). **TWZ-101.72** \$45.75
- Medium Point
- Economy, 16 cm (6¼"). **203** \$9.00
- Dumont, 16 cm (6¼"). **TWZ-101.76** \$43.60
- Snap-Closing, angled, fine points, 15.5 cm (6"). **3610** \$5.80
- serrated tips & handles, guide pin

Tissue



Forceps, stainless steel

- Adson-Brown, 9 x 9 side grasping teeth, 12 cm. **17-574** \$17.40
- Adson Tissue, 1 x 2 teeth, 12 cm (4¾"). **17-560** \$4.90
- Iris Tissue, 1 x 2 teeth
- Straight, economy, 10 cm (4"). **26-295** \$7.70
- Angled, economy, 11.5 cm (4½"). **19-350ANG** \$5.00
- Half Curve, standard, 10 cm (4"). **45-822** \$17.60





- Russian Tissue, circular serrated tips, 15.5 cm (6"). **2275** \$5.25
- Tissue, 1 x 2 teeth, 14 cm (5½"). **17-145** \$4.35

Tongs



- Tongs, Crucible, serrated jaws, stainless steel**
- 24.5 cm (9¾"). **320-101** \$11.60
- 51 cm (20"). **320-120** \$28.10

Hemostats & Clamps




Hemostats, serrated jaws

- Hartman Mosquito, stainless steel
- Straight
- Economy, 9 cm (3½"). **12-180S** \$6.05
- Standard, 9 cm (3½"). **20-352 (525-025)** \$21.30
- Curved
- Economy, 9 cm (3½"). **14-030** \$6.60
- Standard, 10 cm (4"). **20-360** \$21.70
- Halstead Mosquito, 12.5 cm (5"), stainless steel
- Straight
- Economy. **14-065** \$7.75
- Standard. **20-501** \$22.45
- Curved
- Economy. **14-070** \$7.75
- Standard. **20-519** \$22.85
- Kelly, stainless steel
- Straight
- Economy, 14.5 cm (5¾"). **14-225** \$7.75
- Standard, 14.5 cm (5¾"). **525-225** \$22.40
- Curved
- Economy, 14 cm (5½"). **14-230** \$7.75
- Standard, 14.5 cm (5¾"). **20-568** \$22.40
- Rochester-Pean, 20.5 cm (8"), stainless steel
- Straight
- Economy. **1920** \$6.50
- Standard. **525-435** \$48.50
- Curved
- Economy. **1930** \$6.50
- Standard. **525-440** \$48.50
- Polypropylene, 12.5 cm (5"). **HEM-300.00** \$1.50

Continued



Specimen Forceps, 25 cm (9 3/4"), serrated jaws. **300-500** \$17.25



Sponge Forceps, Foerster, serrated jaws, stainless steel
 Straight
 17.5 cm (7"), economy. **14-660** \$16.30
 24 cm (9 1/2"), economy. **14-670** \$19.65
 Curved, 17.5 cm (7"), economy. **14-665** \$16.30

Straight  300-120


Curved  300-121

Clamps

Dieffenbach Serrefine, stainless steel, economy
 Straight, serrated jaws and handles
 48 mm (1 7/8"). **12-247** \$8.90
 60 mm (2 3/8"). **300-120** \$6.80
 Curved, serrated jaws and handles
 48 mm (1 7/8"). **12-248** \$9.15
 60 mm (2 3/8"). **300-121** \$6.80




Schwartz Vessel Clips, 1.7 x 8 mm jaws, 2.5 cm long, Germany
 Straight Jaws, Roboz. **RS-5452** \$21.75
 Slightly Angled Jaws, Roboz. **RS-5457** \$21.75



8.6 cm shaft; 5 mm serrated jaw
 Micro Alligator Forceps, stainless steel. **2769** \$50.00


Needle Holders & Wound Closure



Needle Holders, stainless steel

Baumgartner, cross-serrated jaws, 14 cm (5 1/2"). **16-285S** \$11.05
 Castroviejo, cross-serrated jaws w/lock, 14 cm. **15-302** \$44.45
 Crile-Wood, cross-serrated jaws, 15 cm (6")
 Economy, central groove. **16-293S** \$13.45
 Standard. **32-485** \$33.75
 German. **320-366** \$60.80
 Derf, 12 cm (4 3/4"), economy. **16-238** \$9.70
 serrated jaws with central groove
 Olsen-Hegar, cross-serrated jaws with scissors, 14 cm (5 1/2")
 Economy, central groove. **16-296** \$14.30
 Economy, No Groove. **16-296A** \$14.45
 Standard, central groove. **32-500** (555-160) \$35.25
 German, central groove. **320-556** \$71.50
 Webster, smooth jaws, 12 cm (4 3/4")
 Economy. **15-035** \$14.30
 Standard. **32-101** (555-035) \$30.75
 Tungsten Carbide. **0817** \$45.00
 German. **320-101** \$60.50

Michel Wound Closure System



Michel Clip Applier/Remover, double-ended. **3700** \$13.75
 11.5 cm (4 1/2"), stainless steel



Michel Wound Clips, stainless steel
 11 mm, Bx/1000. **380-519** \$74.00
 14 mm, Bx/1000. **380-543** \$89.00
 16 mm, Bx/1000. **380-568** \$97.00
 18 mm, Bx/1000. **380-584** \$98.00
 22 mm, Bx/1000. **380-626** \$116.00

Reflex Wound Closure System

CellPoint Scientific, USA

Reflex Clip Applier


For 7 mm Clips. **204-1000** \$155.00
 ★ For 9 mm Clips. **202-1000** \$155.00

Reflex Wound Clips, nonsterile, stainless steel

7 mm, 5 packs of 20 clips, Bx/100. **203-1000** \$50.00
 ★ 9 mm, 5 packs of 20 clips, Bx/100. **201-1000** \$40.00

Reflex Clip Remover, fits 7 mm & 9 mm clips. 205-1000 \$99.00

Accessories



Bibulous Paper, 10 x 15 cm, Bk/50 sheets. **BIB-46-50** $\frac{\text{bk}}{\$4.10}$ $\frac{\text{cs/48 bk}}{190.00}$


Forceps

Spade Tips




Angled (Kuehne), 11 cm (4 1/4")
 Stainless Steel, economy. **33-1041** \$4.15
 Teflon® Coated, Bel-Art. **379420000** \$22.50
 Straight, stainless steel, 10.5 cm, economy. **57.926** \$5.85

Square Point Tips, straight



Stainless Steel, 12 cm (4 3/4"), economy. **1244** \$2.25
 Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379430000** \$24.95



Minimize contamination during slide & coverslip handling. Will not scratch; corrosion and heat resistant.

Hemostat, Polypropylene, 12.5 cm. **HEM-300.00** \$1.50

Cytology Funnels

Single,
Metal Clips
Required



Double,
All Plastic



White filter cards will accept up to 0.5 mL of sample, while tan filters will accept volumes up to 0.4 mL.

Cytology Funnels, Simport CytoSep™

Shandon (Thermo) Cytospin® CytoCentrifuges

Single. *Metal slide clips required*

White Filter & Cap, Cs/500. **M964-10FW** \$402.00
 (59 91 040, 110 25 48)

Tan Filter & Cap, Cs/500. **M964-10FT** \$402.00

Single. *All plastic - metal slide clips not required*

White Filter & Cap, Cs/480. **M965-10FW** \$558.00
 (A787 100003)

Tan Filter & Cap, Cs/480. **M965-10FT** \$558.00

Double. *All plastic - metal slide clips not required*

White Filter & Cap, Cs/480. **M965-20FW** \$714.00
 (A787 100005)

Immersion Oils & Bottles



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities. **USA**

Immersion Oil

Type A, low viscosity (150 cSt). *Short W.D. optics*

7.5 mL in LDPE dispensing bottle. **18101** \$7.75
 20 mL in HDPE dispensing bottle. **18102** \$10.85
 480 mL in amber bottle. **18104** \$58.20

Type B, high viscosity (1250 cSt). *Long W.D. optics*

7.5 mL in LDPE dispensing bottle. **18105** \$7.75
 20 mL in HDPE dispensing bottle. **18106** \$10.85
 480 mL in amber bottle. **18108** \$58.20

Type NVH, very high viscosity (21,000 cSt)

480 mL in amber bottle. **18113** \$67.00


Type HF, low fluorescence (700 cSt)

480 mL in amber bottle. **18118** \$119.00

Type FF, fluorescence free (170 cSt)

480 mL in amber bottle. **18124** \$95.60

A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



Bottle, Dispensing, HDPE, 20 mL, each. **18145** \$7.00



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.

Bottle, Balsam, glass, 45 mL (empty), **USA. 616** \$18.65