

Inverted Microscopes



CP-2A1

Jenco

CP-2A1 Inverted, Trinocular, Brightfield/Phase..... \$3,060.00

- Widefield 10X/20 eyepieces
- Trinocular, 30° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
- Six DIN, LWD plan achromat objectives: brightfield 10X, 25X, and 40X; plan phase 10X, 25X, and 40X
- 5-place, reverse-facing nosepiece
- Attachable mechanical stage includes three inserts
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm and filter holder; working distance of 55 mm for phase contrast, or 125 mm if swung out for brightfield
- Centerable, dedicated phase annuli mounted on sliders
- 30 watt halogen Koehler-type illuminator
- 48.25 cm (19")D x 21.6 cm (8½")W x 47 cm (18½")H
- 5-year mechanical & optical/2-year electrical warranty



BC-366

Jenco

BC-366 Inverted, Trinocular, BF/Phase/Fluor \$9,399.00

- Widefield 10X/22 eyepieces
- Binocular, 45° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
- Side camera port
- Six infinity, LWD plan achromat objectives: brightfield/fluorescence 10X, 25X, and 40X; phase 10X, 25X, and 40X
- 6-place, side-facing nosepiece
- Attachable mechanical stage
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm and filter holder; centerable, dedicated phase annuli mounted in turret
- 30 watt halogen transmitted light; Koehler-type illuminator
- 100 watt mercury epi-fluorescence
- Blue (FITC) and green (TRITC) fluorescence filter sets
- 5-year mechanical & optical/2-year electrical warranty

Meiji Access www.meijitechno.com
and call us to place an order

Polarizing Compound Microscopes



Meiji Access www.meijitechno.com
and call us to place an order

Metallurgical Compound Microscopes

Inverted Microscopes



IM7100



IM7200

Meiji

IM7100 Inverted, Binocular, BF/Polarized Light \$6,322.00

- Widefield 10X/22 eyepieces
- Binocular, 30° inclined, Siedentopf design viewing head with depolarizer has single diopter; 53-75 mm interpupillary adjustment
- Infinity-corrected, plan achromat E5X, E10X, E20X, and E50X objectives on 5-place, side-facing nosepiece
- Fixed stage measures 18 cm(W) x 24.5 cm(D)
- Coaxial coarse and fine focus controls; coarse has tension adjustment
- 30 watt halogen incident (reflected) light Koehler illuminator
- Filters provided: polarizer and rotatable analyzer, blue, green, and neutral density
- Camera port is built into front of cast aluminum alloy body
- 55.9 cm (22") x 25.4 cm (10") x 41.3 cm (16.25")H
- Lifetime mechanical & optical/1-year electrical warranty

IM7200 Inverted, Trinocular, BF/Polarized Light \$6,511.00

- Trinocular version of IM7100
- 55.9 cm (22") x 25.4 cm (10") x 47.5 cm (18.75")H

Optional components that can be purchased for the IM7100/IM7200 models include eyepieces, viewing heads, 100X objective, attachable mechanical stage, eyepiece and stage micrometers, and camera adapters.

Jenco Access www.jencointernational.com
and call us to place an order