



**Trocar**



Trocar, 4-in-1 set, 13 cm, stainless steel. **1596**..... \$17.50

**Accessories/Maintenance/Storage**



**Cutting Boards**, white high-density polyethylene  
 10 x 6 x 1/2" (25.4 x 15.2 x 1.3 cm). **1054319**..... \$6.95  
 18 x 12 x 1/2" (45.7 x 30.5 x 1.3 cm). **105CB1218**..... \$16.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1 1/2 x 3" opening. **DEM-345.00**..... \$80.95



Tools not included.

**Dissection Cases**  
 6-Loop with snap fastener, black leatherette. **1257370**..... \$1.30  
 12-Loop with snap fastener, black leatherette. **1257371**..... \$1.40  
 13-Loop with zipper, brown leatherette. **36900**..... \$7.80  
 Features red velveteen interior, protector flap, and ID pouch.

**Dissecting Kits**

*Inquire about customizing a kit for your surgical procedure.  
 Kit components are discounted 10%.*

**Basic kit. 38501**..... \$5.30  
 Contents: Dissecting Case, 6-loop (1257370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

**Mouse kit. 38551**..... \$17.85  
 Contents: Dissecting Case, 6-loop (1257370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

**Dissecting Pans**, stainless steel (without wax)  
 20 7/8 x 12 13/16 x 2 1/2". **AM20029**..... \$21.90  
 25 3/4 x 17 3/4 x 2 1/4". **AM68357**..... \$64.90



**Dissecting Pins**, T-pin, 2", nickel-plated steel  
 1/2 lb Box, ~430 pcs. **214500**..... \$14.85




**Insect Pins**, 38 mm (1 1/2") long, nylon head, Czech Republic  
 Black Enameled, Pk/100..... \$7.00  
 #000, 0.25 mm dia. **355425** #1, 0.40 mm dia. **355440**  
 #00, 0.30 mm dia. **355430** #2, 0.45 mm dia. **355445**  
 #0, 0.35 mm dia. **355435** #3, 0.50 mm dia. **355450**  
 Stainless Steel, Pk/100..... \$8.75  
 #000, 0.25 mm dia. **355525** #1, 0.40 mm dia. **355540**  
 #00, 0.30 mm dia. **355530** #2, 0.45 mm dia. **355545**  
 #0, 0.35 mm dia. **355535** #3, 0.50 mm dia. **355550**



Instrument Cup, silicone, 3 1/2" opening, 2 1/2" deep. **150-1005**..... \$43.00



**Instrument Storage Boxes**, 10 x 4 x 2", 18-8 stainless steel  
 Small, 7 x 3 x 1 1/2". **35-406**..... \$26.60  
 Large, 10 x 4 x 2". **35-408**..... \$35.00



**Rulers**, 15 cm/6" long, scaled in both metric and English units  
 Vinyl, flexible  
 Transparent. **4461**..... \$0.95  
 White. **4460**..... \$0.55  
 Stainless Steel, 12 mm wide, **USA. 605-065**..... \$5.50

*Continued*



**Fuses**

¼ x 1¼"

Fast-Acting, AGC

0.1A/250V. <b>312-100</b> .....	\$2.00
0.2A/250V. <b>312-200</b> .....	\$2.00
1A/250V. <b>312-001</b> .....	\$2.00
2A/250V. <b>312-002</b> .....	\$2.00
3A/250V. <b>312-003</b> .....	\$2.00
4A/250V. <b>312-004</b> .....	\$2.00

Time-Delay (Slow-Blow), MDL

0.7A/250V. <b>313-700</b> .....	\$2.00
1A/250V. <b>313-001</b> .....	\$2.00
2.5A/250V .....	\$2.00

5 x 20 mm

Fast-Acting, GMA


0.5A/250V. <b>235-500</b> .....	\$2.00
1A/250V. <b>235-001</b> .....	\$2.00
3A/250V. <b>235-003</b> .....	\$2.00
5A/125V. <b>235-005</b> .....	\$2.00

Time-Delay (Slow-Blow), Glass tube


0.5A/250V. <b>239-500</b> .....	\$2.00
1A/250V. <b>239-001</b> .....	\$2.00
2A/250V. <b>239-002</b> .....	\$2.00
3A/250V. <b>239-003</b> .....	\$2.00
5A/250V. <b>239-005</b> .....	\$2.00

Time-Delay (Slow-Blow), Ceramic

1A/250V. <b>215-001</b> .....	\$2.00
-------------------------------	--------



Cone Point



Cup Point

**Set Screws, steel**

Cone Point, DIN 914

M4 x 0.7 x 5 mm length, 2 mm hex key size. <b>77702</b> .....	\$0.50
M4 x 0.7 x 10 mm length, 2 mm hex key size. <b>77704</b> .....	\$0.50
M4 x 0.7 x 12 mm length, 2.5 mm hex key size. <b>77706</b> .....	\$0.50
M5 x 0.7 x 10 mm length, 2.5 mm hex key size. <b>77708</b> .....	\$0.50

Cup Point, DIN 916

M3 x 0.5 x 3 mm length, 1.5 mm hex key size. <b>77722</b> .....	\$0.50
M3 x 0.5 x 5 mm length, 1.5 mm hex key size. <b>77724</b> .....	\$0.50

**Tools**




**Cutter, diagonal, 5". PLR-766.00** ..... \$10.95



Demagnetize or magnetize parts and tools.

**Demagnetizer/Magnetizer, 1½ x 3" opening. DEM-345.00** ..... \$80.95


**Gauges**



Vernier Caliper measures in millimeters (0-150 mm with 0.01 mm precision) and inches (0-6"). Provided with a plastic case.

**Vernier Caliper, digital, 3-button. 600VCD001** ..... \$30.40

Battery (replacement). **357** ..... \$4.00



**Dividers, replaceable tips, Germany**

3" Long. **DIV-320.00** ..... \$11.95

6" Long. **DIV-326.00** ..... \$16.95

Tips (replacement), Pair. **DIV-320.05** ..... \$1.95



Reproducibility enhanced by rapid stop. Provided with padded, plastic case.

**Micrometer, 0-25 mm in 0.01 mm div, Germany. 601A** ..... \$38.95



Measure inside and outside diameters of glass and plastic tubing, connecting pieces, rods, test tubes, etc. between 1.5 and 16 mm.

**Tubing Sizer, metric diameter gauge, PP. 133900000** ..... \$37.65

*Continued*