

Slide & Cassette Pens, Markers & Scribes


Write a 0.4 mm stroke on glass slides with this permanent marker.
Pilot pens are ideal for cytologists and pathologists.



Pens, Pilot® , permanent, extra fine point	ea	bx/12
Black. 44102	\$1.70	17.45
Blue. 44103	\$1.70	17.45
Green. 44105	\$1.70	17.45
Red. 44104	—	17.45



Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



Wax Pencil

Markers, USA

Wax Pencil (China Marker), Bx/12. **02089** \$11.50



The 0.75 mm tip is ideal for use on small surfaces, e.g., histology cassettes, microscope slides and cryogenic vials. The permanent ink is waterproof and resistant to most laboratory chemicals and solvents.

Diamond PerfectMark™, permanent, black, Bx/10. **1000** \$21.60 ↓

Fine Tip



Ultra-Fine Tip



Sharpie Retractable, permanent

Fine Tip	ea	bx/12
Black. 36701	—	29.65
Blue. 36703	\$2.85	29.65
Red. 36702	\$2.85	29.65
Ultra-Fine Tip	ea	bx/12
Black. 1735790	\$2.85	29.65
Blue. 1735792	—	29.65

Diamond-Tipped



Carbide-Tipped

Scribes, retractable

Carbide-Tipped (sharpenable). RSC	\$15.30
Carbide Tip (replacement). RCT-1	\$6.20
Diamond-Tipped. RSD	\$22.75
Diamond Tip (replacement). RDT-1	\$12.90




XH-2002

Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.

Slide Warmers

Small, 26 x 18 cm (10 x 7"). XH-2002	\$257.00
Large, 63 x 20 cm (25 x 8"). XH-2001	\$345.00

Staining



Ten slide compartments, separated by 1/2" of space, are additionally isolated by eight barrier dividers when the lid is closed. Individual slides are suspended 1/2" above the chamber floor by four pedestal posts and four corner posts.

Staining/Moisture Chambers

10-Slide, polystyrene, USA	ea	pk/3
Clear. 240-9000-010	\$54.95	154.00
Black. 240-9020-Z10	\$56.10	157.00



M920-1

20-Slide, Simport

Clear Lid (PETG), black base (ABS), each. M920-1	\$169.00
Black Lid (ABS), black base (ABS), each. M920-2	\$169.00



Continued

— Hot Plate, Magnetic Stirrer & Hot Plate Stirrers —




Hot Plates, analog
 6" Diameter Cast Iron Burner, 1000W. **310HOTPLT** \$27.00
 7.5 x 7.5" Ceramic Platform, 80-380°C. **H4000-H**..... \$239.00




Magnetic Stirrers, analog, Benchmark
 3.4 x 3.4" Platform, 100-1500 rpm, MiniMag™. **S1005**..... \$112.00
 (1-year warranty)
 7.1 x 7.1" Ceramic Platform, 60-1500 rpm. **H4000-S**..... \$239.00
 (2-year warranty)

Hot Plate Stirrer, 80-380°C, Benchmark. H4000-HS..... \$358.00



Hydrometers *Inquire*

— Ice Buckets & Pans —



Temperature Range: -196 to 100°C

Ice Buckets with lids, expanded PVC foam, Bel-Art, USA
 2.5 L Capacity (Handy)
 Black. **188482002**..... \$54.55
 Blue. **188482001** \$54.55
 Purple. **188482006** \$54.55
 4 L Capacity (N'Icer)
 Black. **188484002**..... \$66.50
 Blue. **188484001** \$66.50

Ice Pans w/o lids, expanded PVC foam, Bel-Art, USA
 1 L Capacity (Mini), 15.5 x 15.5 x 8.5 cm(H)
 Black. **188481102**..... \$40.65
 Blue. **188481101** \$40.65
 Purple. **188481106** \$40.65
 4 L Capacity (Midi), 33 x 24.5 x 12.5 cm(H)
 Black. **188484102**..... \$64.80
 Blue. **188484101** \$64.80
 Purple. **188484106** \$64.80
 9 L Capacity (Maxi), 44 x 35 x 13 cm(H)
 Black. **188489102**..... \$76.20
 Blue. **188489101** \$76.20
 Purple. **188489106** \$76.20



Incubators, MyTemp™ Mini Digital, Benchmark

- 20 L capacity (0.8 cu ft)
- Digital temperature control with LED display
- Two models: Heat Only, or Heat/Cool
- Accepts bottles and flasks up to 2 L
- Two adjustable/removable shelves
- 9.3"(D) x 10.3"(W) x 12.8"(H) interior;
 14.5"(D) x 13.2"(W) x 18.7"(H) exterior
- 2-year warranty

Heat Only, 1° above ambient to 60°C. **H2200-H** \$474.00
 Heat/Cool, 15° below ambient to 60°C. **H2200-HC** \$713.00
 Heat/Cool/CO₂, 10° below ambient to 55°C. **H2300-HC2** \$3,408.00

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



Markers, USA
 Wax Pencil (China Marker), Bx/12. **02089** \$11.50



Sharpie Retractable, permanent

Fine Tip	ea	bx/12
Black. 36701	---	\$29.65
Blue. 36703	\$2.85	29.65
Red. 36702	\$2.85	29.65
Ultra-Fine Tip	ea	bx/12
Black. 1735790	\$2.85	29.65
Blue. 1735792	---	29.65

Continued