




Simport



Simport  
Macrosette

**Tissue Cassettes, Attached Lids, slotted pores** 1 pk pk/3

Histosette™ I, snap-apart lids, Pk/1500..... \$131.00 371.00  
30° front labeling area.

Blue. **M490-6** Pink. **M490-3**  
Green. **M490-4** White. **M490-2**

1 pk pk/3  
Histosette™ II, integral lids, Pk/1500..... \$175.00 504.00  
45° front labeling area.

Blue. **M492-6** Pink. **M492-3**  
Green. **M492-4** White. **M492-2**

1 pk pk/3  
Uniset™, hinged lids, Pk/1500..... \$141.00 398.00  
35° front labeling area. (Alternative to Sakura's Uni-Cassette®)

Blue. **M505-6** Pink. **M505-3**  
Green. **M505-4** White. **M505-2**

Macrosette™, 2X normal depth, white, Pk/750. **M512**..... \$185.00  
Dimensions: 40.1 x 28.5 x 13 mm(H)

M480-4



Compatible with all  
available metal lids.

**Tissue Cassettes, without Lids, acetyl polymer, premium**

Capsettes™, **USA**, Pk/500..... \$42.40

Aqua. **6039** Lime. **6042**  
Blue. **6034** Peach. **6044**  
Burnt Orange. **6036** Pink. **6032**  
Chartreuse. **6041** Red. **6040**  
Gold. **6038** Rose. **6043**  
Green. **6033** Tan. **6035**  
Lavender. **6045** White. **6030**  
Lilac. **6037** Yellow. **6031**

Simport, 45° front labeling area  
Green, Pk/1500. **M480-4**..... \$97.80

**Cassette Storage**

**Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. 8342**..... \$121.00  
Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks  
in 4 drawers, each divided into 12 rows.  
Cabinets can be stacked up to 8 high.  
Dimensions: 49.5 x 38 x 21 cm.




**Tissue Cassette Cabinet, Steel, 4-drawer. 485400**..... \$250.00



Simport

Store 1500 cassettes or 990 embedding  
rings in 6 stackable drawers.  
High-impact polystyrene resists  
temperatures from -80°C to +80°C.  
Drawer dimensions: 40.5 x 23 x 5.1 cm.

**Tissue Cassette Drawers, Polystyrene, Pk/6. M495-6**..... \$159.00



LB-38



LB-30



LB-31



**Brushes**

Flat Bristles, paint style, camel hair

1/2" Wide, 6 1/4" total length. **LB-36**..... \$3.90

1" Wide, 7" total length. **LB-37**..... \$6.15

2" Wide, 7 3/4" total length. **LB-38**..... \$5.30

Round Bristles, camel hair

3/16" Diameter (artist style), 7 1/2" total length. **LB-31**..... \$1.75

3/4" Diameter, 8 1/2" total length. **LB-30**..... \$10.25

Ultra pure and nonflammable;  
chills to -62°F.

**Freeze Spray, 10 oz Can. 49185**..... \$12.00

