




Simport



Simport
Macrosette

Tissue Cassettes, Attached Lids, slotted pores 1 pk pk/3

Histosette™ I, snap-apart lids, Pk/1500..... \$131.00 371.00
30° front labeling area.

Blue. **M490-6** Pink. **M490-3**
Green. **M490-4** White. **M490-2**

Histosette™ II, integral lids, Pk/1500..... 1 pk pk/3
\$175.00 504.00
45° front labeling area.

Blue. **M492-6** Pink. **M492-3**
Green. **M492-4** White. **M492-2**

Unisette™, hinged lids, Pk/1500..... 1 pk pk/3
\$141.00 398.00
35° front labeling area. (Alternative to Sakura's Uni-Cassette®)

Blue. **M505-6** Pink. **M505-3**
Green. **M505-4** White. **M505-2**

Macrosette™, 2X normal depth, white, Pk/750. **M512**..... \$185.00
Dimensions: 40.1 x 28.5 x 13 mm(H)

M480-4



Compatible with all
available metal lids.

Tissue Cassettes, without Lids, acetyl polymer, premium

Capsettes™, **USA**, Pk/500..... \$42.40

Aqua. **6039** Lime. **6042**
Blue. **6034** Peach. **6044**
Burnt Orange. **6036** Pink. **6032**
Chartreuse. **6041** Red. **6040**
Gold. **6038** Rose. **6043**
Green. **6033** Tan. **6035**
Lavender. **6045** White. **6030**
Lilac. **6037** Yellow. **6031**

Simport, 45° front labeling area
Green, Pk/1500. **M480-4**..... \$97.80

Cassette Storage

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. 8342..... \$121.00
Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks
in 4 drawers, each divided into 12 rows.
Cabinets can be stacked up to 8 high.
Dimensions: 49.5 x 38 x 21 cm.




Tissue Cassette Cabinet, Steel, 4-drawer. 485400..... \$250.00



Simport

Store 1500 cassettes or 990 embedding
rings in 6 stackable drawers.
High-impact polystyrene resists
temperatures from -80°C to +80°C.
Drawer dimensions: 40.5 x 23 x 5.1 cm.

Tissue Cassette Drawers, Polystyrene, Pk/6. M495-6..... \$159.00



LB-38



LB-30



LB-31



Brushes

Flat Bristles, paint style, camel hair

1/2" Wide, 6 1/4" total length. **LB-36**..... \$3.90

1" Wide, 7" total length. **LB-37**..... \$6.15

2" Wide, 7 3/4" total length. **LB-38**..... \$5.30

Round Bristles, camel hair

3/16" Diameter (artist style), 7 1/2" total length. **LB-31**..... \$1.75

3/4" Diameter, 8 1/2" total length. **LB-30**..... \$10.25

Ultra pure and nonflammable;
chills to -62°F.

Freeze Spray, 10 oz Can. 49185..... \$12.00



Borosilicate glass bottle;
blue PP cap (GL45)



Bottles, Media, graduated, autoclavable

| | |
|---|----------|
| 100 mL, Cs/10. 600BM0100 | \$75.85 |
| 250 mL, Cs/10. 600BM0250 | \$85.70 |
| 500 mL, Cs/10. 600BM0500 | \$98.65 |
| 1000 mL, Cs/10. 600BM1000 | \$140.00 |
| Bottle Caps (repl.), GL45, blue, Pk/10. B3000-CAP | \$35.30 |
| Bottle Drip-Free Rings (repl.), GL45, Pk/10. B3000-RIN | \$19.30 |



Bottles, Wash, Safety Labeled, 500 mL (empty), Bel-Art
4-Color diamond, LDPE bottle/PP cap

| | |
|---|---------|
| Soap, Pk/4. 117160014 | \$34.25 |
| Universal (label your own bottle), Pk/4. 117160009 | \$33.50 |



Globally Harmonized System (GHS) wash bottles are vented to prevent solvent drips. Polyethylene bottles (LDPE) with polypropylene (PP) caps.

Bottles, Wash, GHS Labeled, 500 mL (empty), Bel-Art

| | |
|---|---------|
| Acetone, Pk/4. 124160001 | \$56.40 |
| Bleach, Pk/4. 124160015 | \$56.40 |
| Deionized Water, Pk/4. 124160003 | \$56.40 |
| Ethanol, Pk/4. 124160019 | \$56.40 |
| Methanol, Pk/4. 124160011 | \$56.40 |
| Saline, Pk/4. 124160013 | \$56.40 |
| Water, Pk/4. 124160017 | \$56.40 |



Bottles, Wash, Integral Spout, LDPE bottle/PP cap

| | |
|--|---------|
| 250 mL (8 oz), Cs/10. 601633 | \$46.60 |
| 500 mL (16 oz), Cs/10. 601634 | \$56.50 |



Bottles, Narrow-Mouth, HDPE bottle/polypropylene cap

| | |
|--|---------|
| 125 mL, Pk/12. 106200015 | \$21.50 |
| 250 mL, Pk/12. 106200016 | \$29.70 |
| 500 mL, 28 mm cap, Pk/12. 106200007 | \$35.90 |
| 1000 mL, Pk/6. 106200008 | \$28.30 |

Glass



HDPE



Bottles, Wide-Mouth

Glass (flint), straight-sided, PP cap w/polyethylene liner

| | |
|--|--------|
| 250 mL (8 oz), 70-400 cap, each. 17983P-1 | \$2.90 |
| 500 mL (16 oz), 89-400 cap, each. 17984P-1 | \$3.60 |
| 1000 mL (32 oz), 89-400 cap, each. 19316P-1 | \$4.80 |

HDPE bottle/polypropylene cap

| | |
|--|---------|
| 125 mL, 38 mm cap, Pk/12. 106250005 | \$21.50 |
| 250 mL, Pk/12. 106250006 | \$30.95 |
| 500 mL, 53 mm cap, Pk/12. 60033409 | \$35.90 |
| 1000 mL, Pk/6. 106250018 | \$33.35 |

Polypropylene bottle and cap, autoclavable

| | |
|--|---------|
| 125 mL, 38 mm cap, Pk/12. 60033307 | \$21.35 |
| 250 mL, 43 mm cap, Pk/12. 60033308 | \$30.80 |
| 500 mL, 53 mm cap, Pk/12. 106320007 | \$35.80 |
| 1000 mL, 63 mm cap, Pk/6. 60033310 | \$31.45 |



Bottle Top Dispensers

| | |
|---------------------------------|----------|
| 0.5-5 mL. 600BTDR2 | \$251.00 |
| 5-60 mL. 600BTDR5 | \$298.00 |

Boxes



Matchbox Style

Box, Matchbox Style (Pill Box), cardboard, white

| | |
|---|---------|
| 2 1/8 x 1 1/8 x 5/8", Pk/72. 20047 | \$32.95 |
|---|---------|

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. **8342**

Each cardboard box can store approximately 540 paraffin blocks.

Continued