




Simport



Simport  
Macrosette

**Tissue Cassettes, Attached Lids, slotted pores** 1 pk pk/3

Histosette™ I, snap-apart lids, Pk/1500..... \$131.00 371.00  
30° front labeling area.

Blue. **M490-6** Pink. **M490-3**  
Green. **M490-4** White. **M490-2**

Histosette™ II, integral lids, Pk/1500..... \$175.00 504.00  
45° front labeling area.

Blue. **M492-6** Pink. **M492-3**  
Green. **M492-4** White. **M492-2**

Uniset™, hinged lids, Pk/1500..... \$141.00 398.00  
35° front labeling area. (Alternative to Sakura's Uni-Cassette®)

Blue. **M505-6** Pink. **M505-3**  
Green. **M505-4** White. **M505-2**

Macrosette™, 2X normal depth, white, Pk/750. **M512**..... \$185.00  
Dimensions: 40.1 x 28.5 x 13 mm(H)

M480-4



Compatible with all  
available metal lids.

**Tissue Cassettes, without Lids, acetyl polymer, premium**

Capsettes™, **USA**, Pk/500..... \$42.40

Aqua. **6039** Lime. **6042**  
Blue. **6034** Peach. **6044**  
Burnt Orange. **6036** Pink. **6032**  
Chartreuse. **6041** Red. **6040**  
Gold. **6038** Rose. **6043**  
Green. **6033** Tan. **6035**  
Lavender. **6045** White. **6030**  
Lilac. **6037** Yellow. **6031**

Simport, 45° front labeling area  
Green, Pk/1500. **M480-4**..... \$97.80


**Cassette Storage**

**Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. 8342**..... \$121.00  
Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks  
in 4 drawers, each divided into 12 rows.  
Cabinets can be stacked up to 8 high.  
Dimensions: 49.5 x 38 x 21 cm.




**Tissue Cassette Cabinet, Steel, 4-drawer. 485400**..... \$250.00



Simport

Store 1500 cassettes or 990 embedding  
rings in 6 stackable drawers.  
High-impact polystyrene resists  
temperatures from -80°C to +80°C.  
Drawer dimensions: 40.5 x 23 x 5.1 cm.

**Tissue Cassette Drawers, Polystyrene, Pk/6. M495-6**..... \$159.00



LB-38



LB-30



LB-31



**Brushes**

Flat Bristles, paint style, camel hair

1/2" Wide, 6 1/4" total length. **LB-36**..... \$3.90

1" Wide, 7" total length. **LB-37**..... \$6.15

2" Wide, 7 3/4" total length. **LB-38**..... \$5.30

Round Bristles, camel hair

3/16" Diameter (artist style), 7 1/2" total length. **LB-31**..... \$1.75

3/4" Diameter, 8 1/2" total length. **LB-30**..... \$10.25

Ultra pure and nonflammable;  
chills to -62°F.

**Freeze Spray, 10 oz Can. 49185**..... \$12.00





**Boxes, Glass Disposal**, cardboard, polyethylene bag included  
 Bench Model, 8 x 8 x 10", Pk/6. **246530002** ..... \$86.10  
 Floor Model, 12 x 12 x 27", Pk/6. **246530001** ..... \$97.60




**Boxes, Plastic**, (see following chart)

Plastic Boxes				
Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
<b>H-1</b>	1 x 1 x 1/4"	-	Reticles, Pointers	\$0.40
<b>H-1/11B</b>	1 x 1 x 1/4"	1	Reticles, Pointers	\$0.50
<b>H-54</b>	2 9/16 x 2 9/16 x 1/2"	-	Filters	\$1.20
<b>H-54/2916G</b>	2 9/16 x 2 9/16 x 1/2"	1	Filters	\$1.35
<b>H-28/334B2</b>	3 3/4 x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40
<b>H-10</b>	2 7/8 x 1 3/16 x 3/8"	-	Small Parts	\$0.60
<b>H-21/512B</b>	5 1/2 x 1 3/4 x 1/2"	1	Forceps & Scissors	\$1.80



**Brushes**

Bench, horse hair, 10" long. **BRS-969.10** ..... \$9.20  
 Bristle head: 4 3/4"(L) x 1 1/2"(W); bristle length: 1".  
 Brass Wire, 9 3/4" long. **LB-40** ..... \$8.70  
 Bristle head: 5"(L) x 3/4"(W); bristle length: 1 1/8".  
 Flat Bristles, paint style  
 Camel Hair  
 1/2" Wide, 6 3/4" total length. **LB-36** ..... \$3.90  
 1" Wide, 6 1/2" total length. **LB-37** ..... \$6.15  
 2" Wide, 7 3/4" total length. **LB-38** ..... \$5.30  
 Horsehair, heavy-duty bristles, metal handle. **33761880** ..... \$0.85  
 Round Bristles, camel hair, 8 1/2" long. **LB-30** ..... \$10.25  
 Swivel-Head, black nylon bristles. **LB-14** ..... \$4.20  
 Bristle head: 4 1/4"(L) x 1 1/2"(W). Total length: 12".



Pen Style (scratch), retractable, 11.5 cm (4 1/2") long, Germany  
 Brass. **BRS-295.00** ..... \$12.20  
 Nylon. **BRS-293.00** ..... \$9.70  
 Steel. **BRS-292.00** ..... \$12.80

**Burners**



**Burner, Alcohol**, 120 mL, glass body, metal cap. **97-5313** ..... \$6.80  
 Wicks (replacement), 10-14 cm x 8 mm, Pk/12. **97-5310** ..... \$5.00



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing.

H-5550 gridded top, 3.8 cm (1 1/2") diameter, produces a uniform blue cone of light.



**Burners, Bunsen**, natural gas, Humboldt Mfg, **USA**  
 Standard, flame stabilizer. **H-5870** ..... \$29.15  
 High Temperature, Meeker-type. **H-5550** ..... \$98.80