

# GENERAL LAB SUPPLIES & EQUIPMENT



## Applicators, 6" long

Cotton Tipped, wooden shaft  
Nonsterile

Bg/100. <b>95-8702</b> .....	\$1.50
Bx/1000, 10 bg/bx. <b>95-8702</b> .....	\$12.00
Sterile Pk/1, Bx/100. <b>44024-106-1S</b> .....	\$7.80
Sterile Pk/2, Bx/200. <b>44024-106-2S</b> .....	\$8.25
Plain, wood, Bg/1000. <b>95-8701</b> .....	\$7.25



Rubberized Cloth

## Aprons

Disposable, polyethylene, 46"(L) x 28"(W)	
0.8 mil Thick (0.02 mm), Bx/100. <b>105DA2846</b> .....	\$11.85
1.75 mil Thick (0.044 mm), Bx/100. <b>59010485</b> .....	\$38.00
Rubberized Cloth, black	
36"(L) x 27"(W), each. <b>31012760</b> .....	\$7.85
42"(L) x 27"(W), each. <b>31012761</b> .....	\$8.70
★ 42"(L) x 36"(W), each. <b>31012763</b> .....	\$10.65
46"(L) x 36"(W), each. <b>31012762</b> .....	\$11.50



B4000-M

B4000-16



## Autoclaves, benchtop, Benchmark (1-year warranty)

8 L, Bioclave™ Mini. <b>B4000-M</b> .....	\$3,707.00
Chamber Dimensions: 6.7" ID x 12.5" (D)	
16 L, Bioclave™. <b>B4000-16</b> .....	\$5,095.00
Chamber Dimensions: 9" ID x 13.75" (D)	



Maximize DNA visualization with this light blue film that is coated on both sides (double emulsion). Ideal for all autoradiographic and electrophoretic applications. Compatible with automatic processors and manual development methods.


## Autoradiography (Western Blotting) Film

5 x 7", 100 sheets/pk, Cs/5 pk. <b>XC59X</b> .....	\$268.00
8 x 10", 100 sheets/pk, Cs/5 pk. <b>XC6A2</b> .....	\$348.00

## Bags



Collect biohazardous waste in this space-saving, autoclavable bag. Printed with a biohazard symbol and sterilization indicator patch. Measures 20.3 x 30.5 cm (8 x 12") and is 1.57 mil (0.04 mm) thick polypropylene.



HS10044H

<b>Bags, Biohazard, Benchtop, Bx/500. HS10320</b> .....	\$76.00
Biohazard Bag Holder, Benchtop. <b>HS10044H</b> .....	\$23.00
Epoxy-coated wire holder measures 9.1 x 12.8 x 21.4 cm.	



Bags display a biohazard symbol and precautionary procedures. Sterilization indicator patch darkens when autoclaved. Polypropylene bags are bright red. Durable, extra heavy-duty bags provide maximum protection against dangerous tears and leaks.

<b>Bags, Biohazard, Standard, 1.5 mil (0.038 mm) thick</b>	
14 x 19" (36 x 48 cm), Bx/200. <b>131641419</b> .....	\$69.00
19 x 23" (48 x 58 cm), Bx/200. <b>131641923</b> .....	\$99.00
25 x 35" (64 x 89 cm), Bx/200. <b>131642535</b> .....	\$243.00

<b>Bags, Biohazard, Heavy-Duty, 2 mil (0.05 mm) thick</b>	
14 x 19" (36 x 48 cm), Bx/200. <b>131651419</b> .....	\$203.00
19 x 23" (48 x 58 cm), Bx/200. <b>131651923</b> .....	\$296.00
25 x 35" (64 x 89 cm), Bx/200. <b>131652535</b> .....	\$508.00
31 x 38" (79 x 96 cm), Bx/200. <b>131653138</b> .....	\$638.00

Safely transport patient specimens with these OSHA and NCCLS compliant bags. Each resealable, 2 mil thick polyethylene bag has an attached document pouch and imprinted biohazard symbol.



<b>Bag, Specimen Transport, 6 x 9", Bx/1000. 4919</b> .....	\$64.15
---	---------



<b>Bags, Waste Disposal, HDPE, Cleanware™</b>	
8 x 10", Pk/50. <b>131741008</b> .....	\$21.50
12 x 16", Pk/50. <b>131741012</b> .....	\$36.10



Seal polyethylene or polypropylene bags, sheets, or flat tubing from 0.01 to 0.8 mm thick. Sliding blade cuts excess material 19 mm from the seal.

<b>Bag Sealer, electric, 12" (30.5 cm) Wide. HS28272</b> .....	\$192.00
Bag Sealer Replacement Part Kit. <b>HS28280</b> .....	\$31.00