

Balances

Portable



Model	DWP-1001
Capacity (g)	1000
Readability (g)	0.1
Modes	g, oz, ct, gn
Pan	14.6 cm (5.75") stainless steel
Power	4 AA/AC adapter
Price	\$72.00

Warranty: 1 year

Toploading

Benchmark
W3200-1200



Model	W3200-500	W3200-1200
Capacity (g)	500	1200
Readability (g)	0.001	0.01
Modes	g, oz, ct	g, oz, ct
Pan	11.5 cm (4.5")	16 cm (6.3")
Draft Shield	Yes, 3-door glass	No
Price	\$1,174.00	\$950.00

Warranty: 2 years

Analytical

Benchmark



Model	★ W3100-120	W3100A-120
Calibration	External (100 g weight provided)	Internal (Quick-Cal™ software)
Capacity (g)	120	120
Readability (g)	0.0001	0.0001
Modes	g, oz, lb, ct	g, oz, lb, ct
Pan	9 cm (3.5")	9 cm (3.5")
Price	\$1,746.00	\$1,853.00

Warranty: 2 years

Baths



BSH1002

Aluminum
Blocks

Dry Baths, digital, Benchmark (2-year warranty)

1-Block Capacity. BSH1001	\$427.00
2-Block Capacity. BSH1002	\$532.00
4-Block Capacity. BSH1004	\$869.00

Blocks sold separately – please specify configuration

Blocks, aluminum

Holds (48) 0.2 mL Tubes. BSW02	\$118.00
Holds (24) 0.5 mL Tubes. BSW05	\$118.00
Holds (24) 1.5 mL Tubes. BSW1500	\$118.00
Holds (24) 1.5/2.0 mL Tubes. BSW1520	\$118.00
Holds (20) 12/13 mL Tubes. BSW13	\$118.00
Holds (12) 15 mL Tubes. BSW15	\$118.00
Holds (5) 50 mL Tubes. BSW50	\$118.00



40 x 1.5/2.0 mL



15 x 50 mL



41 x 15 mL

Racks are only included with model B2000-4.

Water Baths, digital, MyBath™, Benchmark (2-year warranty)

4 L, 3 tube racks. B2000-4	\$832.00
8 L. B2000-8	\$916.00
12 L. B2000-12	\$1,081.00



Economy

Polypropylene
Beakers



Premiu

Beakers

Glass (borosilicate), Griffin low form, single scale, economy

10 mL, no graduations, each. 6102022910	\$1.85
50 mL, each. 7-440050	\$3.10
100 mL, each. 7-440100	\$3.40
150 mL	
Pk/12. 61020229150	\$46.30
Cs/48. 61020229150	\$156.00
250 mL, each. 61020229250	\$3.50
400 mL, each. 61020229400	\$4.40
600 mL, each. 61020229600	\$5.15
800 mL	
Pk/6. 600BG800	\$63.45
Cs/24. 600BG800	\$234.00
1000 mL, each. 600BG1000	\$8.75
1000 mL, Heavy-Duty, each. 600BG1003	\$14.15
2000 mL, each. 600BG2000	\$16.00
4000 mL, each. 600BG4000	\$44.30

Continued