

## Balances

### Portable



<b>Model</b>	<b>DWP-1001</b>
<b>Capacity (g)</b>	1000
<b>Readability (g)</b>	0.1
<b>Modes</b>	g, oz, ct, gn
<b>Pan</b>	14.6 cm (5.75") stainless steel
<b>Power</b>	4 AA/AC adapter
<b>Price</b>	\$72.00

Warranty: 1 year

### Toploading

Benchmark  
W3200-1200



Model	W3200-500	W3200-1200
<b>Capacity (g)</b>	500	1200
<b>Readability (g)</b>	0.001	0.01
<b>Modes</b>	g, oz, ct	g, oz, ct
<b>Pan</b>	11.5 cm (4.5")	16 cm (6.3")
<b>Draft Shield</b>	Yes, 3-door glass	No
<b>Price</b>	\$1,174.00	\$950.00

Warranty: 2 years

### Analytical


Benchmark



Model	★ W3100-120	W3100A-120
<b>Calibration</b>	External (100 g weight provided)	Internal (Quick-Cal™ software)
<b>Capacity (g)</b>	120	120
<b>Readability (g)</b>	0.0001	0.0001
<b>Modes</b>	g, oz, lb, ct	g, oz, lb, ct
<b>Pan</b>	9 cm (3.5")	9 cm (3.5")
<b>Price</b>	\$1,746.00	\$1,853.00

Warranty: 2 years

## Baths



BSH1002

Aluminum  
Blocks

**Dry Baths**, digital, Benchmark (2-year warranty)

1-Block Capacity. <b>BSH1001</b> .....	\$427.00
2-Block Capacity. <b>BSH1002</b> .....	\$532.00
4-Block Capacity. <b>BSH1004</b> .....	\$869.00

*Blocks sold separately – please specify configuration*

Blocks, aluminum

Holds (48) 0.2 mL Tubes. <b>BSW02</b> .....	\$118.00
Holds (24) 0.5 mL Tubes. <b>BSW05</b> .....	\$118.00
Holds (24) 1.5 mL Tubes. <b>BSW1500</b> .....	\$118.00
Holds (24) 1.5/2.0 mL Tubes. <b>BSW1520</b> .....	\$118.00
Holds (20) 12/13 mL Tubes. <b>BSW13</b> .....	\$118.00
Holds (12) 15 mL Tubes. <b>BSW15</b> .....	\$118.00
Holds (5) 50 mL Tubes. <b>BSW50</b> .....	\$118.00



40 x 1.5/2.0 mL



15 x 50 mL



41 x 15 mL

Racks are only included with model B2000-4.

**Water Baths**, digital, MyBath™, Benchmark (2-year warranty)

4 L, 3 tube racks. <b>B2000-4</b> .....	\$832.00
8 L. <b>B2000-8</b> .....	\$916.00
12 L. <b>B2000-12</b> .....	\$1,081.00



Economy

Polypropylene  
Beakers



Premiu

### Beakers

Glass (borosilicate), Griffin low form, single scale, economy

10 mL, no graduations, each. <b>6102022910</b> .....	\$1.85
50 mL, each. <b>7-440050</b> .....	\$3.10
100 mL, each. <b>7-440100</b> .....	\$3.40
150 mL	
Pk/12. <b>61020229150</b> .....	\$46.30
Cs/48. <b>61020229150</b> .....	\$156.00
250 mL, each. <b>61020229250</b> .....	\$3.50
400 mL, each. <b>61020229400</b> .....	\$4.40
600 mL, each. <b>61020229600</b> .....	\$5.15
800 mL	
Pk/6. <b>600BG800</b> .....	\$63.45
Cs/24. <b>600BG800</b> .....	\$234.00
1000 mL, each. <b>600BG1000</b> .....	\$8.75
1000 mL, Heavy-Duty, each. <b>600BG1003</b> .....	\$14.15
2000 mL, each. <b>600BG2000</b> .....	\$16.00
4000 mL, each. <b>600BG4000</b> .....	\$44.30

*Continued*

**Beakers (cont)**

Polypropylene

Economy, Triple-Spout

- 50 mL, 5 mL graduations, Cs/100. **3640** ..... \$19.50
- 100 mL, 10 mL graduations, Cs/100. **80093** ..... \$29.70
- 250 mL, 20 mL graduations, Cs/100. **3642** ..... \$43.90
- 1000 mL, 50 mL graduations, Cs/100. **80098** ..... \$99.00

Premium, Griffin low form, single scale

- 50 mL, 10 mL graduations, each. **262100000** ..... \$3.05
- 100 mL, 25 mL graduations, each. **262110000** ..... \$2.85
- 250 mL, 50 mL graduations, each. **262130000** ..... \$5.05
- 600 mL, 100 mL graduations, each. **262150000** ..... \$8.05
- 1000 mL, 200 mL graduations, each. **262160000** ..... \$8.15
- 2000 mL, 500 mL graduations, each. **262170000** ..... \$19.70
- 4000 mL, 1000 mL graduations, each. **262190000** ..... \$40.70



**Roll Bench Liner**

Protect your bench surface using this 40 mil absorbent liner backed with 0.5 mil polyethylene film.

**Bench Liners, white**

- Roll, 50' x 20", Labmat™. **246750000** ..... \$67.65
- Sheets,
  - Bel-Art, 20 x 24", Pk/100. **246750020** ..... \$236.00
  - Ahlstrom, 46 x 57 cm (18 x 22½"), Pk/50. **4002-4657** ..... \$82.75

**Bottles**



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.

- Bottle, Balsam, glass, 45 mL (empty), USA. 616** ..... \$18.35



- Bottle Carrier, BOD, holds 8 bottles, polypropylene**
- Bx/2. **600BODCAR** ..... \$37.60




Deliver controlled drops from standard 0.030" tip opening with these LDPE bottles. Sealer cap snaps on easily.

**Bottles, Dispensing, squeeze bottles**

- LDPE (low-density polyethylene), Yorker cap
  - 30 mL (1 oz), Bg/12. **116370001** ..... \$18.40
  - 60 mL (2 oz), Bg/12. **116370002** ..... \$21.40
  - 125 mL (4 oz), Bg/12. **116370004** ..... \$35.60
  - 185 mL (6 oz), Bg/12. **116370006** ..... \$35.60

A long nozzle, textured sides, and slim profile make this chemically inert, HDPE container ideal for dispensing immersion oil.



**HDPE (high-density polyethylene)**

- 20 mL, each. **18145** ..... \$7.00

Amber glass protects light-sensitive liquids. Screw-cap dropper features a black rubber bulb and glass pipet.



**Boston Round**

**Bottles, Dropping, Amber Glass (flint)**

- 30 mL (1 oz), Bx/12. **3104053** ..... \$14.00 ↓
- Bottle Dropper Assembly (replacement)
  - 20/400, 7 x 76 mm Pipet, Pk/12. **3104053D** ..... \$8.25
- 60 mL (2 oz), Bx/12. **3104054** ..... \$16.10 ↓↓
- Bottle Dropper Assembly (replacement)
  - 20/400, 7 x 89 mm Pipet, Pk/12. **12574054D** ..... \$8.80

**French Square**



**Bottles, Dropping, Clear Glass (flint)**

- 15 mL (1/2 oz), Bx/12. **3104062** ..... \$17.10 ↓
- Bottle Dropper Assembly (replacement)
  - 20/400, 7 x 60 mm Pipet, Pk/12. **3104074** ..... \$9.80
- 30 mL (1 oz), Bx/12. **3104063** ..... \$18.50 ↓
- Bottle Dropper Assembly (replacement)
  - 24/400, 7 x 75 mm Pipet, Pk/12. **3104075** ..... \$11.65



*Continued*