



Burner Lighter, replaceable single flint, round file. **SOL-290.00** ... \$3.40
 Flint (replacement), Pk/5. **SOL-290.01** \$1.85

Carboys & Jugs



Carboys with Spigot, polypropylene, round
 10 L (2.5 gal). **118460025** \$70.80
 20 L (5 gal). **60034052** \$99.00



Jug, dispensing, HDPE jug/LDPE spigot
 20 L (5 gal). **118500000** \$46.75



Carts, utility, stainless steel, Gusto
 300 lb Capacity. **CART300S** \$287.00
 26.75"(L) x 15.5"(W) x 36"(H)
 500 lb Capacity. **CART500S** \$428.00
 29.75"(L) x 18"(W) x 36"(H)



Centrifuges, Mini
 5,500 RPM (2000 x g), 12-place, MyFuge™12. **C1012** \$374.00
 (1-year warranty)
 6,000 RPM (2000 x g), 8-place, MyFuge™Mini. **C1008-*** \$283.00
 *Specify clear (-C), green (-G), or purple (-P) for lid color.
 (1-year warranty)
 10,000 RPM (5000 x g), 6-place, blue body. **MLX-110** \$193.00
 Holds (16) 0.2 mL tubes/strips, (6) 0.5 mL tubes,
 or (6) 1.5/2.0 mL tubes with included accessories.
 (2-year warranty)

Clamps & Supports

Rectangular
Steel Base






Triangular
Cast-Iron Base

Support (Ring) Stands, finished metal base/plated support rod

Rectangular Steel Base
 6 x 4" base, 5/16 x 18" rod. **3104SSS66** \$13.15
 8 x 5" base, 3/8 x 20" rod. **3104SSS67** \$16.10
 9 x 6" base, 1/2 x 24" rod. **3104SSS68** \$18.90

Rectangular Cast-Iron Base
 Provides greater stability than steel base.
 6 x 4" base, 5/16 x 18" rod. **3104SSC14** \$14.60
 8 x 5" base, 3/8 x 20" rod. **3104SSC15** \$17.20
 9 x 6" base, 1/2 x 24" rod. **3104SSC16** \$23.40

Triangular Cast-Iron Base
 5" Legs, 1/2 x 24" rod. **1257G2324** \$25.00
 6" Legs, 1/2 x 36" rod. **1257G23** \$30.65



Support Clamps
 Burette, plastic coated round jaws, 360° rotatable. **3104BC37** \$6.90
 Burette, Double, deluxe. **600CLBURDD** \$19.00
 Plastic covered spring jaws on a metal frame securely grip
 burettes up to 30 mm OD. Fits rods up to 3/4" (19 mm) OD.
 Hook Collar, nickel plated. **1257G49** \$5.65
 Fits rods up to 1/2" (13 mm) in diameter.
 Right-Angle, aluminum. **97-4023** \$6.10
 Fits rods up to 5/8" in diameter.
 Thermometer, aluminum, rotates 360°. **600THCL01** \$7.75
 Fits rods up to 1/2" (13 mm) in diameter.



Support Rings, plated cast-iron
 3" OD (2 1/2" ID). **600SRC103** \$5.95
 4" OD (3 1/2" ID). **600SRC104** \$7.80
 5" OD (4 1/2" ID). **600SRC105** \$8.25

Continued