

Stereo Microscopes



420T-430PLL-10

National

420-430PLL-10 Stereoscope, Binocular, 10-40X Zoom \$569.00

- Widefield 10X/20 eyepieces; accept 23 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 10X-40X total magnification; working distance of 79 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual illuminators offer light from top only (LED with rheostat control and directional beam), bottom only (LED), or both top and bottom
- Post stand with locking support collar allows viewing specimens up to 93 mm tall
- Frosted glass and reversible black/white plastic stage plates included
- 24.2 cm(D) x 17.1 cm(W) x 43.3 cm(H)
- Lifetime mechanical and optical/3-year electrical warranty

420T-430PLL-10 Stereoscope, Trinocular, 10-40X Zoom \$659.00

- Trinocular version of 420-430PLL-10
- 2X SLR adapter & 0.4X C-mount adapter included

425-445PLL-10



National

425-445PLL-10 Stereoscope, Binocular, 7.5-50X Zoom \$914.00

- Widefield 10X/23 eyepieces; accept 24.95 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 7.5-50X total magnification; working distance of 113 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual 3W LED illuminators with rheostat control offer cool light from top, bottom, or both top and bottom
- Post stand has a large ergonomic base and locking support collar
- Frosted glass and reversible black/white plastic stage plates included
- Lifetime mechanical and optical/3-year electrical warranty


425T-445PLL-10 Stereoscope, Trinocular, 7.5-50X Zoom ... \$973.00

- Trinocular version of 425-445PLL-10



Meiji Access www.meijitechno.com
and call us to place an order

Stands



Please call regarding availability of items not listed

Stereo Boom Stands, Diagnostic Instruments

All boom stands require a focus mount

- | | |
|---|-----------------------|
| Standard, nylon bearings. SMS16B | \$599.00 |
| Ball Bearing | |
| Medium-Duty. SMS6B | \$759.00 |
| Heavy-Duty. SMS20-__ * | <i>See note below</i> |
| Articulating Arm. SMS25 | \$1,409.00 |

*Access www.spotimaging.com and call us with microscope model and boom stand version to place an order



Meiji PBH Stand
with Adapter Sleeve

Stereo Dual Illumination Stands, 6V/10W, Meiji

- | | |
|---|----------|
| Pole-Type Base, 84 mm focus mount. PBH | \$563.00 |
| <i>Some stereo brands require an adapter sleeve</i> | |
| Stereo Stand Sleeve for 76 mm dia. bodies | \$30.00 |
| Stereo Stage Plates, 94.5 mm | |
| Clear Glass. MA569 | \$45.00 |
| Filter, Blue Frosted, 40 mm. MA563/05 | \$50.00 |
| Black & White Plastic. MA568 | \$45.00 |