

Funnels (cont)

- Standard Stem, Premium, polypropylene, Bel-Art
 Exclusive filter-step internal design provides speedier filtration; external ribs prevent air-lock.
- 55 mm Top ID, takes 9 cm filter, Pk/12. **146760055** \$47.75
 - 65 mm Top ID, takes 11 cm filter, Pk/12. **146760065** \$49.60
 - 75 mm Top ID, takes 12.5 cm filter, Pk/6. **146760075** \$25.65
 - 100 mm Top ID, takes 18.5 cm filter, Pk/4. **146760100** \$32.60
- Standard Stem, Economy, polypropylene
- 50 mm Top ID, each. **310FUNPP50** \$0.80
 - 62 mm Top ID, each. **310FUNPP62** \$0.95
 - 75 mm Top ID, each. **310FUNPP75** \$1.15
 - 100 mm Top ID, each. **310FUNPP100** \$1.65

These 4 x 4" nonsterile gauze sponges are 12-ply.



Gauze Pads, 10 pk/200, Cs/2000. 44122000 \$74.20

Gloves & Hand Protection

Cotton Gloves



Heat-Resistant Gloves



Gloves


- Cotton
 Protect optics and photographs from fingerprints. Can also be used as glove liners.
- Small, Pk/12. **GLV-190.10** \$4.95
 - Large, Pk/12. **GLV-190.20** \$5.25
- Heat-Resistant
 Terry cloth gloves protect hands from high temperatures.
- 33 cm (13") Long, Pair. **132010000** \$35.10
 - 47 cm (18.5") Long, Pair. **132010001** \$47.75



Safely grip hot or cold labware with this silicone rubber hand protector partially lined with mini suction cups. Withstands temperatures from 260°C (500°F) to -57°C (-70°F).

Hand Protector, red. 38000000 \$17.90

Homogenizers & Tissue Grinder



Rapidly homogenize, emulsify, suspend and disrupt 0.1 mL to 250 mL samples. Includes a 5 x 50 mm flat-bottom probe & 7 x 50 mm saw tooth probe for use with micro tubes.

Homogenizer, Handheld. D1000 \$1,020.00



Homogenize a variety of sample types, from soft tissue to bone, in 35 seconds or less.
 Benchmark BeadBlaster™ 24

Homogenizer, Benchtop, Microtube. D2400 \$7,049.00

- Holds (24) 2.0 mL tubes, or (6) 5 mL tubes w/optional tube holder
- Speed adjusts from 4.00 to 7.00 m/s in 0.05 m/s increments
- Stores 50 programs in memory, each with a maximum of 10 cycles
- 2-year warranty

Homogenizer Prefilled Tubes

2.0 mL screw-capped tubes are prefilled with your choice of the following beads:

Glass, acid washed


- 0.1 mm, Pk/50. **D1031-01** \$113.00
- 0.5 mm, Pk/50. **D1031-05** \$113.00

Zirconium Beads, acid washed/heat treated/nuclease-free

- 1.5 mm, Pk/50. **D1032-15** \$158.00
- 3.0 mm, Pk/50. **D1032-30** \$158.00

Stainless Steel Beads, acid washed

- 2.8 mm, Pk/50. **D1033-28** \$158.00



Tissue Grinder, Dual® Glass Pestle

- 1 mL Capacity, **LG10640100** \$85.15
- 3 mL Capacity, **LG10640102** \$95.80
- 5 mL Capacity, **LG10640104** \$109.00