



Dimensions: 29 x 8.8 x 40 cm

Pipet Rack, magnetized, blue, Heathrow. **HS20612M** \$76.10

Pipets, Mohr (Measuring), Class A, individually certified


- 5 mL, (± 0.020 mL), Pk/5. **600PM7058H**..... \$67.05
- 10 mL, (± 0.030 mL), Pk/5. **600PM7058I**..... \$67.05
- 25 mL, (± 0.050 mL), Pk/5. **600PM7058J**..... \$99.90



Dispense fluids accurately and safely with these color-coded and plugged pipets. Negative graduations provide extra volume (not available on 25 mL). RNase-/DNase-/Pyrogen-free, and noncytotoxic. Calibrated "To Deliver" (TD).

Pipets, Serological, polystyrene, UniPlast™

- Nonsterile**
 - 5 mL, blue band, pk/25, Cs/375. **1755**..... \$83.45
 - 10 mL, orange band, pk/25, Cs/250. **1775**..... \$57.15
 - 25 mL, red band, pk/10, Cs/100. **1784**..... \$53.20
- Sterile Pk/1** (individually paper/plastic wrapped)
 - 1 mL, yellow band, Cs/500. **1700**..... \$83.80
 - 2 mL, green band, Cs/500. **1720**..... \$89.75
 - 5 mL, blue band, reverse grad., Cs/250. **1740**..... \$66.30
 - 10 mL, orange band, reverse grad., Cs/200. **1760**..... \$57.80
 - 25 mL, red band, reverse grad., Cs/100. **1780**..... \$62.50
 - 50 mL, black band, reverse grad., Cs/100. **1820**..... \$157.00
- Sterile Bg/25**
 - 1 mL, yellow band, Cs/1000. **1710**..... \$153.00
 - 10 mL, orange band, Cs/250. **1770**..... \$68.00
 - 25 mL, red band, Cs/100. **1782**..... \$55.35



Pipets, Transfer, LDPE (low-density polyethylene)

- Graduated, Nonsterile, 15 cm long**
 - 5 mL capacity, 1 mL graduated, 25 drops/mL
 - #1** Bx/500. **PTP-02**..... \$17.50
 - Cs/10 Bx. **PTP-02**..... \$168.00
 - 5 mL capacity, thick-walled
 - Cs/10 Bx. **1881200**..... \$140.00
 - 7 mL capacity, 3 mL graduated, 14 drops/mL
 - Bx/500. **PTP-01**..... \$17.50
 - Cs/10 Bx. **PTP-01**..... \$168.00
- Graduated, Sterile Pk/1** (individually wrapped), 15 cm long
 - 5 mL capacity, 1 mL graduated, Simport
 - Bx/400. **P200-521S**..... \$43.30
 - Cs/10 Bx. **P200-521S**..... \$398.00
- Graduated, Sterile Pk/20**, 15 cm long
 - 5 mL capacity, 1 mL graduated, Bx/500. **HS206371B**..... \$37.75

- Long Stem, Sterile Pk/1 (individually wrapped)
 - 22.5 cm (9") long, 6 mL capacity, Cs/400. **139038**..... \$58.20
- Long Stem, Nonsterile
 - 30 cm (12") long, 23 mL capacity, Bx/100. **139050**..... \$21.30

Made of borosilicate glass in compliance with ASTM specification E-969, Class A requirements. Individually serialized and certified. Calibrated "To Deliver" (TD).



Pipets, Volumetric, Class A, individually certified

- 2 mL (± 0.006 mL), orange band, Pk/5. **600PT7100B**..... \$68.30
- 5 mL (± 0.010 mL), white band, Pk/5. **600PT7100C**..... \$71.75
- 10 mL (± 0.020 mL), red band, Pk/5. **600PT7100D**..... \$81.95
- 20 mL (± 0.030 mL), yellow band, Pk/5. **600PT7100E**..... \$116.00



These mini, air displacement pipettors provide an economical and accurate method for pipetting fixed volumes. Pipettors feature a two-step plunger and are fully autoclavable.

Pipettors, Single Channel

- Fixed Volume, pocket-size, Diamond Jr.**
 - 5 μ L, red. **3322-5**..... \$22.50
 - 10 μ L, pink. **3322-10**..... \$22.50
 - 20 μ L, green. **3322-20**..... \$22.50
 - 50 μ L, purple. **3322-50**..... \$22.50
 - 100 μ L, orange. **3322-100**..... \$22.50



- Adjustable Volume, Diamond Advance™** (3-year warranty)
 - 0.5-10 μ L, 0.02 μ L increments, red. **3341-10**..... \$225.00
 - 5-50 μ L, 0.1 μ L increments, yellow. **3341-50**..... \$225.00
 - 10-100 μ L, 0.1 μ L increments, blue. **3341-100**..... \$225.00
 - 100-1000 μ L, 1 μ L increments, orange. **3341-1000**..... \$225.00



Pipettor Holders, Carousel, Flip & Grip™, Heathrow

- Blue. **HS120026**..... \$80.50
- Lime. **HS120056**..... \$80.50

Continued

Continued