





Pipettor Tips, Filtered (Barrier), sterile *Inquire*
 1000 μ L, 96/rack, Cs/35 rks (3360). **BP71008B** **NEW** \$223.00

Pipettor Tips, Standard, graduated

Nonsterile, natural
 0.1-20 μ L, 45 mm, Bg/1000. **151145** \$20.30
 Sterile, blue
 100-1250 μ L, 6 racks of 96, Bx/576. **151153BRS-96** \$45.00



Reagent Reservoir, 50 mL, v-bottom, Heathrow
 Sterile Bg/1, white polystyrene, Pk/100. **120644** \$94.90



Polypropylene



Stainless Steel

Pitchers

Polypropylene, graduated (mL), tall form, tapered
 1000 mL, 50 mL graduations, each. **289910000** \$6.20
 2000 mL, 100 mL graduations, each. **289920000** \$7.10
 3000 mL, 100 mL graduations, each. **289930000** \$13.80
 Stainless Steel, graduated (mL/oz), tapered, Polar Ware
 1000 mL/32 oz, each. **T1063** \$73.20

Rockers, Shakers & Mixers



Rockers, Fixed Speed, Benchmark (2-year warranty)

2-D (seesaw), 12 x 12" platform w/flat mat
 BenchBlotter™. **BR1000** \$460.00

- Fixed speed of 12 rocks/minute
- Platform angle is $\pm 8^\circ$
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)

3-D (nutating)
 Small, 10.5 x 7.5" platform w/dimpled mat
 MiniMixer™. **B3D1020** \$393.00

Large, 12 x 12" platform w/dimpled mat
 BioMixer™. **B3D1320** \$474.00

- Fixed speed of 24 rocks/minute
- Platform angle is $\pm 20^\circ$
- Cold room & incubator safe (4-45°C)
- Load capacity is 1.8 kg (4 lbs)



Rockers, Variable Speed

Benchmark (2-year warranty)

2-D (seesaw), 14 x 12" platform w/flat mat
 BenchRocker™ 2D. **BR2000** \$732.00

Stacking Platform. **BR2000-STACK** \$121.00
 14 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP** \$43.00
 Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 3.6 kg (8 lbs)

3-D (nutating), 12 x 12" platform w/dimpled and flat mats
 BenchRocker™ 3D. **B3D2300** \$860.00

Stacking Platform. **BR1000-STACK** \$111.00
 12 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP** \$43.00
 Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)

Continued