





Pipettor Tips, Filtered (Barrier), sterile *Inquire*
 1000 µL, 96/rack, Cs/35 rks (3360). **BP71008B** **NEW** \$223.00

Pipettor Tips, Standard, graduated

Nonsterile, natural
 0.1-20 µL, 45 mm, Bg/1000. **151145** \$20.30
 Sterile, blue
 100-1250 µL, 6 racks of 96, Bx/576. **151153BRS-96** \$45.00



Reagent Reservoir, 50 mL, v-bottom, Heathrow
 Sterile Bg/1, white polystyrene, Pk/100. **120644** \$94.90



Polypropylene




Stainless Steel

Pitchers

Polypropylene, graduated (mL), tall form, tapered
 1000 mL, 50 mL graduations, each. **289910000** \$6.20
 2000 mL, 100 mL graduations, each. **289920000** \$7.10
 3000 mL, 100 mL graduations, each. **289930000** \$13.80
 Stainless Steel, graduated (mL/oz), tapered, Polar Ware
 1000 mL/32 oz, each. **T1063** \$73.20

Rockers, Shakers & Mixers



Rockers, Fixed Speed, Benchmark (2-year warranty)

2-D (seesaw), 12 x 12" platform w/flat mat
 BenchBlotter™. **BR1000** \$460.00
 • Fixed speed of 12 rocks/minute
 • Platform angle is ±8°
 • Cold room & incubator safe (4-45°C)
 • Load capacity is 2.3 kg (5 lbs)
 3-D (nutating)
 Small, 10.5 x 7.5" platform w/dimpled mat
 MiniMixer™. **B3D1020** \$393.00
 Large, 12 x 12" platform w/dimpled mat
 BioMixer™. **B3D1320** \$474.00
 • Fixed speed of 24 rocks/minute
 • Platform angle is ±20°
 • Cold room & incubator safe (4-45°C)
 • Load capacity is 1.8 kg (4 lbs)



BR2000



B3D2300

Rockers, Variable Speed

Benchmark (2-year warranty)
 2-D (seesaw), 14 x 12" platform w/flat mat
 BenchRocker™ 2D. **BR2000** \$732.00
 Stacking Platform. **BR2000-STACK** \$121.00
 14 x 12" platform w/flat mat has 3" of clearance.
 Stacking Platform Extenders, Pk/4. **BR2000-SP** \$43.00
 Adds 1.25" to platform separation.
 • Variable speed of 2-30 rocks/minute
 • Platform tilt angle adjusts from 0-30°
 • Cold room & incubator safe (4-45°C)
 • Load capacity is 3.6 kg (8 lbs)
 3-D (nutating), 12 x 12" platform w/dimpled and flat mats
 BenchRocker™ 3D. **B3D2300** \$860.00
 Stacking Platform. **BR1000-STACK** \$111.00
 12 x 12" platform w/flat mat has 3" of clearance.
 Stacking Platform Extenders, Pk/4. **BR2000-SP** \$43.00
 Adds 1.25" to platform separation.
 • Variable speed of 2-30 rocks/minute
 • Platform tilt angle adjusts from 0-30°
 • Cold room & incubator safe (4-45°C)
 • Load capacity is 2.3 kg (5 lbs)

Continued



Rockers, Variable Speed

Reliable Scientific (5-year warranty), **USA**

2-D (seesaw)

Standard, 11 x 14" platform w/rubber mat

Single Platform. **55S1114**.....\$603.00

Double Platform, 3/4" apart. **55D1114**.....\$665.00

Large, 12 x 16" platform w/rubber mat

Single Platform. **55S1216**.....\$647.00

Double Platform, 3/4" apart. **55D1216**.....\$704.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



Rocker, Tube, 12 x 4.5" grooved mat, Benchmark. **M2100** \$348.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is ±30°
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs)
- 2-year warranty

Orbital Shakers



BT30

Aerate and mix with these low to high capacity orbital shakers.



BT302

Orbital Shakers, Benchmark (2-year warranty)

Low Speed, 14 x 12" platform w/flat rubber mat

Orbi-Blotter™. **BT30**.....\$786.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

Small, 9.5 x 8" platform w/flat rubber mat

Orbi-Shaker™ JR. **BT302**\$909.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min



BT3001

Large, 13 x 12" platform w/flat rubber mat

Orbi-Shaker™. **BT3001**.....\$1,499.00

- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L, or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR**...\$255.00

250 mL Flask Magnetic Clamp. **H1000-MR-250** \$36.00

500 mL Flask Magnetic Clamp. **H1000-MR-500** \$36.00

1000 mL Flask Magnetic Clamp. **H1000-MR-1000** \$55.00



R4045

R4040

Tube Rotator, programmable, Benchmark (2-year warranty)

Roto-Bot Mini. **R4040**\$550.00

Roto-Bot. **R4045**.....\$827.00

- Speed adjusts from 1-80 rpm
- Tube capacity (starter set)
 - R4040: 12 x 1.5 mL, 2 x 15 mL & 2 x 50 mL
 - R4045: 30 x 1.5 mL, 8 x 15 mL & 2 x 50 mL
- Timer is programmable from 1 sec – 99 hrs, 59 min / Cont.


Vortexers



VM-03UA



BV101



BV1000

Vortex Mixers, Benchmark (2-year warranty)

MiniMixer, fixed 3000 rpm, blue cup. **VM-03UA**.....\$136.00

Vornado™, mini vortexer, fixed 2800 rpm. **BV101**__*\$199.00

**Specify Gray, Blue or Purple cup head*

BenchMixer™, single cup head, 200-3200 rpm. **BV1000**\$254.00

Flat Head Attachment, 3". **BV1000-FLAT**\$36.50

Continued