



**Rockers, Variable Speed**

Reliable Scientific (5-year warranty), **USA**

2-D (seesaw)

Standard, 11 x 14" platform w/rubber mat

Single Platform. **55S1114**.....\$603.00

Double Platform, 3/4" apart. **55D1114**.....\$665.00

Large, 12 x 16" platform w/rubber mat

Single Platform. **55S1216**.....\$647.00

Double Platform, 3/4" apart. **55D1216**.....\$704.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



**Rocker, Tube**, 12 x 4.5" grooved mat, Benchmark. **M2100** ..... \$348.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is ±30°
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs)
- 2-year warranty

**Orbital Shakers**



BT30

Aerate and mix with these low to high capacity orbital shakers.



BT302

**Orbital Shakers**, Benchmark (2-year warranty)

Low Speed, 14 x 12" platform w/flat rubber mat

Orbi-Blotter™. **BT30**.....\$786.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

Small, 9.5 x 8" platform w/flat rubber mat

Orbi-Shaker™ JR. **BT302** .....\$909.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min



BT3001

Large, 13 x 12" platform w/flat rubber mat

Orbi-Shaker™. **BT3001**.....\$1,499.00

- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L, or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

*Optional Magic Clamp™ platform & accessories sold separately - specify configuration*

Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR**...\$255.00

250 mL Flask Magnetic Clamp. **H1000-MR-250** ..... \$36.00

500 mL Flask Magnetic Clamp. **H1000-MR-500** ..... \$36.00

1000 mL Flask Magnetic Clamp. **H1000-MR-1000** ..... \$55.00



R4045

R4040

**Tube Rotator**, programmable, Benchmark (2-year warranty)

Roto-Bot Mini. **R4040** .....\$550.00

Roto-Bot. **R4045**.....\$827.00

- Speed adjusts from 1-80 rpm
- Tube capacity (starter set)
  - R4040: 12 x 1.5 mL, 2 x 15 mL & 2 x 50 mL
  - R4045: 30 x 1.5 mL, 8 x 15 mL & 2 x 50 mL
- Timer is programmable from 1 sec – 99 hrs, 59 min / Cont.


**Vortexers**



VM-03UA



BV101



BV1000

**Vortex Mixers**, Benchmark (2-year warranty)

MiniMixer, fixed 3000 rpm, blue cup. **VM-03UA**.....\$136.00

Vornado™, mini vortexer, fixed 2800 rpm. **BV101**\_\_\* .....\$199.00

*\*Specify Gray, Blue or Purple cup head*

BenchMixer™, single cup head, 200-3200 rpm. **BV1000** .....\$254.00

Flat Head Attachment, 3". **BV1000-FLAT** .....\$36.50

*Continued*