


## Rulers, Meter Sticks & Tape Measure

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. **4461** ..... \$0.95

White. **4460** ..... \$0.55

Stainless Steel, 12 mm wide, 0.5 mm thick, **USA. 605-065** ..... \$5.60

Meter Sticks, lacquered wood

One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

Plain Ends. **605-050** ..... \$4.90

Metal Ends. **605-051** ..... \$6.15

Tape Measure, 152.4 cm/60", reinforced fiberglass. **11036** ..... \$3.70

One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.



Sharps Containers, red polypropylene

1 Quart, translucent locking lid. **440065** ..... \$3.80  
Dimensions: 4.25" x 4.25" x 6.25"(H)

5.4 Quart, horizontal entry lid, Prevent™, each. **4402269** ..... \$8.35  
Dimensions: 12" x 4.75" x 11"(H)

12 Gallon, hinged vertical entry lid, each. **440321391** ..... \$28.40  
Dimensions: 18.25" x 12.75" x 18.75"(H)

Screw-Cap, 90 mL



Screw-Cap,  
480 & 960 mL



Snap-Cap

Specimen Containers, polypropylene (PP)

Screw-Cap

90 mL, sterile, Cs/400. **840175** ..... \$110.00

480 mL (16 oz), Cs/24. **59014510** ..... \$71.95

83 mm (3 1/4") opening

960 mL (32 oz), Cs/12. **59014511** ..... \$59.05

111 mm (4 3/8") opening

Snap-Cap, white HDPE bottom/white LDPE lid

480 mL (16 oz), Cs/100. **271016W** ..... \$62.85

9.9 cm opening x 11.1 cm high

1000 mL (34 oz), Cs/100. **271032W** ..... \$92.30

11.7 cm opening x 13.7 cm high.



Spot Plates, glazed white porcelain, 21 mm diameter wells

6-Well, each. **LPC-161** ..... \$2.35

12-Well, 115 x 95 x 7 mm, each. **3104187701** ..... \$3.85

## Stirring Bars



FDA approved Teflon® (PTFE) coating is a minimum of 1.6 mm (0.062") thick, except for micro bars which have a 1.32 mm (0.050") thick cover.

Stirring Bars, Teflon®, Bel-Art

Centered Sphere, Saturn Spinbar®

40 mm (1 1/2") Long. **371160002** ..... \$26.00

50 mm (2") Long. **37116003** ..... \$38.80

Micro, Spinbar®

White

5 x 2 mm. **371190005** ..... \$5.15

7 x 2 mm. **371190007** ..... \$5.15

10 x 3 mm. **371190010** ..... \$5.75

12.7 x 3 mm. **371190127** ..... \$5.75

15 x 1.5 mm. **371190015** ..... \$5.05

Red

7 x 2 mm. **371210013** ..... \$5.75

8 x 1.5 mm. **371210001** ..... \$5.75

Octagonal, molded-on pivot ring, Spinbar®

White

12.7 x 8 mm (1/2 x 5/16"). **371101256** ..... \$4.85

15.9 x 8 mm (5/8 x 5/16"). **371100058** ..... \$4.85

25.4 x 9.5 mm (1 x 3/8"). **371100138** ..... \$6.40

38.1 x 9.5 mm (1 1/2 x 3/8"). **371101128** ..... \$7.60

50.8 x 9.5 mm (2 x 3/8"). **371100238** ..... \$10.10

63.5 x 9.5 mm (2 1/2 x 3/8"). **371102128** ..... \$11.70

76.2 x 12.7 mm (3 x 1/2"). **371100003** ..... \$21.00

Red, 12.7 x 3.2 mm (1/2 x 1/8"), no pivot ring. **371090034** ..... \$5.95

Stirring Bar Assortment, octagonal, Pk/6. **371600000** ..... \$37.40

Includes two of each size: 25.4 x 8 mm (1 x 5/16"),

38.1 x 8 mm (1 1/2 x 5/16"), and 50.8 x 8 mm (2 x 5/16").

Polygonal, molded-on pivot ring, Spinbar®

12 x 4.5 mm (7/16 x 3/16"). **371220012** ..... \$5.10

Triangular, vanes w/cross-mounted magnets, Spinvane®

Fits 1 mL Vials. **371340000** ..... \$10.05

Fits 3/5 mL Vials. **371350000** ..... \$10.05

