


Rulers, Meter Sticks & Tape Measure

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

- Transparent. **4461** \$0.95
- White. **4460** \$0.55
- Stainless Steel, 12 mm wide, 0.5 mm thick, **USA. 605-065** \$5.60

Meter Sticks, lacquered wood

One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

- Plain Ends. **605-050** \$4.90
- Metal Ends. **605-051** \$6.15

Tape Measure, 152.4 cm/60", reinforced fiberglass. **11036** \$3.70


One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.



Sharps Containers, red polypropylene

- 1 Quart, translucent locking lid. **440065** \$3.80
Dimensions: 4.25" x 4.25" x 6.25"(H)
- 5.4 Quart, horizontal entry lid, Prevent™, each. **4402269** \$8.35
Dimensions: 12" x 4.75" x 11"(H)
- 12 Gallon, hinged vertical entry lid, each. **440321391** \$28.40
Dimensions: 18.25" x 12.75" x 18.75"(H)

Screw-Cap, 90 mL



Specimen Containers, polypropylene (PP)

Screw-Cap

- 90 mL, sterile, Cs/400. **840175** \$110.00
- 480 mL (16 oz), Cs/24. **59014510** \$71.95
83 mm (3 1/4") opening
- 960 mL (32 oz), Cs/12. **59014511** \$59.05
111 mm (4 3/8") opening
- Snap-Cap, white HDPE bottom/white LDPE lid
- 480 mL (16 oz), Cs/100. **271016W** \$62.85
9.9 cm opening x 11.1 cm high
- 1000 mL (34 oz), Cs/100. **271032W** \$92.30
11.7 cm opening x 13.7 cm high.



Spot Plates, glazed white porcelain, 21 mm diameter wells

- 6-Well, each. **LPC-161** \$2.35
- 12-Well, 115 x 95 x 7 mm, each. **3104187701** \$3.85

Stirring Bars



FDA approved Teflon® (PTFE) coating is a minimum of 1.6 mm (0.062") thick, except for micro bars which have a 1.32 mm (0.050") thick cover.

Stirring Bars, Teflon®, Bel-Art

Centered Sphere, Saturn Spinbar®

- 40 mm (1 1/2") Long. **371160002** \$26.00
- 50 mm (2") Long. **37116003** \$38.80

Micro, Spinbar®

- White
- 5 x 2 mm. **371190005** \$5.15
- 7 x 2 mm. **371190007** \$5.15
- 10 x 3 mm. **371190010** \$5.75
- 12.7 x 3 mm. **371190127** \$5.75
- 15 x 1.5 mm. **371190015** \$5.05

Red

- 7 x 2 mm. **371210013** \$5.75
- 8 x 1.5 mm. **371210001** \$5.75

Octagonal, molded-on pivot ring, Spinbar®

- White
- 12.7 x 8 mm (1/2 x 5/16"). **371101256** \$4.85
- 15.9 x 8 mm (5/8 x 5/16"). **371100058** \$4.85
- 25.4 x 9.5 mm (1 x 3/8"). **371100138** \$6.40
- 38.1 x 9.5 mm (1 1/2 x 3/8"). **371101128** \$7.60
- 50.8 x 9.5 mm (2 x 3/8"). **371100238** \$10.10
- 63.5 x 9.5 mm (2 1/2 x 3/8"). **371102128** \$11.70
- 76.2 x 12.7 mm (3 x 1/2"). **371100003** \$21.00

Red, 12.7 x 3.2 mm (1/2 x 1/8"), no pivot ring. **371090034** \$5.95

Stirring Bar Assortment, octagonal, Pk/6. **371600000** \$37.40

Includes two of each size: 25.4 x 8 mm (1 x 5/16"), 38.1 x 8 mm (1 1/2 x 5/16"), and 50.8 x 8 mm (2 x 5/16").

Polygonal, molded-on pivot ring, Spinbar®

- 12 x 4.5 mm (7/16 x 3/16"). **371220012** \$5.10

Triangular, vanes w/cross-mounted magnets, Spinvane®

- Fits 1 mL Vials. **371340000** \$10.05
- Fits 3/5 mL Vials. **371350000** \$10.05





Stirring Bar Retainer, Spinbar® Restrainer. 377861000..... \$14.25


Stirring Bar Retrievers, magnetic, Spinbar®

30.5 cm (12") Long. **377720000** \$15.75
 45.7 cm (18") Long. **377730000** \$18.20



Prevent stirring bars from going down the sink. Strainer's four strong magnets capture stirring bars, but fluids drain easily. Blue.

Stirring Bar Sink Strainer, magnetic, Spinbar®. 377870000..... \$30.60



Stirring Rods, glass, annealed ends

150 x 5 mm (6 x 3/16"), Pk/12. **600GSR006** \$4.45
 Rubber Policeman, fits 5 mm rod, Pk/12. **201-BB**..... \$6.00
 254 x 6 mm (10 x 1/4"), Pk/12. **600GSR010** \$6.30
 305 x 10 mm (12 x 3/8"), Pk/12. **600GSR012**..... \$23.60

Stoppers



Stoppers, rubber, black

Stopper Sizing Guide

Size #	Top (mm)	Bottom (mm)	Fits ID Openings (mm)
1	19	14	15-17
2	20	16	16-18.5
3	24	18	18-21
4	26	20	20-23
5	27	23	23-25
6	32	26	26-27
6½	34	27	27-31.5
7	37	30	30-34
8	41	33	33-37
9	45	37	37-41
10	50	42	42-46
12	64	54	54-59

All stoppers are 25 mm long. Stopper size #1 has a 4 mm hole; sizes #2 and larger have 5 mm holes.

Solid	<u>pk/5</u>	<u>lb</u>
Size 2, ~39 pcs/lb. 50020	\$3.15	13.45
Size 3, ~29 pcs/lb. 50030	\$4.25	13.45
Size 4, ~23 pcs/lb. 50040	\$5.35	13.45
Size 5, ~20 pcs/lb. 50050	—	13.45
Size 6, ~15 pcs/lb. 50060	—	13.45
Size 6.5, ~14 pcs/lb. 500650	—	13.45
Size 7, ~11 pcs/lb. 50070	—	13.45
Size 8, ~9 pcs/lb. 50080	—	13.45
1-Hole	<u>pk/5</u>	<u>lb</u>
Size 1, ~51 pcs/lb. 50011	\$2.40	13.45
Size 2, ~43 pcs/lb. 50021	\$2.85	13.45
★ Size 3, ~30 pcs/lb. 50031	\$4.10	13.45
Size 4, ~24 pcs/lb. 50041	\$5.15	13.45
Size 5, ~21 pcs/lb. 50051	—	13.45
Size 6, ~18 pcs/lb. 50061	—	13.45
Size 9, ~8 pcs/lb. 50091	—	13.45
Size 12, ~3 pcs/lb. 50121	—	13.45
2-Hole	<u>lb</u>	
Size 6, ~17 pcs/lb. 50062	\$13.45	
Size 7, ~12 pcs/lb. 50072	\$13.45	
Size 8, ~9 pcs/lb. 50082	\$13.45	
Size 9, ~7 pcs/lb. 50092	\$13.45	
Size 10, ~5 pcs/lb. 50102	\$13.45	
Assortment, 1 lb. 50012	\$17.45	

Solid, 1-hole, and 2-hole stoppers in sizes 1-6.