



Stirring Bar Retainer, Spinbar® Restrainer. 377861000..... \$14.25


Stirring Bar Retrievers, magnetic, Spinbar®

30.5 cm (12") Long. **377720000** \$15.75
 45.7 cm (18") Long. **377730000** \$18.20



Prevent stirring bars from going down the sink. Strainer's four strong magnets capture stirring bars, but fluids drain easily. Blue.

Stirring Bar Sink Strainer, magnetic, Spinbar®. 377870000 \$30.60



Stirring Rods, glass, annealed ends

150 x 5 mm (6 x 3/16"), Pk/12. **600GSR006** \$4.45
 Rubber Policeman, fits 5 mm rod, Pk/12. **201-BB** \$6.00
 254 x 6 mm (10 x 1/4"), Pk/12. **600GSR010** \$6.30
 305 x 10 mm (12 x 3/8"), Pk/12. **600GSR012** \$23.60

Stoppers



Stoppers, rubber, black

Stopper Sizing Guide

Size #	Top (mm)	Bottom (mm)	Fits ID Openings (mm)
1	19	14	15-17
2	20	16	16-18.5
3	24	18	18-21
4	26	20	20-23
5	27	23	23-25
6	32	26	26-27
6½	34	27	27-31.5
7	37	30	30-34
8	41	33	33-37
9	45	37	37-41
10	50	42	42-46
12	64	54	54-59

All stoppers are 25 mm long. Stopper size #1 has a 4 mm hole; sizes #2 and larger have 5 mm holes.

Solid	pk/5	lb
Size 2, ~39 pcs/lb. 50020	\$3.15	13.45
Size 3, ~29 pcs/lb. 50030	\$4.25	13.45
Size 4, ~23 pcs/lb. 50040	\$5.35	13.45
Size 5, ~20 pcs/lb. 50050	—	13.45
Size 6, ~15 pcs/lb. 50060	—	13.45
Size 6.5, ~14 pcs/lb. 500650	—	13.45
Size 7, ~11 pcs/lb. 50070	—	13.45
Size 8, ~9 pcs/lb. 50080	—	13.45
1-Hole	pk/5	lb
Size 1, ~51 pcs/lb. 50011	\$2.40	13.45
Size 2, ~43 pcs/lb. 50021	\$2.85	13.45
★ Size 3, ~30 pcs/lb. 50031	\$4.10	13.45
Size 4, ~24 pcs/lb. 50041	\$5.15	13.45
Size 5, ~21 pcs/lb. 50051	—	13.45
Size 6, ~18 pcs/lb. 50061	—	13.45
Size 9, ~8 pcs/lb. 50091	—	13.45
Size 12, ~3 pcs/lb. 50121	—	13.45
2-Hole	lb	
Size 6, ~17 pcs/lb. 50062	\$13.45	
Size 7, ~12 pcs/lb. 50072	\$13.45	
Size 8, ~9 pcs/lb. 50082	\$13.45	
Size 9, ~7 pcs/lb. 50092	\$13.45	
Size 10, ~5 pcs/lb. 50102	\$13.45	
Assortment, 1 lb. 50012	\$17.45	

Solid, 1-hole, and 2-hole stoppers in sizes 1-6.



Accurately determine the correct size of rubber stopper, size 00-12, or cork, size 0-26, for any vessel using this durable plastic tool. Reverse side has a circular scale for sizing 4-52 mm OD tubing.

Stopper/Tubing Sizer. 133920000 \$14.25

Non-fenestrated towels are sterile, latex-free and poly-lined.



Surgical Drape, 18 x 26", Bx/50. 44025-515..... \$25.50

Tape

Kraft paper with rubber resin adhesive.



Tape, Label, Write-On™, Bel-Art
1" Core, 15 yd (540") roll

½" Wide, each \$5.80

Green. **134623005** Yellow. **134622005**

¾" Wide, Pk/4 \$25.00

Green. **134820075** Yellow. **134810075**

White. **134800075**

3" Core, 40 yd roll

½" Wide, each \$10.55

Green. **134633005** White. **134630005**


¾" Wide, each \$14.00

Green. **134633075** Yellow. **134632075**

¾" Wide, White, Pk/4. **134830075** \$56.10

¾" Wide, Mixed, Pk/6. **134630600** \$85.35

blue, green, orange, red, white & yellow



Label Tape Dispenser, multi-roll. 134730000..... \$44.40



Versatile multi-format block holds (96) 0.2 mL tubes/strips, (39) 0.5 mL tubes, or (1) nonskirted 96-well PCR plate. In-situ adapter sold separately. Touchscreen control and height adjustable heated lid. 2-year warranty.




0.2 mL Tubes



96-Well Plate



0.5 mL Tubes



In-Situ Adapter

Thermal Cycler, multi-format block, Benchmark. T5000-96 .. \$4,471.00
In-Situ Adapter, holds 4 slides. **T5000-A-IS** \$439.00