



Accurately determine the correct size of rubber stopper, size 00-12, or cork, size 0-26, for any vessel using this durable plastic tool. Reverse side has a circular scale for sizing 4-52 mm OD tubing.

**Stopper/Tubing Sizer. 133920000 ..... \$14.25**

Non-fenestrated towels are sterile, latex-free and poly-lined.



**Surgical Drape, 18 x 26", Bx/50. 44025-515..... \$25.50**

**Tape**

Kraft paper with rubber resin adhesive.



**Tape, Label, Write-On™, Bel-Art**  
1" Core, 15 yd (540") roll

1/2" Wide, each ..... \$5.80

Green. **134623005**      Yellow. **134622005**

3/4" Wide, Pk/4 ..... \$25.00

Green. **134820075**      Yellow. **134810075**

White. **134800075**

3" Core, 40 yd roll

1/2" Wide, each ..... \$10.55

Green. **134633005**      White. **134630005**


3/4" Wide, each ..... \$14.00

Green. **134633075**      Yellow. **134632075**

3/4" Wide, White, Pk/4. **134830075** ..... \$56.10

3/4" Wide, Mixed, Pk/6. **134630600** ..... \$85.35


blue, green, orange, red, white & yellow



**Label Tape Dispenser, multi-roll. 134730000..... \$44.40**



Versatile multi-format block holds (96) 0.2 mL tubes/strips, (39) 0.5 mL tubes, or (1) nonskirted 96-well PCR plate. In-situ adapter sold separately. Touchscreen control and height adjustable heated lid. 2-year warranty.




0.2 mL Tubes



96-Well Plate



0.5 mL Tubes



In-Situ Adapter

**Thermal Cycler, multi-format block, Benchmark. T5000-96 .. \$4,471.00**  
In-Situ Adapter, holds 4 slides. **T5000-A-IS** ..... \$439.00