



**Thermometers, Glass, student grade**

Model	Range	Divisions	Length	Features	Price
97-3104	-20 to 230°F/ -30 to 110°C	2°F/1°C	5 3/4"	Red Alcohol/Aluminum V-Back.	\$2.20
97-3110	-10 to 102°C	1°C	6"	Red Alcohol/White Background. Partial immersion.	\$3.25
97-3108	0-230°F/ -20 to 110°C	2°F/1°C	12"	Red Alcohol/White Background. Partial immersion.	\$3.25

**Trays**



**Trays**

- Cafeteria Style, polypropylene, 12 x 16 1/4", NSF, **USA**
- Orange. **105CT1216O** ..... \$5.95
  - Red. **105CT1216R** ..... \$5.95
  - Black. **105CT1216B** ..... \$5.95



Contain spills or transport lab products with this rimmed, durable, large white tray. Holds 5 gallons (18.9 L). *Not autoclavable.*

- HDPE (high-density polyethylene), **USA. 1053506** ..... \$33.95  
 Dimensions: 26 x 18 x 3 1/2" (66 x 45.7 x 8.9 cm)



- Polypropylene, autoclavable
- 12 x 7 3/4 x 5 1/8", Bel-Art, each. **162600000** ..... \$32.35
  - 15 x 14 x 5", Pk/6. **60081721** ..... \$133.00
  - 20 x 17 x 5", Pk/6. **60081722** ..... \$140.00

**Tubes & Vials**



**Tubes, Centrifuge**

- Conical bottom, sterile, white writing patch  
RNase-/DNase-/Pyrogen-free
- 15 mL, natural, polypropylene (PP) tube/HDPE cap  
Sterile Bag/25, Cs/500. **BP1502** ..... \$79.20  
Maximum RCF is 12,000 x g
- Sterile Rack/25, Cs/500. **BP1501** ..... \$89.75  
Maximum RCF is 12,000 x g
- 50 mL, polypropylene (PP) tube/HDPE cap  
Natural, 20 sterile bags/25, Cs/500. **BP5002** ..... \$102.00  
Maximum RCF is 12,000 x g
- Natural, 20 sterile bags/25, Cs/500. **6297** ..... ~~\$176.00~~ \$176.00  
Maximum RCF is 20,000 x g
- Black, 10 sterile bags/50, Cs/500. **HS4429** ..... \$218.00  
Maximum RCF is 9,500 x g
- Self-Standing (skirted)
- 30 mL  
Nonsterile, PP tube/PP cap, Cs/500. **109220** ..... \$128.00  
Dimensions: 23(OD) x 90 mm
- Sterile, PP tube/HDPE cap, Cs/500. **C572-1** ..... \$138.00  
Molded-in graduations from 10 to 30 mL in 5 mL increments. Dimensions: 25.3(OD) x 111 mm

*Continued*