

Simport plastic culture tubes with dual-position caps.



Tubes, Culture,


- 10 x 75 mm (3 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1503** \$47.60
- 12 x 75 mm (5 mL)
Polystyrene, Sterile, dual-position caps, Simport
Cs/500, 20 sterile bags/25. **T415-2** \$73.10
Cs/1000, 8 sterile bags/125. **T415-3** \$104.00
- Polystyrene, Nonsterile, Cs/1000. **110410** \$37.70
- Polypropylene, Sterile, dual-position caps, Simport
Cs/500, 20 sterile bags/25. **T415-2A** \$73.10
- 13 x 100 mm (7 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1510** \$71.60
- 17 x 95 mm (14 mL)
Polypropylene, Sterile, dual-position caps, Simport
Cs/500, 20 sterile bags/25. **T416-2A** \$104.00
- 16 x 125 mm (12 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1515** \$122.00
- 16 x 150 mm (15 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1517** \$116.00



RNase/DNase-free

Tubes, Microcentrifuge, Economy, polypropylene

- 0.5 mL, 0.1 mL graduations, natural, Pk/1000. **HS4422** \$35.00
- 1.5 mL, 0.25 mL graduations, Heathrow
Natural, max RCF is 16,400 x g
Pk/500. **HS4323** \$17.00
Pk/1000. **HS4325** \$33.00
- Black, max RCF is 18,000 x g, Pk/500. **HS4323K** \$19.00
- 5 mL, natural, max RCF is 25,000 x g, Benchmark
Nonsterile, 4 bags/50, Pk/200. **C1005T5** \$57.25
Sterile, 10 bags/20, Pk/200. **C1005T5ST** \$63.65



Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable. Free of RNase/DNase.

Tubes, Microcentrifuge, Premium, PP, G-Tube™, Bio Plas, USA

Amber color protects light sensitive specimens.

- 0.6 mL, Pk/500 \$22.50
Amber. **4048** Lavender. **4047** Natural. **4040**
Mixed (blue, green, orange, pink, & yellow). **4049** \$27.05
- 1.5 mL, 0.5 mL graduations, frosted area on side, Pk/500 \$22.50
Amber. **4038** Blue. **4033** Lavender. **4037** Natural. **4030**
Mixed (blue, green, orange, pink, & yellow). **4039** \$27.05
- 2.0 mL, 0.5 mL graduations, frosted area on side, Pk/500 \$24.15
Amber. **4058** Blue. **4053** Natural. **4050**
Mixed (blue, green, orange, pink, & yellow). **4059** \$28.25

Tubes, Microcentrifuge, Siliconized, G-Tube™, Bio Plas, USA

- 1.5 mL, natural color, Pk/500. **4165SL** \$45.10
Siliconized products should NOT be autoclaved



Tubes, PCR, 0.2 mL thin-wall, polypropylene

- Individual, flat caps, Simport
Blue, Pk/1000. **T325-2VB** \$46.20
Natural, Pk/1000. **T325-2VN** \$46.20
- Strips, Amplitube™, Low Profile, individually attached caps
8-Tube, Pk/125 (1000 tubes). **T320-2LPN** \$94.25
Holds 0.1 mL. Tubes have flat caps & natural color.
- Strips, Bio Plas, *Caps are sold separately; available in natural only.*
8-Tube, Pk/125 (1000 tubes) \$69.10 ↓
Blue. **5010-4** Natural. **5010-1**
8-Tube Cap Strips, domed, Pk/300 (2400 caps). **5030-1**.. \$75.95
12-Tube, Pk/100 (1200 tubes) \$81.35
Blue. **5020-4** Natural. **5020-1**
12-Tube Cap Strips, domed, Pk/200 (2400 caps). **5040-1** \$75.95



Vials, Cryogenic, round bottom, internal thread

- 4 mL, Cs/1000. **T301-4** \$273.00



Vial Decapping Pliers, vinyl grip, USA

- For 13 mm Crimp Seals. **D-13** \$72.00
- For 20 mm Crimp Seals. **D-20** \$72.00
- For 30 mm Crimp Seals. **D-30** \$72.00