


## Illuminators

*Please call regarding availability of items not listed*



Fiber optic light source with dual gooseneck and two focusing lenses is shown.

- Fiber Optics**, Fiber-Lite™, Dolan-Jenner
- Light Source, 150W EKE lamp, rheostat. **FL150/115** ..... \$392.00  
(MI-150)
  - Gooseneck, Dual, chrome plated, 23" long. **FL150/05**..... \$230.00  
(EEG2823M)
  - Focusing Lens, each. **FL150/50** ..... \$39.00  
(LH-759)

## Accessories & Parts

Cushion and support your elbows/forearms with a pair of these ergonomic armrests. Pliable, latex bag stays firmly in position on the workbench, where it conforms to surface and your arm.



- Armrests**, latex rubber bag, black, **USA**
- ★ Low-Profile, 1.25 lbs each, Pair. **14146**..... \$25.00  
*Recommended for newer 'scopes with low controls*
  - High-Profile, 2 lbs each, Pair. **14141** ..... \$27.00




**Boxes**, Plastic, hinged polystyrene, **USA**

Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
<b>H-1</b>	1 x 1 x ¼"	-	Reticles, Pointers	\$0.40
<b>H-1/11B</b>	1 x 1 x ¼"	1	Reticles, Pointers	\$0.50
<b>H-54</b>	2 9/16 x 2 9/16 x ½"	-	Filters	\$1.20
<b>H-54/2916G</b>	2 9/16 x 2 9/16 x ½"	1	Filters	\$1.35
<b>H-28/334B2</b>	3 ¾ x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40

## Counting Chambers

Bright-Line®  
Hemacytometer



- Hemacytometer**, bright-lined, V-slash, Germany  
Chamber with 2 Coverslips, Marienfeld. **0650030** ..... \$134.00
- Hemacytometers**, Bright-Line®, molded borosilicate glass, **USA**  
Standard Thickness (5.5 mm thick)
- Chamber with 2 Coverslips, Hauser. **1492** ..... \$286.00
  - Chamber only, Hauser. **1490** ..... \$258.00
  - Coverslips, 20 x 26 x 0.4 mm, Erie, Pk/2. **2026H** ..... \$11.05
- Reduced Thickness/phase counting chamber (3.5 mm thick)
- Chamber with 2 Coverslips, Hauser. **1475** ..... \$305.00
  - Coverslips, 20 x 26 x 0.4 mm, Erie, Pk/2. **2026H** ..... \$11.05
- Howard Mold Chamber**, Hauser, **USA**  
*Requires Howard Disc - see RETICLES & STAGE MICROMETERS*
- Chamber with 2 Coverslips (5080 & 5090). **3820** ..... \$383.00
  - Coverslip only, 28 x 33 x 0.5 mm (thin), each. **5080** ..... \$46.80
- Petroff-Hausser Bacteria Counter**, Hauser, **USA**  
Coverslip only, 20 x 26 x 0.2 mm, each. **5051** ..... \$46.80

- England Field Finder**, 1 mm thick, **USA**. **KR-834** ..... \$212.00

- Eyepiece Pointer**, wire. **14301** ..... \$2.00




- Eyepiece Sleeve**, 23.2 mm increased to 30 mm. **14320** ..... \$35.00

# SLIDE & CASSETTE PROCESSING

## Microscope Slides

### Slides, general, 75 x 25 x 1 mm

	bx/72	cs/20 bx
Plain, ground edges		
Premier®, safety corners. <b>9101</b> ..... 	\$3.90	76.00
Globe, 90° corners, "Globe". <b>1301</b> .....	—	110.00
Frosted, beveled edges		
Premier®, sandblasted, safety corners. <b>8205</b> ....	\$6.05	115.00
Globe, etched, clipped corners, "Globe". <b>1334</b> ...	—	137.00
White End, clipped corners		
Premier®. <b>9108W-E</b> .....	—	115.00
Superfrost®, "SPECIMEN". <b>4951-001T</b> .....	—	266.00





### Slides, specialty


Charged (silanated), 75 x 25 x 1 mm		
Recommended for adhering tissue sections and cell preparations.		
Diamond™, white end, Globe		
Cs/20 Bx. <b>1358W</b> .....		\$302.00
Superfrost® Plus, white end, Erie (Epredia)		
Cs/20 Bx. <b>4951PLUS-001</b> .....		\$617.00
Preferred slide for Ventana instrumentation. (12-550-15, 48311-703)		
Concavity		
75 x 25 x 1.3 mm, 16 mm dia., 0.8 mm deep, economy		
1-Well, Bx/12. <b>600CS3X11</b> .....		\$6.55
2-Well, Bx/12. <b>600CS3X12</b> .....		\$9.00
76 x 26 x 3.2 mm, 18 mm dia., 1.75 mm deep, Erie (Epredia)		
1-Well, Bx/12. <b>1526</b> .....		\$97.35
2-Well, Bx/12. <b>1527</b> .....		\$109.00
Fused Quartz, NSG Type OZ®, 3" x 1" x 1 mm, <b>USA</b>		
<b>QMS3X1X1</b> , each.....		\$20.00
Fused Silica, Corning 7980, 3" x 1" x 1 mm, <b>USA</b>		
<b>FSMS3X1X1</b> , each.....		\$29.00


## Coverslips (Cover Glasses)


### Coverslips, premium, borosilicate glass, 1 oz/pk, Cs/10 oz

#0 Thickness (0.08-0.12 mm)		
Circles, Germany* 		
25 mm, ~222 pcs/oz, Cs/10 oz. <b>0-25</b> .....		\$214.00

#1 Thickness (0.13-0.17 mm)		
Squares & Rectangles, Erie (Epredia) 	1 oz	cs/10 oz
22 x 32 mm, ~162 pcs/oz. <b>1-2222</b> .....	—	\$64.95
22 x 30 mm, ~119 pcs/oz. <b>1-2230</b> .....	—	\$58.80
22 x 50 mm, ~71 pcs/oz. <b>1-2250</b> .....	—	\$58.80
24 x 50 mm, ~65 pcs/oz. <b>1-2450</b> .....	—	\$58.80
24 x 60 mm, ~54 pcs/oz. <b>1-2460</b> .....	—	\$58.80

#1 Thickness (0.13-0.17 mm)		
Circles, Germany* 	cs/5 oz	cs/10 oz
12 mm, ~663 pcs/oz. <b>1-12</b> .....		\$268.00
18 mm, ~296 pcs/oz. <b>1-18</b> .....		\$150.00

#1.5 Thickness (0.16-0.19 mm). Recommended for most users		
Squares & Rectangles, Erie (Epredia) 	1 oz	cs/10 oz
18 x 18 mm, ~209 pcs/oz. <b>1.5-1818</b> .....	—	\$62.95
22 x 22 mm, ~140 pcs/oz. <b>1.5-2222</b> .....	\$7.70	\$64.40
22 x 40 mm, ~77 pcs/oz. <b>1.5-2240</b> .....	—	\$58.80
22 x 50 mm, ~62 pcs/oz. <b>1.5-2250</b> .....	—	\$58.80
24 x 30 mm, ~94 pcs/oz. <b>1.5-2430</b> .....	—	\$58.80
24 x 40 mm, ~71 pcs/oz. <b>1.5-2440</b> .....	—	\$58.80
24 x 50 mm, ~56 pcs/oz. <b>1.5-2450</b> .....	\$7.00	\$58.80
24 x 60 mm, ~47 pcs/oz. <b>1.5-2460</b> .....	—	\$58.80

#1.5 Thickness (0.16-0.19 mm). Recommended for most users		
Circles, Germany* 	cs/5 oz	cs/10 oz
12 mm, ~552 pcs/oz. <b>1.5-12</b> .....	\$223.00	—
18 mm, ~245 pcs/oz. <b>1.5-18</b> .....	\$250.00	—
22 mm, ~164 pcs/oz. <b>1.5-22</b> .....	—	\$252.00
25 mm, ~127 pcs/oz. <b>1.5-25</b> .....	—	\$252.00

\*German glass is preferred for cell growth/adhesion

### Coverslips, specialty

Counting Chamber - see Counting Chamber Coverslips in MICROSCOPE Accessories & Parts section

#### Fused Quartz, **USA**

0.15-0.25 mm Thick, NSG Type OZ®

1" Square, each. <b>QCS1X1</b> .....	\$19.00
1" Circle, each. <b>QCS1X0</b> .....	\$19.00

0.16 ± 0.01 mm Thick, A1

25 mm Square, each. <b>R525025</b> .....	\$32.40
25 mm Circle, each. <b>R525000</b> .....	\$23.40

#### Fused Silica, **USA**

0.15-0.25 mm Thick, Corning 7980

1" Square, each. <b>FSCS1X1</b> .....	\$21.50
1" Circle, each. <b>FSCS1X0</b> .....	\$21.50

0.16 ± 0.01 mm Thick, S1-UV

25 mm Square, each. <b>R425025</b> .....	\$32.00
25 mm Circle, each. <b>R425000</b> .....	\$28.80