


Tube Storage



Baskets, test tube, aluminum, Heathrow

- Small, 15.2 x 15.2 x 15.2 cm(H). **HS20341A**..... \$27.00
 Medium, 25.6 x 15.2 x 15.2 cm(H). **HS20341B**..... \$30.00



Each cube holds (36) 1.5/2.0 mL tubes, or (96) 0.2 mL tubes, or (12) 0.2 mL strips, or (1) 96-well PCR plate and (10) 1.5/2.0 mL tubes. Maintains 0°C for 4 hours.

Tube Cooler, benchtop, 0°C, CoolCube™, Benchmark

- Microtube and PCR Plate. **R1000** \$121.00



SB2CC-64 Cryogenic
Tube Storage Box

Tube Storage Boxes, Cardboard, vented

- 2" Box Only, 5¼ x 5¼ x 2", each. **SB2CC**..... \$1.65
 2" Box w/dividers, white
 100-Cell (10 x 10), Pk/12. **SB2CC-100** \$26.95
 ★ 81-Cell (9 x 9), Pk/12. **SB2CC-81** \$26.95
 64-Cell (8 x 8), Pk/12. **SB2CC-64** \$26.95
 2" Mini box w/dividers, 3 x 3 x 2"
 25-Cell (5 x 5), Pk/12. **MB2C-25** \$29.55
 Holds 13 mm vials and 4.0-5.0 mL cryovials.
 3" Box Only, 5¼ x 5¼ x 3", each. **SB3C** \$1.75
 3" Box w/dividers
 100-Cell (10 x 10), Pk/12. **SB3C-100** \$28.00
 81-Cell (9 x 9), Pk/12. **SB3C-81** \$28.00
 64-Cell (8 x 8), Pk/12. **SB3C-64** \$28.00
 Dividers only, Fit 2 & 3" Boxes
 100-Cell (10 x 10), Pk/12. **CD-100** \$11.55
 Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.
 81-Cell (9 x 9), Pk/12. **CD-81** \$11.55
 Holds 13 mm vials and 4.0-5.0 mL cryovials.
 64-Cell (8 x 8), Pk/12. **CD-64** \$11.55
 Holds 14 mm vials.
 49-Cell (7 x 7), Pk/12. **CD-49** \$11.55
 Holds 14 mm vials.
 25-Cell (5 x 5), 25 mm cell, Pk/12. **CD-25** \$11.55
 16-Cell (4 x 4), 31 mm cell, Pk/12. **CD-16** \$11.55

- 5" Box Only, 5¼ x 5¼ x 5", each. **SB4C** \$3.20
 5" Box w/dividers
 64-Cell (8 x 8), Pk/12. **SB4C-64** \$46.35
 36-Cell (6 x 6), 15 mL tubes, Pk/12. **SB4C-36** \$46.35
 16-Cell (4 x 4), 50 mL tubes, Pk/12. **SB4C-16** \$46.35
 Dividers only, Fit 5" Box
 64-Cell (8 x 8), Pk/12. **CD4-64** \$11.55
 36-Cell (6 x 6), Pk/12. **CD4-36** \$11.55
 16-Cell (4 x 4), Pk/12. **CD4-16** \$11.55



Holds 1.0-2.0 mL tubes and cryovials.
 Features a locking hinged lid, and base with molded alphanumeric coordinates.
 Dimensions: 12.8 x 12.8 x 5.4 cm.



Tube Storage Boxes, Polypropylene

- 81-Well (9x9), 2 blue/2 yellow boxes, Pk/4. **60066501** \$27.00



Holds 1.5-2.0 mL tubes
in a two-piece rack.

Tube Storage Boxes, Polystyrene Foam (styrofoam), USA

- 50-Well, Pk/10. **HS4283P** \$42.85
 100-Well, Pk/10. **HS4282** \$53.60




Tube Racks, Grip, polypropylene, No-Wire™, Bel-Art


- 15-16 mm Tubes, 60-place, each. **187490001** \$27.35
 18-20 mm Tubes, 40-place, each. **187490002** \$27.35

Continued

Continued



- Tube Racks, Interlocking**, polypropylene, Heathrow
- Cube, mixed colors, Pk/5. **HS29050A** \$42.70
Holds (4) 50 mL tubes, (10) 15 mL tubes,
(12) 12 x 75 or 12 x 100 tubes, or (16) 1.5-2.0 mL tubes.
 - 4-Way, mixed colors, Pk/5. **HS29022G** \$39.35
Holds (4) 50 mL tubes, (12) 15 mL tubes,
(32) 1.5 mL tubes, or (32) 0.5 mL tubes.
 - 4-Way Rotatable, Rota-Rack™, four connected modules
Large Tubes, 22.4 x 11 x 11 cm, each. **HS29040A** \$41.50
Each module holds (2) 50 mL tubes, (8) 15 mL tubes,
(6) 20 mm tubes, or (10) 10-13 mm tubes. Mixed colors.



- Tube Racks, Mega Rack**, polypropylene, Heathrow
- 10-13 mm Tubes
 - 216-Place, 33.7 x 26.7 x 9.1 cm(H)
 - Blue, each. **HS24322A** \$38.00
 - White, each. **HS24322E** \$38.00
 - 432-Place, 67.4 x 26.7 x 9.1 cm(H)
 - Blue, each. **HS120113** \$70.00
 - White, each. **HS120112** \$70.00




- Tube Racks, Microtube**, polypropylene, Heathrow
- 80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm
Holds 1.5-2.0 mL tubes in conical wells.
 - Blue, Pk/5. **HS29025B** \$25.35
 - Pink, Pk/5. **HS29025D** \$25.35
 - Mixed, Pk/5. **HS29025G** \$25.35
 - 96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.
 - Blue, Pk/5. **HS2345C** \$49.25
 - Green, Pk/5. **HS2345D** \$49.25
 - Pink, Pk/5. **HS2345E** \$49.25


Available in blue,
green, orange,
purple & yellow.



- Tube Rack, PCR**, polypropylene, mixed colors, Heathrow
Reversible, 96-Well, Pk/5. **HS2344A** \$53.30
Holds 0.2, 0.5 and 1.5 mL tubes.



- Tube Racks, Peg**, 14-17 mm Tubes, 50-place, polypropylene
- Blue, Pk/2. **HS24311B** \$22.75
 - White, Pk/2. **HS24311A** \$22.75



Tilt your test tube rack 25°, 35°, or 45°. Create agar slant tubes, or adjust the angle of tube racks when pipetting to reduce hand strain.

- Tube Rack Slanter**, Bel-Art. **188900100** \$22.05
Platform size: 19.6 cm(W) x 12.7 cm(D)




- Tube Racks, Teardrop, OneRack™**, POM, Heathrow
- Half-Size
 - 16 mm Tubes, 36-place, 12.7 x 12.7 x 7 cm(H)
 - Magenta, Pk/8. **HS27502D** \$100.00
 - Full-Size
 - 13 mm Tubes, 72-place, 20.2 x 10.4 x 5.9 cm(H)
 - Blue, Pk/8. **HS27511B** \$130.00
 - Magenta, Pk/8. **HS27511D** \$130.00
 - 16 mm Tubes, 72-place, 25 x 12.7 x 7 cm(H)
 - Blue, Pk/8. **HS27512B** \$130.00
 - Magenta, Pk/8. **HS27512D** \$130.00




Polypropylene cube and lid has two snap-in, gridded inserts to hold either (16) 15 mL tubes or (9) 50 mL tubes.

- Tube Racks, Transformer Cube**, Pk/5. **HS23271A** \$53.85

Continued



Tube Racks, Wire, epoxy-coated, blue
 10-13 mm Tubes, 108-place (9 x 12), each. **345120759** \$32.00
 13-16 mm Tubes
 108-Place (9 x 12), each. **345120763** \$32.00
 180-Place (9 x 20), each. **187860770** \$51.45
 Dimensions: 38.6 x 18.3 x 6.4 cm(H)



17515ML-25

Tube Racks, 15 mL
 25-Place (5 x 5), stainless steel, **USA**, ea. **17515ML-25** \$14.80
 80-Place (5 x 16), stainless steel, **USA**, ea. **17515ML-80** \$33.40
 Dimensions: 40.7 x 12.7 x 10.2 cm(H)



17550ML-9

9-Place (3 x 3), stainless steel, **USA**, ea. **17550ML-9** \$14.80
 16-Place (4 x 4), epoxy-coated wire, ea. **187940001** \$50.55

Weighing Dishes



Hexagonal



Standard



Pour

Weigh Boats, polystyrene, Heathrow

Standard

- ★ Small, 4.6 x 4.6 x 0.8 cm(D), Bx/500. **HS1420A** \$21.45
- ★ Medium, 8.5 x 8.5 x 2.4 cm(D), Bx/500. **HS1420B** \$35.05
- ★ Large, 14 x 14 x 2.2 cm(D), Bx/500. **HS1420C** \$57.25

Hexagonal, anti-static

- Small, 3.8 cm top dia., Bx/500. **HS14251A** \$15.40
- Medium, 6.4 cm top dia., Bx/500. **HS14251B** \$27.35
- Large, 12 x 14 cm top OD, Bx/500. **HS14251C** \$71.25
(25433-104)
- X-Large, 14 x 16 cm top OD, Bx/500. **HS14251D** \$90.95
(25433-106)

Pour (Canoes)

- Small, 5.8 x 4.3 x 1.3 cm(D), Bx/500. **HS1419A** \$35.75
- Medium, 13.2 x 8.3 x 2.6 cm(D), Bx/500. **HS1419B** \$70.45
- Large, 18.3 x 10.8 x 2.6 cm(D), Bx/500. **HS1419C** \$101.00

Wipes



Extra low-lint, nonabrasive, highly absorbent, 1-ply wipes are perfect for microscope maintenance.

Kimwipes® bx/280 cs/60 bx
 4½ x 8½" Wipes. **34155** \$3.55 190.00



WypAll® L40

Lab Wipes, ¼ fold, low lint pk/56 cs/18 pk
 WypAll® L40, limited use. **43005701** \$6.00 96.00
 Measures 12 x 12.5"
 WypAll® X60, pk/100, cs/9 pk, Lot/9 cs. **43034770** \$2,610.00
 Measures 11 x 23"