


- Tube Racks, Interlocking**, polypropylene, Heathrow
- Cube, mixed colors, Pk/5. **HS29050A** ..... \$42.70  
Holds (4) 50 mL tubes, (10) 15 mL tubes,  
(12) 12 x 75 or 12 x 100 tubes, or (16) 1.5-2.0 mL tubes.
  - 4-Way, mixed colors, Pk/5. **HS29022G** ..... \$39.35  
Holds (4) 50 mL tubes, (12) 15 mL tubes,  
(32) 1.5 mL tubes, or (32) 0.5 mL tubes.
  - 4-Way Rotatable, Rota-Rack™, four connected modules  
Large Tubes, 22.4 x 11 x 11 cm, each. **HS29040A** ..... \$41.50  
Each module holds (2) 50 mL tubes, (8) 15 mL tubes,  
(6) 20 mm tubes, or (10) 10-13 mm tubes. Mixed colors.



- Tube Racks, Mega Rack**, polypropylene, Heathrow
- 10-13 mm Tubes
  - 216-Place, 33.7 x 26.7 x 9.1 cm(H)
  - Blue, each. **HS24322A** ..... \$38.00
  - White, each. **HS24322E** ..... \$38.00
  - 432-Place, 67.4 x 26.7 x 9.1 cm(H)
  - Blue, each. **HS120113** ..... \$70.00
  - White, each. **HS120112** ..... \$70.00




- Tube Racks, Microtube**, polypropylene, Heathrow
- 80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm  
Holds 1.5-2.0 mL tubes in conical wells.
  - Blue, Pk/5. **HS29025B** ..... \$25.35
  - Pink, Pk/5. **HS29025D** ..... \$25.35
  - Mixed, Pk/5. **HS29025G** ..... \$25.35
  - 96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm  
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.
  - Blue, Pk/5. **HS2345C** ..... \$49.25
  - Green, Pk/5. **HS2345D** ..... \$49.25
  - Pink, Pk/5. **HS2345E** ..... \$49.25

Available in blue,  
green, orange,  
purple & yellow.



- Tube Rack, PCR**, polypropylene, mixed colors, Heathrow  
Reversible, 96-Well, Pk/5. **HS2344A** ..... \$53.30  
Holds 0.2, 0.5 and 1.5 mL tubes.



- Tube Racks, Peg**, 14-17 mm Tubes, 50-place, polypropylene
- Blue, Pk/2. **HS24311B** ..... \$22.75
  - White, Pk/2. **HS24311A** ..... \$22.75



Tilt your test tube rack 25°, 35°, or 45°. Create agar slant tubes, or adjust the angle of tube racks when pipetting to reduce hand strain.

- Tube Rack Slanter**, Bel-Art. **188900100** ..... \$22.05  
Platform size: 19.6 cm(W) x 12.7 cm(D)



- Tube Racks, Teardrop, OneRack™**, POM, Heathrow
- Half-Size
  - 16 mm Tubes, 36-place, 12.7 x 12.7 x 7 cm(H)
  - Magenta, Pk/8. **HS27502D** ..... \$100.00
  - Full-Size
  - 13 mm Tubes, 72-place, 20.2 x 10.4 x 5.9 cm(H)
  - Blue, Pk/8. **HS27511B** ..... \$130.00
  - Magenta, Pk/8. **HS27511D** ..... \$130.00
  - 16 mm Tubes, 72-place, 25 x 12.7 x 7 cm(H)
  - Blue, Pk/8. **HS27512B** ..... \$130.00
  - Magenta, Pk/8. **HS27512D** ..... \$130.00



Polypropylene cube and lid has two snap-in, gridded inserts to hold either (16) 15 mL tubes or (9) 50 mL tubes.

- Tube Racks, Transformer Cube**, Pk/5. **HS23271A** ..... \$53.85  
*Continued*