

MAINTENANCE & REPLACEMENT PARTS

Bulbs (Lamps)



Please call regarding availability of bulbs not listed

ANSI codes and manufacturer's part numbers are cross-referenced (in italics and parenthesis) when available.

Tungsten Halogen & Incandescent

* Recommend Bulb Removal/Installation Tool (shown after listing)

2.5V/0.35A Bullet. 252	\$4.00
<i>(634, 8000155)</i>	
3.5V/0.72A Bullet, Boehm. 03100	\$10.60
<i>(3287, 8000017)</i>	
6V/0.2A Bullet (for teaching head). 328	\$1.50
<i>(79098, 31-31-50, 8000168)</i>	
6V/0.2A Bullet w/Lens (for multi-TH), Nikon. MXA20101	\$11.00
<i>(79186)</i>	
6V/1.2A Horizontal Filament, Hosobuchi. OP2101	\$15.45
<i>(MA560, 8000369)</i>	
6V/1.2A Vertical Filament, Hosobuchi. OP2101/U	\$15.45
<i>(MA561, 8000370)</i>	
6V/2A SC Bayonet, Frosted, round, Ushio. 8-C201	\$7.10
<i>(8000301)</i>	
6V/3A (6V/18W) Frosted, ½ Silver, Swift. MA737	\$12.00
<i>(8000350)</i>	
6V/4.5A DC Prefocus*, Ushio. 1468	\$5.00
<i>(581, 8000067)</i>	
6.4V/2.63A DC Bayonet. S8.1130	\$2.30
6.5V/2.75A (6.5V/15W) DC Bayonet. 1493	\$2.75
<i>(360, 700-034, 8-B102, MA402, 21138, 8000070)</i>	
6.5V/2.75A DC Prefocus*, Ushio. 1460X	\$5.50
<i>(31-31-75, 31-31-76, 31-01-05, 8000066)</i>	
6.5V/2.75A DC Prefocus*. GE1630	\$9.30
<i>(8000092)</i>	
6.5V/2.75A DC Prefocus*. 1631X	\$8.00
<i>(31-31-85, 31-31-86, 31-31-79, 8000093)</i>	
6.5V/2.75A DC Prefocus*. GE1631X	\$13.50
<i>(31-31-85, 31-31-86, 31-31-79, 8000093)</i>	
6.5V/2.75A DC Prefocus*, Ushio. 1649	\$9.90
<i>(31-31-42, 31-31-43, 8000100)</i>	
6V/5W Screw, Ushio. 8000302	\$10.00
<i>(8-G101)</i>	
6V/9W MR11, Philips. 13529	\$21.50
6V/10W (6V/1.67A) Halogen, G4. GE787	\$7.50
<i>(31-31-77-01, 31-31-77-10, MA570, 76342, 8000254)</i>	
6V/10W Halogen, G4, Osram. 64223	\$17.00
<i>(38-61-08, 91-01-54, 5972, 1000857)</i>	
6V/10W Halogen, G4, Osram. 64225	\$4.50
<i>(ESA/FHD, 76342, 050-550, 050-556, 700-065, 8-C403, 907-090, 624, 655C, 7387, 3369, 1000531)</i>	

6V/10W Halogen, G4, Philips. 5972	\$7.50
<i>(38-61-08, 91-01-54, 64223, 1000857)</i>	
6V/10W Halogen, G4, Philips. 7387	\$5.20
<i>(ESA/FHD, 76342, 050-550, 050-556, 700-065, 8-C403, 907-090, 624, 655C, 7387, 3369, 1000531)</i>	
6V/10W SC Bayonet, round. 5001	\$1.83
<i>(3040-0260, 6814P)</i>	
6V/10W Screw, Hosobuchi. OP2103	\$9.75
<i>(8-G102, 8000303)</i>	
6V/12W Halogen, G4. GE786	\$8.20
<i>(8000252)</i>	
6V/15W Bulb and Socket, Ushio. 8-C103	\$22.50
6V/15W MR11, Philips. 13528	\$22.00
<i>(8-C408, 630)</i>	
6V/15W DC Bayonet, Flanged, Oblong. 77903	\$10.00
<i>(OP2130, 8000240)</i>	
6V/15W DC Bayonet, Flanged, Nikon. 77910	\$13.50
<i>(8000242)</i>	
6V/15W DC Bayonet, Oblong, Ushio. 8018	\$7.40
<i>(3800-29-9020, 76906, 166-324, 500-012, 8000269)</i>	
6V/15W DC Flanged. 177Z	\$9.25
<i>(3800-18-1740, 38-01-77, OP2131, 8000172)</i>	
6V/15W (6V/3A) DC Bayonet, Round. OP2114	\$6.25
<i>(76907, 8000233)</i>	
6V/20W Halogen, G4. GE778	\$8.90
<i>(31-31-81-01, 31-31-81-10, 569, 31-31-21, 8000239)</i>	
6V/20W (6V/3.33A) Halogen, G4. GE788	\$7.90
<i>(MA260/05, 1121)</i>	
6V/20W Halogen, G4, Osram. 64250 HLX	\$4.40
<i>(ESB/FHE, 1120, 362-658, 3801-20-1740D, 11143, 500-245, 89500, 79099, 8-C405, MA780, 7388, M420, MA07, 687.95, 3256, 12-561-BLB1, MA311, 1000532)</i>	
6V/20W Halogen with Lead, Narva. 55133	\$33.00
6V/20W MR11, Ushio. JCR/M6V20W	\$18.40
<i>(83314)</i>	
6V/25W Flanged, G4. 39-01-53	\$19.00
<i>(8000174)</i>	
6V/25W Halogen with Lead, Narva. 55147	\$33.00
<i>(33-01-53)</i>	
6V/30W (6V/5A) Bayonet, Oblong, Hosobuchi. OP2118	\$6.70
<i>(MA738, 77904, 8-C101, 8000241)</i>	
6V/30W Bulb and Socket, Ushio. 8-C102	\$22.50
<i>(8000299)</i>	
6V/30W (6V/5A) Candelabra, Osram. 8100	\$18.75
<i>(500-011, 8000277)</i>	
6V/30W DC Flanged, Ushio. 8000175	\$20.00
<i>(39-01-58)</i>	
6V/30W Halogen, G4, Philips. 5761	\$7.40
<i>(84246, 8-C410, MA326, 3256-30W, LB-30WH, SE06.0215, MINJ-F-B2, 1000867)</i>	
6V/30W Halogen, G4, Ushio. 1000867	\$11.90
<i>(84246, 8-C410, MA326, 3256-30W, LB-30WH, SE06.0215, 5761, MINJ-F-B2)</i>	
6V/30W Halogen with Lead, Narva. 55168	\$39.00
6V/35W Halogen, G4, vertical filament, Osram. 64275	\$7.50
8V/20W MR11, Osram. 64255	\$18.20
<i>(Y99702)</i>	

Continued

10V/70W Horizontal Filament, Nikon. 71256	\$16.00	21V/150W MR16, 40 hrs, Ushio. EJA	\$11.50
(8000226)		(A08110, 76223, 8-C404, 902-063, N650013, N650012, 1000297)	
10.8V/30W MR16, Ushio. EKZ	\$16.50	21V/150W MR16, 40 hrs. EJV	\$13.20
(1000315)		(1000302)	
12V/10W DC Bayonet. GE1004	\$2.45	21V/150W MR16, 200 Hrs, Ushio. EKE	\$11.00
(800-411)		(902-058, 31-31-52, A08120, FL180/70, FL150/70, FT150/70, FT190/70, 3378, 1000306)	
12V/10W SC Bayonet, round, horizontal filament. S-5203Z ...	\$33.30	21V/150W MR16, 1000 Hrs, Osram. EKE/X	\$13.80
(800-421)		(EKE/L, 1003370)	
12V/10W Halogen, G4. 64415	\$3.40	24V/150W Halogen, High Intensity, Osram. 64640 HLX	\$3.40
12V/10W Halogen, G4. 800-170	\$13.50	(FCS, 71826, 799-892, MA786, 7158, 037-717, 799-892, 58101, 501-109, 1000492)	
12V/10W Tubular. 800-400	\$9.90	24V/150W Halogen, Long Life, Osram. 64642 HLX	\$8.00
12V/15W Halogen, G4. 800-422	\$15.30	(FDV)	
12V/15W MR11C, faceted. JCR12V15W/G4	\$5.90	24V/250W Halogen, G6.35, Osram. 64657	\$5.30
12V/20W Halogen, G4, 2000 hrs, Osram. 64425	\$2.85	(EVC/FGX, 6958)	
(58661, 800-160, 3253, HT-59-1152, SE06.0098)		24V/250W MR16, Osram. 64653 HLX	\$12.50
12V/22W MR11. JCR12V22WA/3 , KLS.....	\$34.85	(ELC, 13163, Y99704, 1000318)	
(8-S110)		82V/250W MR16. EVW	\$12.75
12V/25W SC Bayonet. 93	\$1.25	82V/300W MR13, Eiko. FHS	\$12.25
12V/30W Halogen, G6.35, Osram. 64261	\$14.50	120V/6W DC Bayonet. 6S6DC120V	\$2.50
(500-317, 115000317)		(10407, 1005, 8000220)	
12V/30W Halogen, PG22, Osram. 64260	\$17.60	120V/15W Candelabra, Ushio. 15S11/13CL	\$5.20
12V/35W MR11C, Philips. 14626	\$13.25	(626, 8000079)	
(FTE-CG)		120V/15W DC Bayonet. 15T7DC	\$3.00
12V/45W SC Bayonet, clear. BA15s	\$15.00	(77906, 12610, 12603, 800-101, 8000083)	
12V/50W Halogen, G6.35, Osram. 64610 HLX	\$4.50	120V/15W Medium Screw, Eiko. 15S11/102CL	\$3.00
(BRL, 76516, 8-C407, MA784, 500-182, 7027, M650, MA275/05, 1000071)		(31-31-15, 611, 8-B150, 8000077)	
12V/50W Halogen, Special Base, Dr Fischer. 384-643	\$59.00	120V/20W ½ Silver, DC Bayonet, Oblong. 77427	\$13.00
(4643 HLX, 90-1300, M-01020)		(MA361/10)	
12V/60W DC Bayonet. 3800-18-2520	\$20.50	120V/20W ½ Silver, DC Bayonet, Oblong, Ushio. 78595	\$11.80
(38-02-16, 500-010, 8000173)		(MA361/10, MXA20163, 8000385)	
12V/75W MR16, Ushio. EFN	\$14.00	120V/20W ½ Silver, DC Bayonet, Round. MA886	\$12.20
(3800-79-9030, 64615 HLX, 6853, 1000270)		120V/20W DC Bayonet, Oblong. 800-102	\$13.50
12V/100W Halogen, High Intensity, Osram. 64625 HLX	\$4.75	120V/20W DC Bayonet, Round, Eiko. MA2201	\$7.00
(FCR, 76513, 500-974, 8-C401, 193-049, 3800-59-1660, MA782, 7023, 9-8600-15, M32-05-45, 2052, 1000490)		(MA2251B, 8000344)	
12V/100W Halogen, Long Life, Osram. 64623 HLX	\$5.70	120V/20W Frosted/Silver, Round. MA736	\$14.00
(EVA, 84125, 8-C406, 043-144, 7724, 3800-79-9540, 1000381)		120V/30W DC Bayonet, Round. MA722A	\$7.00
12V/100W Halogen, 2-pin prefocus, Osram. 64628	\$11.40	120V/30W DC Bayonet, S11. BLC	\$4.80
(FDT, 31-31-58, 5973, 1000504)		(31-31-71-03, 31-31-71-05, 8-B152, 1000060)	
12V/100W MR16, Osram. EPF	\$9.70	120V/30W Miniature Flood. MA263/05	\$20.00
(64627, 098-145, 6834, 3800-79-9040, 1000271)		120V/30W ½ Silver, DC Bayonet, Ushio. 8000165	\$18.00
13.8V/85W MR16, Osram. 64618	\$12.00	(78508, MXA20450, SM-31-74-28)	
(DED, 13288, 13194)		120V/150W Medium Screw, Higuchi. 150P25/2SB	\$20.00
14.5V/90W MR16, Ushio. EPV	\$13.75	120V/150W MR16, faceted, 12 hrs, Ushio. ESD	\$26.55
15V/150W Halogen, G6.35, EVB/BRJ	\$5.25	(906117, 1000356)	
(64633 HLX, 6550, 100069)		120V/250W MR16, faceted, Ushio. ENH	\$17.75
15V/150W MR16, Osram. 64635 HLX	\$42.00	(54986, 1006926, 1000333)	
15V/150W MR16, Ushio. EFR	\$15.50	120V/300W GY9.5, Ushio. JCD120V300WC	\$14.50
(997-053, 8-S401, 31-01-98, 64634 HLX, 6423, 58102, 10001-111, Y99701, 1000272)		(1000896)	
20V/1A DC Prefocus*. 1634	\$9.00	120V/300W MR16, stippled, oval 2-pin, Osram. ELH	\$11.75
(31-31-37, 050-540, 506-014, 8000095)		(64660, 1000321)	
20V/150W MR16, faceted, 500 hrs, Ushio. DDL	\$11.00	125V/300W MR16, Long Life, faceted, Eiko. ELH-5	\$17.00
(A08130, 82931, 3378, 23500-014, 1000173)		(1003020)	
21V/150W MR16, 2-pin vented, Ushio. DNF	\$20.00	130V/50W Halogen, G8, Hiraki. JCD130V50W/G8	\$5.50
(61351, 001160, 243054, SY0023, 4000-0996-0019, 6000-150, 52118, 5212-11, 5212-51, 5212-81, 1000207)		(25088A)	

Continued

LED

3.4V/0.06W LED, Wedge. 800-001	\$13.50
(MA2215)	
3.4V/3W LED Assembly. 800-500	\$35.10
LED Assembly. 800-425	\$35.10
LED Assembly. 800-426	\$36.90
LED Assembly. 800-427	\$45.00
LED Assembly. 800-428	\$13.50

Metal Halide

60V/24W MR11, 300 hrs, Welch Allyn. 09500-U	\$270.00
(AL-1824, 8-61257)	

Mercury Short Arc

50W Mercury Arc, 100 Hrs, Osram. HBO50W/AC/L1	\$115.00
(78564, 78588, 8-B199, 38-16-19, 414140, 2055, 3040-0090, 500-137, SHL-50AC1, BA-001, 5000487)	
50W Mercury Arc, 100 Hrs, Osram. HBO50W/AC/L2	\$115.00
(31-31-65, SHL-50AC2, BA-002, 5000857)	
50W Mercury Arc, 200 Hrs, Osram. HBO50W/3	\$140.00
(BA-003)	
100W Mercury Arc, 200 Hrs, Ushio. USH-102D	\$109.00
(78589, 8-B192U, 500-138, 3800-18-4060/D, 2087-001, 3040-0380, BA-005A, HBO100W/2, M/5000/133, 5000274)	
100W Mercury Arc, 300 Hrs, Ushio. USH-103D	\$123.00
(78590, 8-B1920, 380301-9350, USH-103OL, 5001329)	
100W Mercury Arc, 200 Hrs, Osram. HBO100W/2	\$109.00
(78589, 8-B192U, 500-138, 3800-18-4060/D, M/5000/133 2087-001, 3040-0380, BA-005, 3254, USH-102D)	
100W Mercury Arc, 300 Hrs, Osram. HBO103W/2	\$123.00
(Not compatible with all Hg power supplies)	
(78590, 8-B1920, 380301-9350, 500-321)	
200W Mercury Arc, 1000 hrs, Ushio. USH-200DP	\$135.80
(HBO200W/2, 31-01-01, 649, 860014, 78580, 8-B190, 3800-18-4070, 38-16-15, 050-320, MA885, 91-00-50, 505-620, 78590, 5000276)	

Xenon Short Arc

75W Xenon Arc, 400 hrs, Ushio. UXL-75XE	\$195.00
(XBO75W/2, 8-B198, 910881-2013, 38-16-02, 500-139, 91-01-42, 3800-79-9190, 78682, 620-125, 903-050, BA-013, 5000371)	
75W Xenon Arc, 400 hrs, Osram. XBO75W/2	\$184.00
(UXL-75XE, 8-B198, 910881-2013, 38-16-02, 500-139, 91-01-42, 3800-79-9190, 78682, 620-125, 903-050, BA-013, 5000371)	
150W Xenon Arc, Short Gap. XBO150W/CROFR	\$260.00
(8-B195, BA-018, UXL-150S-0, 3801-22-2450, 500-008, 5000326)	
180W Xenon Arc, Osram. XBOR180W45COFR	\$450.00
(69183, 407.309)	

Fluorescent

4W Fluorescent, Tubular, Eiko. F4T5/D	\$7.00
(31-31-03, MA309/05, 3000108)	
5W Fluorescent, G23, Sylvania. CF5DS/841	\$5.00
(21278, 800-138, MA68, 48920-81, 20303)	
9W Fluorescent, G23, Sylvania. CF9DS/841	\$4.95
(20305)	
22W Fluorescent, 8" OD Circle, Sylvania. FC8T9/D/RS	\$19.00
(29980)	
32W Fluorescent, 12" OD Circle, Sylvania. FC12T9/D/RS	\$15.90
(20030)	

Germicidal

4W Germicidal (UV-C), G5, Osram. G4T5/OF	\$10.35
(23397)	
15W Germicidal (UV-C), G13, Sylvania. G15T8	\$17.00
(21612, 3000007)	
25W Germicidal (UV-C), G13, Ushio. G25T8	\$19.80
(3000008)	
30W Germicidal (UV-C), G13, Sylvania. G30T8	\$20.00
(23112, 00204, 11080, 3000009)	
39W Germicidal, single-pin, Ushio. G36T5L	\$35.80
(23443, 00206, 3000312)	
65W Germicidal, single-pin, 155.4 cm, Higuchi. G64T5L	\$83.00

Heat (Infrared)

250W, R40, shatter resistant, GE. 250BR40/HI/TC	\$14.00
--	---------

* Recommend Bulb Removal/Installation Tool (shown below)



Remove and install prefocused/flanged bulbs safely with this simple tool. Push rubber-lined tube directly onto top of bulb and turn clockwise to install, or counterclockwise to remove.

Bulb Removal/Installation Tool. 3283	\$16.00
---	---------

Compressed Air, 10 oz Can (285 g). 77001	\$14.00
Ultra pure duster is nonflammable.	



Continued