



Reticles & Stage Micrometers



Reticles, eyepiece, 1.5 mm thick, Klarmann Rulings, USA

*Specify reticle diameter using 2-digit suffix,
for example, use **KR-207-21** for 21 mm diameter reticle

★ 5 mm/100 Parts. KR-209-__ *	\$80.00
10 mm/100 Parts. KR-207-__ *	\$80.00
Crossline 90°. KR-301-__ *	\$80.00
Crossline 90° with 10 mm/100 Parts Scale. KR-257-__ *	\$80.00
Howard Disc. KR-832-__ *	\$80.00
Indexed Square, 100 sqs, each 1 x 1 mm. KR-406B-__ *	\$80.00
1-10 on x-axis; A-J on y-axis	
Miller Disc	
★ Corner Square. KR-802-__ *	\$80.00
Central Square. KR-802A-__ *	\$80.00
Walton Beckett, G-22. KR-827-1-__ *	\$197.00
Whipple Disc. KR-413-__ *	\$80.00



Stage Micrometers, Transmitted Light

★ 1 mm/100 Parts Scale (0.01 mm), Meiji. MA285	\$105.00
1 mm-0.01 mm/div on 1 mm thick slide, protected with microcoverslip. Preferred for compound microscopes.	
10 mm/100 Parts Scale (0.1 mm), Meiji. MA187	\$105.00
10 mm-0.1 mm/div on 1 mm thick slide, protected with 1 mm thick coverslip. Designed for stereo microscopes.	
25 mm Scale, USA. KR-812	\$180.00
25mm-1 mm/div, 5 mm-0.1 mm/div, 2 mm-0.02 mm/div on 3 mm thick slide. Preferred at magnifications below 400X. Ideal for large field stereo microscopes.	
0.2 mm Square/400 Squares, Meiji. MA287	\$105.00
20 vertical lines, 1 mm long and spaced 0.01 mm apart, intersect 20 similar horizontal lines to create a central 0.2 mm square consisting of 400 squares, each 0.01 mm. 1 mm thick slide is protected with microcoverslip. Useful for x-y calibration.	
0.04"/40 Parts Scale (0.001"), Meiji. MA286	\$105.00
0.04"-0.001"/div on 1 mm thick slide, protected with microcoverslip.	

Stage Micrometer, Reflected Light

1 mm/100 Parts Scale (0.01 mm), Meiji. MA292	\$210.00
1 mm-0.01 mm/div on 3 mm thick metal slide	

Cameras

Digital Cameras



National

D-Moticam 2 Digital Camera	\$299.00 ↓
• 2.0 megapixel (1600 x 1200 pixels) USB 2.0 output with card	
• 1/3" CMOS (Complementary Metal Oxide Semiconductor) chip	
• Adapters for attaching camera to microscope include a C-mount ring and couplers to fit 30-38 mm OD eyepieces	
• 12 mm focusable lens and viewing tube for macro viewing	
• 4-dot calibration slide	
• Motic Images Plus 3.0 imaging software	



National

D-Moticam 5+ Digital Camera	\$586.00
• 5.0 megapixel (2592 x 1944 pixels) USB 3.0 output with card	
• 1/3" CMOS (Complementary Metal Oxide Semiconductor) chip	
• Adapters for attaching camera to microscope include a C-mount ring and couplers to fit 30-38 mm OD eyepieces	
• 12 mm focusable lens and viewing tube for macro viewing	
• 4-dot calibration slide	
• Motic Images Plus 3.0 imaging software	

Digital SLR



Digital SLR Camera, 24.2 MP, Canon. EOS Rebel T8i	\$936.00
Includes SD-16 GB card	
AC Adapter Kit. DRE18	\$67.00