



## Reticles & Stage Micrometers



### Reticles, eyepiece, 1.5 mm thick, Klarmann Rulings, USA

\*Specify reticle diameter using 2-digit suffix,  
for example, use **KR-207-21** for 21 mm diameter reticle

★ 5 mm/100 Parts. <b>KR-209-__</b> *	\$80.00
10 mm/100 Parts. <b>KR-207-__</b> *	\$80.00
Crossline 90°. <b>KR-301-__</b> *	\$80.00
Crossline 90° with 10 mm/100 Parts Scale. <b>KR-257-__</b> *	\$80.00
Howard Disc. <b>KR-832-__</b> *	\$80.00
Indexed Square, 100 sqs, each 1 x 1 mm. <b>KR-406B-__</b> *	\$80.00
1-10 on x-axis; A-J on y-axis	
Miller Disc	
★ Corner Square. <b>KR-802-__</b> *	\$80.00
Central Square. <b>KR-802A-__</b> *	\$80.00
Walton Beckett, G-22. <b>KR-827-1-__</b> *	\$197.00
Whipple Disc. <b>KR-413-__</b> *	\$80.00



### Stage Micrometers, Transmitted Light

★ 1 mm/100 Parts Scale (0.01 mm), Meiji. <b>MA285</b>	\$105.00
1 mm-0.01 mm/div on 1 mm thick slide, protected with microcoverslip. Preferred for compound microscopes.	
10 mm/100 Parts Scale (0.1 mm), Meiji. <b>MA187</b>	\$105.00
10 mm-0.1 mm/div on 1 mm thick slide, protected with 1 mm thick coverslip. Designed for stereo microscopes.	
25 mm Scale, USA. <b>KR-812</b>	\$180.00
25mm-1 mm/div, 5 mm-0.1 mm/div, 2 mm-0.02 mm/div on 3 mm thick slide. Preferred at magnifications below 400X. Ideal for large field stereo microscopes.	
0.2 mm Square/400 Squares, Meiji. <b>MA287</b>	\$105.00
20 vertical lines, 1 mm long and spaced 0.01 mm apart, intersect 20 similar horizontal lines to create a central 0.2 mm square consisting of 400 squares, each 0.01 mm. 1 mm thick slide is protected with microcoverslip. Useful for x-y calibration.	
0.04"/40 Parts Scale (0.001"), Meiji. <b>MA286</b>	\$105.00
0.04"-0.001"/div on 1 mm thick slide, protected with microcoverslip.	

### Stage Micrometer, Reflected Light

1 mm/100 Parts Scale (0.01 mm), Meiji. <b>MA292</b>	\$210.00
1 mm-0.01 mm/div on 3 mm thick metal slide	

## Cameras

### Digital Cameras



#### National

<b>D-Moticam 2</b> Digital Camera	\$299.00 ↓
<ul style="list-style-type: none"> <li>• 2.0 megapixel (1600 x 1200 pixels) USB 2.0 output with card</li> <li>• 1/3" CMOS (Complementary Metal Oxide Semiconductor) chip</li> <li>• Adapters for attaching camera to microscope include a C-mount ring and couplers to fit 30-38 mm OD eyepieces</li> <li>• 12 mm focusable lens and viewing tube for macro viewing</li> <li>• 4-dot calibration slide</li> <li>• Motic Images Plus 3.0 imaging software</li> </ul>	



#### National

<b>D-Moticam 5+</b> Digital Camera	\$586.00
<ul style="list-style-type: none"> <li>• 5.0 megapixel (2592 x 1944 pixels) USB 3.0 output with card</li> <li>• 1/3" CMOS (Complementary Metal Oxide Semiconductor) chip</li> <li>• Adapters for attaching camera to microscope include a C-mount ring and couplers to fit 30-38 mm OD eyepieces</li> <li>• 12 mm focusable lens and viewing tube for macro viewing</li> <li>• 4-dot calibration slide</li> <li>• Motic Images Plus 3.0 imaging software</li> </ul>	

### Digital SLR



<b>Digital SLR Camera, 24.2 MP, Canon. EOS Rebel T8i</b>	\$936.00
Includes SD-16 GB card	
AC Adapter Kit. <b>DRE18</b>	\$67.00

## Accessories



**Adapter, CS-Mount to C-Mount, 5 mm spacer** ..... \$15.00

*Do not use this adapter ring on an existing C-mount camera, as it will negatively impact the quality and parfocality of the imaging system.*



**Couplers, Diagnostic Instruments, USA**

1X C-Mount for 38 mm ISO Photoport. **D10NLC** ..... \$135.00

*Fits newer Nikon and Leica microscopes*

**Couplers, National**

1X C-Mount for 23.2 mm ID Eytube. **930-023** ..... \$29.00

*Lens free - this coupler is not parfocal*

## Magnifiers



A 4X magnifying lens is attached to fine-tipped, stainless steel forceps. Forceps double as a lens support, freeing viewer's hands. Folds for storage.

**Magnifier, Forceps. M15** ..... \$4.00



Lenses are lightweight, water-clear, and shatterproof. Molded handles protect the lens and offer comfort and balance.

**Magnifier, Handle, Ultra-Optix**

2" Lens, 3X. **SV-2P** ..... \$5.45

3" Lens, 2X with 5X bifocal inset. **SV-3P** ..... \$7.45

4" Lens, 2X with 4X bifocal inset. **SV-4P** ..... \$12.55



Provides precise binocular viewing. Visor tilts out of the way when optical glass magnifiers are not needed. Adjustable headband is leather lined. Other magnifications and replacement lenses are available.

**Magnifier, Headband, Optivisor®, USA**

1.75X at 14", DA3. **ELP-503.00** ..... \$51.95

2.5X at 8", DA5. **ELP-505.00** ..... \$51.95

2.75X at 6", DA7. **ELP-507.00** ..... \$53.95

Light Attachment, Quasar™ LED. **ELP-558.00** ..... \$29.95



View small details in low light with this compact, illuminated magnifier. Protective case keeps the lens scratch-free and serves as a handle when lens is extended.

**Magnifier, Lighted, 2 x 2" lens, 2X. SV-BPL** ..... \$7.45