


Vienna Nasal Speculum




Specula

- Graves Vaginal, medium, economy. **3310** \$18.00
- Vienna Nasal, child, 13 cm (5 1/8"). **2722** \$10.50
lightweight, stainless steel

Rongeurs & Bone Instruments

Blumenthal Rongeur



Rongeurs, stainless steel


- Blumenthal, 15 cm (5 3/4")
 - Economy. **410-571** \$42.10
 - Standard. **650-571** \$117.00
- Blumenthal, Slightly Angled Jaws, 15 cm (5 3/4"). **9571** \$37.50
- Mini-Blumenthal, 11.5 cm (4 1/2 "). **410-560** \$54.85
- Friedman, 14 cm (5 1/2")
 - Economy. **410-574** \$42.10
 - Standard. **650-574** \$117.00

Liston Bone Cutter




Bone Cutters, stainless steel, economy

- Liston
 - 14 cm (5 1/2"), 20 mm blades. **320-086** \$24.30
 - 21.5 cm (8 1/2"), 30 mm blades. **320-087 (1095)** \$29.80
- Scissor Style, 10 cm (4"), 18 mm blades. **320-025** \$8.95
straight, strong tips




Bone Nipper, 10 cm (4") long, 6 mm cutting edge. B286..... \$12.00




Bone Shears, Constable, 23 cm (9") long, 50 mm (2") cutting jaws

- German. **965-620** \$160.00



Postmortem Hammer, hook handle, 25 cm (10") long

- Standard. **74-043** \$80.00



16 mm (5/8") cutting blade;
9 cm (3 1/2") long

Skull Breaker, Virchow, economy. 33-1011 \$42.40

Scalpels



Economy

Scalpels

Disposable, plastic handle, sterile

Economy, gray handle	bx/10	cs/50 bx
★ #10 carbon steel blade. 9410	\$9.95	431.00
#11 carbon steel blade. 9411	\$9.95	431.00
#22 carbon steel blade. 9422	\$9.95	431.00

Lance, green handle, England

#10 carbon steel blade, Bx/10. 92410	\$14.50
#11 carbon steel blade, Bx/10. 92411	\$14.50
#21 carbon steel blade, Bx/10. 92421	\$15.50
#22 carbon steel blade, Bx/10. 92422	\$15.50



Miniature, plastic thumb handle, nonsterile, England



Smooth Handle




Fluted Handle

Reusable, one-piece stainless steel, 1 1/2" blade

(1445)

Fluted Handle, narrow blade, 15.5 cm (6"). **14-360**..... \$3.15

(1460)



Scalpel Blades, individually wrapped

Carbon Steel, Economy, Bx/100	\$19.70
#10 blade. 9310	#21 blade. 9321
#11 blade. 9311	#22 blade. 9322
#12 blade. 9312	#23 blade. 9323
#15 blade. 9315	#24 blade. 9324
#20 blade. 9320	
Carbon Steel, Germany, nonsterile	
#60 blade, 25 strips of 4, Bx/100. 2200260	\$55.00
#70 blade, 25 strips of 4, Bx/100. 2200270	\$69.50
Carbon Steel, Swann-Morton, nonsterile	
#60 blade, Bx/50. 22002SM60	\$33.20
Stainless Steel, Economy, Bx/100	\$21.95
#10 blade. 36-010S	#21 blade. 052-610
#11 blade. 052-602	#22 blade. 052-612
#20 blade. 36-020S	
Stainless Steel, Lance, Bx/100	\$45.50
#11 blade. 92111	#15 blade. 92115
Stainless steel, Swann-Morton, sterile	
#22 blade, Bx/100. 22001SM22	\$51.00
Stainless Steel, Economy	
#60 blade, Bx/100. 052-620	\$45.00
#70 blade, Bx/100. 052-622	\$45.00



Safely remove and store up to 100 used blades in this disposable box.

Scalpel Blade Remover. 22048..... \$3.90



Scalpel Handles, reusable

Stainless Steel, solid	
#3 Style, accepts blades #9-16	
Economy. 36-405	\$4.35
Standard. 10-305	\$11.50
#4 Style, accepts blades #20-26	
Economy. 36-410	\$4.35
Standard. 10-310 (565-410)	\$11.50

#7 Style, accepts blades #9-16

Economy. **36-415** \$4.35

Standard. **10-315** (565-415) \$18.50

Stainless Steel, Hollow

 #8 Style, accepts blades #60 & #70. **HAN8** \$8.00

PVC Coated Stainless Steel. *Not autoclavable*

 #3 Style, accepts blades #10-17. **KNF-225.00**..... \$5.35

Plastic, contoured to fit palm of hand

 #5 Style, accepts blades #10-17. **KNF-221.00**..... \$5.00

 #6 Style, accepts blades #20-26. **1605** \$5.00

 #8 Style, accepts blades #60 & #70

 Economy. **HAN8P**
 \$6.00 |

Bard-Parker™, **USA. 371080**..... \$38.25

Scissors




Barraquer

Scissors, Spring, ambidextrous, stainless steel

 Barraquer (butterfly), 6 cm long, economy. **3737** \$34.50

 10 mm angled-on-flat blades



36205/5-254/5-

M3720

 Castroviejo, straight, sharp tips, economy

 Effective cutting edge is 3 mm; tip-to-tip is 3 mm.

 8 cm (3"), 7 mm blades. **5-254** \$45.00

 Effective cutting edge is 6 mm; tip-to-tip is 4 mm.

 8.5 cm (3 1/4"), 5 mm blades. **5-855** \$45.00

 Effective cutting edge is 4 mm; tip-to-tip is 3 mm.

 11.5 cm (4 1/2"), 17 mm blades. **M3720**..... \$26.00


 Effective cutting edge is 12 mm.



Uniband™

Lifetime guarantee against rust and defects.

Featherlite™



 Uniband™, sharp tips, Laschal Surgical, **USA**

 Straight, 12 mm Blades. **LA-4** \$32.00

 Straight, 20 mm Blades. **LA-1** \$32.00

 Straight, 20 mm Serrated Blades. **LA-2** \$32.00

 Curved, 12 mm Blades. **LA-4C** \$32.00

 Curved, 20 mm Blades. **LA-1C** \$32.00

 Featherlite™, Castroviejo, sharp tips, Laschal Surgical, **USA**

 Greater tactile sensitivity than Uniband™ scissors.

 Straight, 12 mm Blades. **MPF-4** \$80.80

 Straight, 20 mm Blades. **MPF-1** \$80.80

 Curved, 12 mm blades. **MPF-4C** \$80.80

Continued

Continued