



Trocar



Trocar, 4-in-1 set, 13 cm, stainless steel. **1596**..... \$17.50

Accessories/Maintenance/Storage



Cutting Boards, white high-density polyethylene
 10 x 6 x 1/2" (25.4 x 15.2 x 1.3 cm). **1054319**..... \$6.95
 18 x 12 x 1/2" (45.7 x 30.5 x 1.3 cm). **105CB1218**..... \$16.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1 1/2 x 3" opening. **DEM-345.00**..... \$80.95



Tools not included.

Dissecting Cases
 6-Loop with snap fastener, black leatherette. **1257370**..... \$1.30
 12-Loop with snap fastener, black leatherette. **1257371**..... \$1.40
 13-Loop with zipper, brown leatherette. **36900**..... \$7.80
 Features red velveteen interior, protector flap, and ID pouch.

Dissecting Kits

*Inquire about customizing a kit for your surgical procedure.
 Kit components are discounted 10%.*

Basic kit. 38501..... \$5.30
 Contents: Dissecting Case, 6-loop (1257370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

Mouse kit. 38551..... \$17.85
 Contents: Dissecting Case, 6-loop (1257370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

Dissecting Pans, stainless steel (without wax)
 20 7/8 x 12 13/16 x 2 1/2". **AM20029**..... \$21.90
 25 3/4 x 17 3/4 x 2 1/4". **AM68357**..... \$64.90



Insect Pins, 38 mm (1 1/2") long, nylon head, Czech Republic


#000, 0.25 mm dia. 355425	#1, 0.40 mm dia. 355440
#00, 0.30 mm dia. 355430	#2, 0.45 mm dia. 355445
#0, 0.35 mm dia. 355435	#3, 0.50 mm dia. 355450
Stainless Steel, Pk/100 \$8.75	
#000, 0.25 mm dia. 355525	#1, 0.40 mm dia. 355540
#00, 0.30 mm dia. 355530	#2, 0.45 mm dia. 355545
#0, 0.35 mm dia. 355535	#3, 0.50 mm dia. 355550



Instrument Cup, silicone, 3 1/2" opening, 2 1/2" deep. **150-1005**..... \$43.00



Instrument Storage Boxes, 10 x 4 x 2", 18-8 stainless steel
 Large, 10 x 4 x 2". **35-408**..... \$35.00



Rulers, 15 cm/6" long, scaled in both metric and English units

Transparent. 4461	\$0.95
White. 4460	\$0.55
Stainless Steel, 12 mm wide, USA. 605-065	\$5.50



Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

Sterilizer, Hot Bead, Germinator™ 500, USA. 5-1450..... \$419.00
 Glass Beads (replacement), 1.5 mm, Bg/10 oz. **5-1455**..... \$21.00

Continued