



**Trocar**



Trocar, 4-in-1 set, 13 cm, stainless steel. **1596**..... \$17.50

**Accessories/Maintenance/Storage**



**Cutting Boards**, white high-density polyethylene  
 10 x 6 x 1/2" (25.4 x 15.2 x 1.3 cm). **1054319**..... \$6.95  
 18 x 12 x 1/2" (45.7 x 30.5 x 1.3 cm). **105CB1218**..... \$16.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1 1/2 x 3" opening. **DEM-345.00**..... \$80.95



Tools not included.

**Dissecting Cases**  
 6-Loop with snap fastener, black leatherette. **1257370**..... \$1.30  
 12-Loop with snap fastener, black leatherette. **1257371**..... \$1.40  
 13-Loop with zipper, brown leatherette. **36900**..... \$7.80  
 Features red velveteen interior, protector flap, and ID pouch.

**Dissecting Kits**

*Inquire about customizing a kit for your surgical procedure.  
 Kit components are discounted 10%.*

**Basic kit. 38501**..... \$5.30  
 Contents: Dissecting Case, 6-loop (1257370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

**Mouse kit. 38551**..... \$17.85  
 Contents: Dissecting Case, 6-loop (1257370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

**Dissecting Pans**, stainless steel (without wax)  
 20 7/8 x 12 13/16 x 2 1/2". **AM20029**..... \$21.90  
 25 3/4 x 17 3/4 x 2 1/4". **AM68357**..... \$64.90



**Insect Pins**, 38 mm (1 1/2") long, nylon head, Czech Republic


|                                      |                                |
|--------------------------------------|--------------------------------|
| #000, 0.25 mm dia. <b>355425</b>     | #1, 0.40 mm dia. <b>355440</b> |
| #00, 0.30 mm dia. <b>355430</b>      | #2, 0.45 mm dia. <b>355445</b> |
| #0, 0.35 mm dia. <b>355435</b>       | #3, 0.50 mm dia. <b>355450</b> |
| Stainless Steel, Pk/100 ..... \$8.75 |                                |
| #000, 0.25 mm dia. <b>355525</b>     | #1, 0.40 mm dia. <b>355540</b> |
| #00, 0.30 mm dia. <b>355530</b>      | #2, 0.45 mm dia. <b>355545</b> |
| #0, 0.35 mm dia. <b>355535</b>       | #3, 0.50 mm dia. <b>355550</b> |



**Instrument Cup**, silicone, 3 1/2" opening, 2 1/2" deep. **150-1005**..... \$43.00



**Instrument Storage Boxes**, 10 x 4 x 2", 18-8 stainless steel  
 Large, 10 x 4 x 2". **35-408**..... \$35.00



**Rulers**, 15 cm/6" long, scaled in both metric and English units  
 Transparent. **4461**..... \$0.95  
 White. **4460**..... \$0.55  
 Stainless Steel, 12 mm wide, **USA. 605-065**..... \$5.50



Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

**Sterilizer, Hot Bead, Germinator™ 500, USA. 5-1450**..... \$419.00  
 Glass Beads (replacement), 1.5 mm, Bg/10 oz. **5-1455**..... \$21.00

*Continued*



Safely incinerate microorganisms in 5-7 seconds at 815°C (1500°F) within this ceramic, stainless steel-encased chamber. Compact and flameless operation allows this unit to be used in anaerobic and aerobic chambers. Ideal for sterilizing loops and needles, glass tube/pipette mouths, and various glass and metal instruments. 1-year warranty.


**Sterilizer, Infrared Heat, BactiZapper™. B1000** ..... \$389.00




- Blue film/white paper
- EO gas/steam
- Latex-free
- Self-sealing

**Sterilization Pouches**

3½ x 9" (8.9 x 22.9 cm), Bx/200. **440166423** ..... \$25.40  
 5¼ x 10" (13.3 x 25.4 cm), Bx/200. **440166424** ..... \$35.95  
 7½ x 13" (19.1 x 33 cm), Bx/200. **440166425** ..... \$62.60



6100A



Tools not included.

4600A

**Sterilizing Trays, plastic, autoclavable, PST, USA**

B/L/M (base, lid & mat) included.

6 x 2½ x ¾" (15.2 x 6.4 x 1.9 cm). **2600A** ..... \$23.00  
 (500252)  
 ★ 6½ x 4 x ¾" (16.2 x 10.2 x 1.9 cm). **4600A** ..... \$33.00  
 10 x 6 x ¾" (25.4 x 15.2 x 1.9 cm). **6100A** ..... \$67.00  
 (500254)




60081736

**Instrument Trays**

Polypropylene, Deep w/Flat Cover, autoclavable

8½ x 6 x 2¼" (22 x 15 x 7 cm). **60081736** ..... \$10.90



Shallow Tray



Deep Trays

**Stainless Steel**

Shallow (Mayo), **USA**

15½ x 10½ x ¾" (38.4 x 27 x 1.9 cm). **15F** ..... \$51.65

Deep w/Dome Cover

8¾ x 3¼ x 1½" (22.2 x 8.4 x 3.8 cm). **35-309** ..... \$31.25

Deep - Order trays and covers separately, **USA**

8¾ x 5 x 2" (22.7 x 12.9 x 5.1 cm). **952** ..... \$28.40

Cover, Fits 952 Tray, dome w/handle. **952D-2** ..... \$44.90



This natural Arkansas stone outlasts and produces a better finish than a man-made stone. Mounted in a wooden box. Dimensions: 4 x 2 x ½"

**Sharpening Stone, fine grit (hard). ABR-800.10** ..... \$24.95

Includes a plastic cleaning basket. Requires only 42 kHz of power. Dimensions: 6 x 3½ x 2"(D).



**Ultrasonic Cleaner, 600 mL, mechanical timer. CLN-520.20** ..... \$45.95  
 Accessories

Cleaning Basket, fine mesh, 5 x 4 x 1¼". **CLN-650.10** ..... \$14.95

Cleaning Solution, 1 Quart. **CLN-850.10** ..... \$9.95  
 Biodegradable concentrate mixes 1:30.