



Trocar



Trocar, 4-in-1 set, 13 cm, stainless steel. **1596**..... \$17.50

Accessories/Maintenance/Storage



Cutting Boards, white high-density polyethylene
 10 x 6 x 1/2" (25.4 x 15.2 x 1.3 cm). **1054319**..... \$6.95
 18 x 12 x 1/2" (45.7 x 30.5 x 1.3 cm). **105CB1218**..... \$16.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1 1/2 x 3" opening. **DEM-345.00**..... \$80.95



Tools not included.

Dissecting Cases
 6-Loop with snap fastener, black leatherette. **1257370**..... \$1.30
 12-Loop with snap fastener, black leatherette. **1257371**..... \$1.40
 13-Loop with zipper, brown leatherette. **36900**..... \$7.80
 Features red velveteen interior, protector flap, and ID pouch.

Dissecting Kits

*Inquire about customizing a kit for your surgical procedure.
 Kit components are discounted 10%.*

Basic kit. 38501..... \$5.30
 Contents: Dissecting Case, 6-loop (1257370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

Mouse kit. 38551..... \$17.85
 Contents: Dissecting Case, 6-loop (1257370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

Dissecting Pans, stainless steel (without wax)
 20 7/8 x 12 13/16 x 2 1/2". **AM20029**..... \$21.90
 25 3/4 x 17 3/4 x 2 1/4". **AM68357**..... \$64.90




Insect Pins, 38 mm (1 1/2") long, nylon head, Czech Republic
 #000, 0.25 mm dia. **355425** #1, 0.40 mm dia. **355440**
 #00, 0.30 mm dia. **355430** #2, 0.45 mm dia. **355445**
 #0, 0.35 mm dia. **355435** #3, 0.50 mm dia. **355450**
 Stainless Steel, Pk/100 \$8.75
 #000, 0.25 mm dia. **355525** #1, 0.40 mm dia. **355540**
 #00, 0.30 mm dia. **355530** #2, 0.45 mm dia. **355545**
 #0, 0.35 mm dia. **355535** #3, 0.50 mm dia. **355550**



Instrument Cup, silicone, 3 1/2" opening, 2 1/2" deep. **150-1005**..... \$43.00



Instrument Storage Boxes, 10 x 4 x 2", 18-8 stainless steel
 Large, 10 x 4 x 2". **35-408**..... \$35.00




Rulers, 15 cm/6" long, scaled in both metric and English units
 Transparent. **4461**..... \$0.95
 White. **4460**..... \$0.55
 Stainless Steel, 12 mm wide, **USA. 605-065**..... \$5.50



Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

Sterilizer, Hot Bead, Germinator™ 500, USA. 5-1450..... \$419.00
 Glass Beads (replacement), 1.5 mm, Bg/10 oz. **5-1455**..... \$21.00

Continued




- Tube Rotator**, programmable, Benchmark (2-year warranty)
- Roto-Bot Mini. **R4040** \$599.00
- Roto-Bot. **R4045** \$899.00
- Speed adjusts from 1-80 rpm
 - Tube capacity (starter set)
 - R4040: 12 x 1.5 mL, 2 x 15 mL & 2 x 50 mL
 - R4045: 30 x 1.5 mL, 8 x 15 mL & 2 x 50 mL
 - Timer is programmable from 1 sec – 99 hrs, 59 min / Cont.

Vortexers



- Vortex Mixers**, Benchmark (2-year warranty)
- MiniMixer, fixed 3000 rpm, blue cup. **VM-03UA** \$138.00
- Vornado™, mini vortexer, fixed 2800 rpm. **BV101**__* \$209.00
- *Specify Gray, Blue or Purple cup head
- BenchMixer™, single cup head, 200-3200 rpm. **BV1000** \$269.00
- Flat Head Attachment, 3". **BV1000-FLAT** \$39.40

Rulers, Meter Stick & Tape Measure



- Rulers**, 15 cm/6" long, scaled in both metric and English units
- Vinyl, flexible
- Transparent. **4461** \$0.95
- White. **4460** \$0.55
- Stainless Steel, 12 mm wide, 0.5 mm thick, **USA. 605-065** \$5.70

- Meter Stick**, lacquered wood, plain ends. **605-050** \$4.90
- One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

- Tape Measure**, 152.4 cm/60", reinforced fiberglass. **11036** \$3.70
- One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.

Sharps Containers



- Sharps Containers**, red polypropylene
- 1 Quart, translucent locking lid. **440065** \$3.80
- Dimensions: 4.25" x 4.25" x 6.25"(H)
- 5.4 Quart, horizontal entry lid, Prevent™, each. **4402269** \$9.90
- Dimensions: 12" x 4.75" x 11"(H)
- 12 Gallon, hinged vertical entry lid, each. **440321391** \$28.95
- Dimensions: 18.25" x 12.75" x 18.75"(H)

Specimen Containers



- Specimen Containers**, polypropylene (PP)
- Screw-Cap
- 90 mL, sterile, Cs/400. **840175** \$110.00
- 480 mL (16 oz), Cs/24. **59014510** \$71.95
- 83 mm (3 1/4") opening
- 960 mL (32 oz), Cs/12. **59014511** \$59.05
- 111 mm (4 3/8") opening
- Snap-Cap, white HDPE bottom/white LDPE lid
- 480 mL (16 oz), Cs/100. **271016W** \$72.50
- 9.9 cm opening x 11.1 cm high
- 1000 mL (34 oz), Cs/100. **271032W** \$109.00
- 11.7 cm opening x 13.7 cm high



- Spot Plates**, glazed white porcelain, 21 mm diameter wells
- 6-Well, each. **LPC-161** \$2.35
- 12-Well, 115 x 95 x 7 mm, each. **3104187701** \$4.15

Novelty Items



Jewelry, sterling silver

Charms

- Caduceus, 7/8" long. **88501** \$11.00
- Chimpanzee, 1 1/16" long. **88503** \$12.00
- Microscope, 3/4" long. **88505** \$11.00
- Mouse, 1" wide. **88507** \$12.00

Earrings, pair

- Stud, microscope, 3/4" long. **88521** \$34.00

Wire Hoops

- Microscope, 3/4" long. **88531** \$25.00
- Mouse, 1" wide. **88533** \$27.00

- Pin, Microscope. **88541** \$17.00

Micro Bible. 89585 \$15.00

View with a 10X objective all 1,245 pages of the King James version of the Holy Bible on this approximately 1" square microfilm.



Ruler, 1 Decimeter, stainless steel, Germany. 605-060 \$7.40 ↓

Ruler is 1 dm long, 1 cm wide, and 1 mm thick. Its volume is 1 cm³ and its weight in grams is the specific gravity. One edge of one side is graduated in centimeters with 1 mm divisions.

Reduce waste & save money with this heavy-duty wringer by completely extracting contents from even the most rigid tube. Features 3" aluminum rollers with steel handle.



Tube-Wringer®, USA. 295401 \$27.45