

**Inverted Microscopes**



CP-5

**Jenco**

**CP-5** Inverted, Trinocular, Brightfield/Phase ..... \$3,389.00

- Widefield 10X/22 eyepieces
- Trinocular, 45° inclined, Siedentopf design viewing head with single diopter; 48-75 mm interpupillary adjustment
- Seven infinity-corrected, plan achromat objectives: brightfield 10X, 20X, and 40X; phase 5X, 10X, 20X, and 40X
- 5-place reverse-facing nosepiece
- Attachable mechanical stage
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- Centered condenser with working distance of 70 mm; iris diaphragm and filter holder; centerable, dedicated phase annuli mounted on sliders
- 9 watt variable LED illuminator
- 56 cm (22")D x 22 cm (8.7")W x 51 cm (20")H
- 5-year mechanical & optical/2-year electrical warranty

**CP-5S** Inverted, Trinocular, Phase..... \$2,889.00

- Widefield 10X/22 eyepieces
- Trinocular, 45° inclined, Siedentopf design viewing head with single diopter; 48-75 mm interpupillary adjustment
- Infinity-corrected, phase plan achromat 5X, 10X, and 20X objectives on 5-place, reverse-facing nosepiece
- Fixed stage measures 22.9 cm (9") x 22.9 cm (9")
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- Centered condenser with working distance of 70 mm; iris diaphragm and filter holder; centerable, dedicated phase annuli mounted on single slider
- 9 watt variable LED illuminator
- 56 cm (22")D x 22 cm (8.7")W x 51 cm (20")H
- 5-year mechanical & optical/2-year electrical warranty



CP-2A1

**Jenco**

**CP-2A1** Inverted, Trinocular, Brightfield/Phase..... \$3,274.00

- Widefield 10X/20 eyepieces
- Trinocular, 30° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
- Six DIN, LWD plan achromat objectives: brightfield 10X, 25X, and 40X; phase 10X, 25X, and 40X
- 5-place, reverse-facing nosepiece
- Attachable mechanical stage includes three inserts
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm and filter holder; working distance of 55 mm for phase contrast, or 125 mm if swung out for brightfield
- Centerable dedicated phase annuli mounted on sliders
- Centerable 30 watt halogen bulb; Koehler-type illuminator
- 48.25 cm (19")D x 21.6 cm (8½")W x 47 cm (18½")H
- 5-year mechanical & optical/2-year electrical warranty



BC-366

**Jenco**

**BC-366** Inverted, Trinocular, BF/Phase/Fluor ..... \$10,349.00

- Widefield 10X/22 eyepieces
- Binocular, 45° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
- Side camera port
- Six infinity, LWD plan achromat objectives: brightfield/fluorescence 10X, 25X, and 40X; phase 10X, 25X, and 40X
- 6-place, side-facing nosepiece
- Attachable mechanical stage
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm and filter holder; centerable, dedicated phase annuli mounted in turret
- 30 watt halogen transmitted light; Koehler-type illuminator
- 100 watt mercury epi-fluorescence
- Blue (FITC) and green (TRITC) fluorescence filter sets
- 5-year mechanical & optical/2-year electrical warranty

**Meiji** ..... Access [www.meijitechno.com](http://www.meijitechno.com) and call us to place an order