



Safely incinerate microorganisms in 5-7 seconds at 815°C (1500°F) within this ceramic, stainless steel-encased chamber. Compact and flameless operation allows this unit to be used in anaerobic and aerobic chambers. Ideal for sterilizing loops and needles, glass tube/pipette mouths, and various glass and metal instruments. 1-year warranty.


Sterilizer, Infrared Heat, BactiZapper™. B1000 \$389.00




- Blue film/white paper
- EO gas/steam
- Latex-free
- Self-sealing

Sterilization Pouches

3½ x 9" (8.9 x 22.9 cm), Bx/200. **440166423** \$25.40
 5¼ x 10" (13.3 x 25.4 cm), Bx/200. **440166424** \$35.95
 7½ x 13" (19.1 x 33 cm), Bx/200. **440166425** \$62.60



6100A



Tools not included.

4600A

Sterilizing Trays, plastic, autoclavable, PST, USA

B/L/M (base, lid & mat) included.

6 x 2½ x ¾" (15.2 x 6.4 x 1.9 cm). **2600A** \$23.00
 (500252)
 ★ 6½ x 4 x ¾" (16.2 x 10.2 x 1.9 cm). **4600A** \$33.00
 10 x 6 x ¾" (25.4 x 15.2 x 1.9 cm). **6100A** \$67.00
 (500254)




60081736

Instrument Trays

Polypropylene, Deep w/Flat Cover, autoclavable

8½ x 6 x 2¼" (22 x 15 x 7 cm). **60081736** \$10.90



Shallow Tray



Deep Trays

Stainless Steel

Shallow (Mayo), **USA**

15½ x 10½ x ¾" (38.4 x 27 x 1.9 cm). **15F** \$51.65

Deep w/Dome Cover

8¾ x 3¼ x 1½" (22.2 x 8.4 x 3.8 cm). **35-309** \$31.25

Deep - Order trays and covers separately, **USA**

8¾ x 5 x 2" (22.7 x 12.9 x 5.1 cm). **952** \$28.40

Cover, Fits 952 Tray, dome w/handle. **952D-2** \$44.90



This natural Arkansas stone outlasts and produces a better finish than a man-made stone. Mounted in a wooden box. Dimensions: 4 x 2 x ½"

Sharpening Stone, fine grit (hard). ABR-800.10 \$24.95

Includes a plastic cleaning basket. Requires only 42 kHz of power. Dimensions: 6 x 3½ x 2"(D).



Ultrasonic Cleaner, 600 mL, mechanical timer. CLN-520.20 \$45.95
 Accessories

Cleaning Basket, fine mesh, 5 x 4 x 1¼". **CLN-650.10** \$14.95

Cleaning Solution, 1 Quart. **CLN-850.10** \$9.95

Biodegradable concentrate mixes 1:30.