


Accessories



Bibulous Paper, 10 x 15 cm, Bk/50 sheets. **BIB-46-50** bk \$4.10 cs/48 bk 190.00


Forceps

Spade Tips




Angled (Kuehne), 11 cm (4 1/4")
 Stainless Steel, economy. **33-1041** \$4.15
 Teflon® Coated, Bel-Art. **379420000** \$22.50
 Straight, stainless steel, 10.5 cm, economy. **57.926** \$5.85

Square Point Tips, straight



Stainless Steel, 12 cm (4 3/4"), economy. **1244** \$2.25
 Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379430000** \$24.95



Minimize contamination during slide & coverslip handling. Will not scratch; corrosion and heat resistant.

Hemostat, Polypropylene, 12.5 cm. **HEM-300.00** \$1.50

Cytology Funnels

Single,
Metal Clips
Required



Double,
All Plastic



White filter cards will accept up to 0.5 mL of sample, while tan filters will accept volumes up to 0.4 mL.

Cytology Funnels, Simport CytoSep™

Shandon (Thermo) Cytospin® CytoCentrifuges

Single. *Metal slide clips required*

★ White Filter & Cap, tray/50, Cs/500. **M964-10FW** \$402.00
 (59 91 040, 110 25 48)

White Filter & Cap, Pk/1, Cs/500. **M964-10FW1** ★ \$530.00

Tan Filter & Cap, tray/50, Cs/500. **M964-10FT** \$402.00

Single. *All plastic - metal slide clips not required*

White Filter & Cap, Cs/480. **M965-10FW** \$583.00
 (A787 100003)

Tan Filter & Cap, Cs/480. **M965-10FT** \$583.00

Double. *All plastic - metal slide clips not required*

White Filter & Cap, Cs/480. **M965-20FW** \$762.00
 (A787 100005)

Immersion Oil & Bottles



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities.
USA

Immersion Oil

Type A, low viscosity (150 cSt). *Short W.D. optics*

7.5 mL in LDPE dispensing bottle. **18101** \$10.00

20 mL in HDPE dispensing bottle. **18102** \$11.55

480 mL in amber bottle. **18104** \$68.25

Type B, high viscosity (1250 cSt). *Long W.D. optics*

7.5 mL in LDPE dispensing bottle. **18105** \$10.00

20 mL in HDPE dispensing bottle. **18106** \$11.55

480 mL in amber bottle. **18108** \$68.25

Type NVH, very high viscosity (21,000 cSt)

480 mL in amber bottle. **18113** \$67.00


Type HF, low fluorescence (700 cSt)

480 mL in amber bottle. **18118** \$127.00

Type FF, fluorescence free (170 cSt)

480 mL in amber bottle. **18124** \$105.50

A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



Bottle, Dispensing, HDPE, 20 mL, each. **18145** \$7.15



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.

Bottle, Balsam, glass, 45 mL (empty), **USA. 616** \$18.75

Beakers



Polypropylene Beakers



Premiu

Beakers

Glass (borosilicate), Griffin low form, single scale, economy

10 mL, no graduations, each. 6102022910	\$1.85
50 mL, each. 7-440050	\$3.10
100 mL, each. 61020229100	\$3.40
150 mL	
Pk/12. 61020229150	\$46.30
Cs/48. 61020229150	\$156.00
250 mL, each. 61020229250	\$3.50
400 mL, each. 61020229400	\$4.40
600 mL, each. 61020229600	\$5.15
800 mL	
Pk/6. 600BG800	\$65.80
Cs/24. 600BG800	\$236.00
1000 mL, each. 600BG1000	\$9.10
1000 mL, Heavy-Duty, each. 600BG1003	\$14.75
2000 mL, each. 600BG2000	\$16.65
4000 mL, each. 600BG4000	\$46.25

Polypropylene

Economy, Triple-Spout

50 mL, 5 mL graduations, Cs/100. 3640	\$25.90
100 mL, 10 mL graduations, Cs/100. 80093	\$39.40
250 mL, 20 mL graduations, Cs/100. 3642	\$58.20
1000 mL, 50 mL graduations, Cs/100. 80098	\$129.00

Premium, Griffin low form, single scale

50 mL, 10 mL graduations, each. 262100000	\$3.20
100 mL, 25 mL graduations, each. 262110000	\$3.05
250 mL, 50 mL graduations, each. 262130000	\$5.40
600 mL, 100 mL graduations, each. 262150000	\$8.50
1000 mL, 200 mL graduations, each. 262160000	\$8.60
2000 mL, 500 mL graduations, each. 262170000	\$23.15
4000 mL, 1000 mL graduations, each. 262190000	\$45.10

Bench Protectors



Bench Liner, Roll

Protect your bench surface using this 40 mil absorbent liner backed with 0.5 mil polyethylene film.

Bench Liners, white

Roll, 50' x 20", Labmat™. 246750000	\$70.40
Sheet, Ahlstrom, 46 x 57 cm (18 x 22½"), Pk/50. 4002-4657	\$86.90

Bottles



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.

Bottle, Balsam, glass, 45 mL (empty), **USA. 616**\$18.75




Deliver controlled drops from standard 0.030" tip opening with these LDPE bottles. Sealer cap snaps on easily.

Bottles, Dispensing, squeeze bottles

LDPE (low-density polyethylene), Yorker cap

30 mL (1 oz), Bg/12. 116370001	\$19.35
60 mL (2 oz), Bg/12. 116370002	\$22.50
125 mL (4 oz), Bg/12. 116370004	\$37.40
185 mL (6 oz), Bg/12. 116370006	\$37.40

A long nozzle, textured sides, and slim profile make this chemically inert, HDPE container ideal for dispensing immersion oil.



HDPE (high-density polyethylene)

20 mL, each. 18145	\$7.15
---------------------------------	--------

Amber glass protects light-sensitive liquids. Screw-cap dropper features a black rubber bulb and glass pipet.



Boston Round

Bottles, Dropping, Amber Glass (flint)

30 mL (1 oz), Bx/12. 600BDPSC30	\$14.50
Bottle Dropper Assembly (replacement)	
20/400, 7 x 76 mm Pipet, Pk/12. 3104053D	\$9.95
60 mL (2 oz), Bx/12. 600BDPSC60	\$16.70
Bottle Dropper Assembly (replacement)	
20/400, 7 x 89 mm Pipet, Pk/12. 3104054D	\$9.85

French Square



Bottles, Dropping, Clear Glass (flint)

15 mL (1/2 oz), Bx/12. 600BDFSQ15	\$17.55
Bottle Dropper Assembly (replacement)	
20/400, 7 x 60 mm Pipet, Pk/12. 3104074	\$10.85
30 mL (1 oz), Bx/12. 600BDFSQ30	\$19.00
Bottle Dropper Assembly (replacement)	
24/400, 7 x 75 mm Pipet, Pk/12. 3104075	\$13.70

Continued