

Slide & Cassette Pens, Marker & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker.
Pilot pens are ideal for cytologists and pathologists.



| Pens, Pilot®, permanent, extra fine point | ea | bx/12 |
|---|--------|-------|
| Black. 44102 | \$1.70 | 17.45 |
| Blue. 44103 | \$1.70 | 17.45 |
| Green. 44105 | \$1.70 | 17.45 |
| Red. 44104 | — | 17.45 |



The 0.75 mm tip is ideal for use on small surfaces, e.g., histology cassettes, microscope slides and cryogenic vials. The permanent ink is waterproof and resistant to most laboratory chemicals and solvents.

Marker, Diamond PerfectMark™, black, Bx/10. **1000**..... \$24.90





| Scribes, retractable | | |
|--|---------|--|
| Carbide-Tipped (sharpenable). RSC | \$15.30 | |
| Carbide Tip (replacement). RCT-1 | \$6.20 | |
| Diamond-Tipped. RSD | \$22.75 | |
| Diamond Tip (replacement). RDT-1 | \$12.90 | |



Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.

| Slide Warmers | |
|---|----------|
| Small, 26 x 18 cm (10 x 7"). XH-2002 | \$273.00 |
| Large, 63 x 20 cm (25 x 8"). XH-2001 | \$364.00 |

Staining



Ten slide compartments, separated by ½" of space, are additionally isolated by eight barrier dividers when the lid is closed. Individual slides are suspended ½" above the chamber floor by four pedestal posts and four corner posts.

Staining/Moisture Chambers

| 10-Slide, polystyrene, USA | ea | pk/3 |
|-----------------------------------|---------|--------|
| Clear. 240-9000-010 | \$54.95 | 154.00 |
| Black. 240-9020-Z10 | \$56.10 | 157.00 |



| 20-Slide, Simport | |
|---|----------|
| Clear Lid (PETG), black base (ABS), each. M920-1 | \$173.00 |
| Black Lid (ABS), black base (ABS), each. M920-2 | \$173.00 |



Staining Dishes

| Glass (flint) | |
|--|---------|
| 20-Slide, holds slides horizontally in metal rack Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H) Set (dish, rack & lid). SDS-20 | \$27.25 |
| Dish & Lid. SD-20 | \$16.60 |
| Rack only, stainless steel. SR-20 | \$13.65 |
| 30-Slide, holds slides vertically in metal rack Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H) Set (dish, rack & lid). SDS-30 | \$37.20 |
| Dish & Lid. SD-30 | \$18.50 |
| Rack only, stainless steel. SR-30 | \$22.85 |
| 60-Slide, holds slides vertically in metal rack Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H) Set (dish, rack & lid). SDS-60 | \$52.00 |
| Dish & Lid. SD-60 | \$28.45 |
| Rack only, stainless steel. SR-60 | \$29.35 |

Continued