

Slide & Cassette Pens, Marker & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker.
Pilot pens are ideal for cytologists and pathologists.



Pens, Pilot®, permanent, extra fine point	ea	bx/12
Black. 44102	\$1.70	17.45
Blue. 44103	\$1.70	17.45
Green. 44105	\$1.70	17.45
Red. 44104	—	17.45



The 0.75 mm tip is ideal for use on small surfaces, e.g., histology cassettes, microscope slides and cryogenic vials. The permanent ink is waterproof and resistant to most laboratory chemicals and solvents.

Marker, Diamond PerfectMark™, black, Bx/10. **1000**..... \$24.90





Scribes, retractable		
Carbide-Tipped (sharpenable). RSC	\$15.30	
Carbide Tip (replacement). RCT-1	\$6.20	
Diamond-Tipped. RSD	\$22.75	
Diamond Tip (replacement). RDT-1	\$12.90	



Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.

Slide Warmers		
Small, 26 x 18 cm (10 x 7"). XH-2002	\$273.00	
Large, 63 x 20 cm (25 x 8"). XH-2001	\$364.00	

Staining



Ten slide compartments, separated by ½" of space, are additionally isolated by eight barrier dividers when the lid is closed. Individual slides are suspended ½" above the chamber floor by four pedestal posts and four corner posts.

Staining/Moisture Chambers

10-Slide, polystyrene, USA	ea	pk/3
Clear. 240-9000-010	\$54.95	154.00
Black. 240-9020-Z10	\$56.10	157.00




20-Slide, Simport		
Clear Lid (PETG), black base (ABS), each. M920-1	\$173.00	
Black Lid (ABS), black base (ABS), each. M920-2	\$173.00	



Staining Dishes


Glass (flint)		
20-Slide, holds slides horizontally in metal rack		
Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H)		
Set (dish, rack & lid). SDS-20	\$27.25	
Dish & Lid. SD-20	\$16.60	
Rack only, stainless steel. SR-20	\$13.65	
30-Slide, holds slides vertically in metal rack		
Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H)		
Set (dish, rack & lid). SDS-30	\$37.20	
Dish & Lid. SD-30	\$18.50	
Rack only, stainless steel. SR-30	\$22.85	
60-Slide, holds slides vertically in metal rack		
Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H)		
Set (dish, rack & lid). SDS-60	\$52.00	
Dish & Lid. SD-60	\$28.45	
Rack only, stainless steel. SR-60	\$29.35	

Continued



Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

Plastic (polyoxymethylene [POM]), holds 25 slides
Set (dish, rack & lid). **0990-99** \$50.90



Staining Jars, hold slides (up to 76 x 26 mm) vertically

Coplin, holds 10 slides back-to-back	ea	cs/12
Glass w/Ground Glass Cover, USA. 673	\$11.85	134.00
Glass Cover, Fits 673 Jar. 673-LID	\$4.40	—
Glass w/Screw-On Phenolic Cap, USA 674	\$11.85	134.00
Phenolic Cap, Fits 674 Jar, black. 674-LID ...	\$4.40	—
PP w/Screw-On Cap. 60062101	\$4.75	47.05
Polypropylene jar is unbreakable, microwave proof, stain resistant, and chemically inert. Cap is domed.		
Wide Top, Glass w/Square Lid, USA 675	\$15.95	179.00
Holds 16 slides back-to-back		
Glass Lid, Fits 675 Jar. 675-LID	\$4.85	—



Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

Wide Top, Plastic (POM)

12-Slide, EasyDip™, Simport	ea	cs/6
Jar only, Blue. M900-12B	\$11.20	61.75
Jar only, Green. M900-12G	\$11.20	61.75
Jar only, Pink. M900-12P	\$11.20	61.75
Jar only, White. M900-12W	\$11.20	61.75
Jar only, Yellow. M900-12Y	\$11.20	61.75
Rack/Lid only, gray. M905-12DGY	\$12.45	68.25





Anodized aluminum staining jar holder

Staining Jar Holder, 6-Place. **M906** \$75.20



24-Slide, light green, microwaveable, acetal Jar & Lid. **4501030** \$11.55
Rack only, gray. **4501042** \$18.90
Staining Jar Station, 12-Place. **4501011** \$221.00
Includes metal frame, (12) 24-slide dishes w/ lids (4501030) & (1) rack (4501042).
Dimensions: 29(L) x 5(W) x 5"(H)



Tissue (Embedding) & Biopsy Cassettes



Base Molds

15 x 15 x 5 mm, Cs/1000. M475-2	\$140.00
24 x 24 x 5 mm, Cs/1000. M475-3	\$140.00
★ 30 x 24 x 5 mm, Cs/1000. M475-4	\$140.00
37 x 24 x 5 mm, Cs/1000. M475-5	\$140.00
37 x 24 x 10 mm, Cs/500. M475-10	\$72.65
Designed for Simport Macrosette™.	
Metal, reusable	
15 x 15 x 5 mm, Pk/12. HB-15	\$24.60
24 x 24 x 5 mm, Pk/12. HB-24	\$24.60
30 x 24 x 5 mm, Pk/12. HB-30	\$24.60
37 x 24 x 5 mm, Pk/12. HB-37	\$24.60

Solvent resistant, 0.2 x 0.2 mm mesh.




Biopsy Bags

Nylon, 30 x 45 mm, 5 bags/200, Pk/1000. 4502011	\$235.00
Paper (Tea Bags), 2¼ x 2¼", Pk/500. 56100	\$40.55

M476-1

Made of 60 PPI cellular polyester urethane foam. Use rectangles with cassettes, circles with capsules.



M476-2

Biopsy Pads

Rectangular	
Black Foam, 30.2 x 25.4 mm, Bg/1000. 339014	\$42.00
Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. M476-1	\$201.00
Circular, Simport	
Blue Foam, 25.4 x 2.7 mm, Cs/10,000. M476-2	\$244.00