


Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

Plastic (polyoxymethylene [POM]), holds 25 slides
Set (dish, rack & lid). **0990-99** \$50.90



Staining Jars, hold slides (up to 76 x 26 mm) vertically

| | | |
|--|---------|--------|
| Coplin, holds 10 slides back-to-back | ea | cs/12 |
| Glass w/Ground Glass Cover, USA. 673 | \$11.85 | 134.00 |
| Glass Cover, Fits 673 Jar. 673-LID | \$4.40 | — |
| Glass w/Screw-On Phenolic Cap, USA 674 | \$11.85 | 134.00 |
| Phenolic Cap, Fits 674 Jar, black. 674-LID ... | \$4.40 | — |
| PP w/Screw-On Cap. 60062101 | \$4.75 | 47.05 |
| Polypropylene jar is unbreakable, microwave proof, stain resistant, and chemically inert. Cap is domed. | ea | cs/12 |
| Wide Top, Glass w/Square Lid, USA 675 | \$15.95 | 179.00 |
| Holds 16 slides back-to-back | | |
| Glass Lid, Fits 675 Jar. 675-LID | \$4.85 | — |



Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

Wide Top, Plastic (POM)

| | | |
|--|---------|-------|
| 12-Slide, EasyDip™, Simport | ea | cs/6 |
| Jar only, Blue. M900-12B | \$11.20 | 61.75 |
| Jar only, Green. M900-12G | \$11.20 | 61.75 |
| Jar only, Pink. M900-12P | \$11.20 | 61.75 |
| Jar only, White. M900-12W | \$11.20 | 61.75 |
| Jar only, Yellow. M900-12Y | \$11.20 | 61.75 |
| Rack/Lid only, gray. M905-12DGY | \$12.45 | 68.25 |





Anodized aluminum staining jar holder

Staining Jar Holder, 6-Place. **M906** \$75.20



24-Slide, light green, microwaveable, acetal
Jar & Lid. **4501030** \$11.55
Rack only, gray. **4501042** \$18.90
Staining Jar Station, 12-Place. **4501011** \$221.00
Includes metal frame, (12) 24-slide dishes
w/ lids (4501030) & (1) rack (4501042).
Dimensions: 29(L) x 5(W) x 5"(H)



Tissue (Embedding) & Biopsy Cassettes



Base Molds

| | |
|--|----------|
| 15 x 15 x 5 mm, Cs/1000. M475-2 | \$140.00 |
| 24 x 24 x 5 mm, Cs/1000. M475-3 | \$140.00 |
| ★ 30 x 24 x 5 mm, Cs/1000. M475-4 | \$140.00 |
| 37 x 24 x 5 mm, Cs/1000. M475-5 | \$140.00 |
| 37 x 24 x 10 mm, Cs/500. M475-10 | \$72.65 |
| Designed for Simport Macrosette™. | |
| Metal, reusable | |
| 15 x 15 x 5 mm, Pk/12. HB-15 | \$24.60 |
| 24 x 24 x 5 mm, Pk/12. HB-24 | \$24.60 |
| 30 x 24 x 5 mm, Pk/12. HB-30 | \$24.60 |
| 37 x 24 x 5 mm, Pk/12. HB-37 | \$24.60 |

Solvent resistant,
0.2 x 0.2 mm mesh.




Biopsy Bags

| | |
|--|----------|
| Nylon, 30 x 45 mm, 5 bags/200, Pk/1000. 4502011 | \$235.00 |
| Paper (Tea Bags), 2¼ x 2¼", Pk/500. 56100 | \$40.55 |

M476-1

Made of 60 PPI cellular polyester urethane foam. Use rectangles with cassettes, circles with capsules.



M476-2

Biopsy Pads

| | |
|---|----------|
| Rectangular | |
| Black Foam, 30.2 x 25.4 mm, Bg/1000. 339014 | \$42.00 |
| Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. M476-1 | \$201.00 |
| Circular, Simport | |
| Blue Foam, 25.4 x 2.7 mm, Cs/10,000. M476-2 | \$244.00 |